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SURGICAL
ACCESSORY
English

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Imagem ilustrativa	Código	Modelo	Dimensões Principais (mm)
	018-15-F40-00-13	STOP DRILL 3MM	Ø 0,5 x 3,0
	018-15-F40-00-14	STOP DRILL 5MM	Ø 0,5 x 5,0
	018-15-F40-00-15	STOP DRILL 7MM	Ø 0,5 x 7,0
	018-15-F70-00-01	STOP DRILL Ø3,5X5MM	Ø 3,5 x 5,0
	018-15-F70-00-02	STOP DRILL Ø3,5X6MM	Ø 3,5 x 6,0
	018-15-F70-00-03	STOP DRILL Ø3,5X7MM	Ø 3,5 x 7,0
	018-15-F70-00-04	STOP DRILL Ø3,5X8MM	Ø 3,5 x 8,0
	018-15-F70-00-05	STOP DRILL Ø3,5X10MM	Ø 3,5 x 10,0
	018-15-F70-00-06	STOP DRILL Ø3,5X11,5MM	Ø 3,5 x 11,5
	018-15-F70-00-07	STOP DRILL Ø3,5X13MM	Ø 3,5 x 13,0
	018-15-F70-00-08	STOP DRILL Ø3,5X15MM	Ø 3,5 x 15,0
	018-15-F70-00-09	STOP DRILL Ø4,5X5MM	Ø 4,5 x 5,0
	018-15-F70-00-10	STOP DRILL Ø4,5X6MM	Ø 4,5 x 6,0
	018-15-F70-00-11	STOP DRILL Ø4,5X7MM	Ø 4,5 x 7,0
	018-15-F70-00-12	STOP DRILL Ø4,5X8MM	Ø 4,5 x 8,0
	018-15-F70-00-13	STOP DRILL Ø4,5X10MM	Ø 4,5 x 10,0
	018-15-F70-00-14	STOP DRILL Ø4,5X11,5MM	Ø 4,5 x 11,5
	018-15-F70-00-15	STOP DRILL Ø4,5X13MM	Ø 4,5 x 13,0
	018-15-F70-00-16	STOP DRILL Ø4,5X15MM	Ø 4,5 x 15,0
	018-15-F70-00-17	STOP DRILL Ø6,5X5MM	Ø 6,5 x 5,0
	018-15-F70-00-18	STOP DRILL Ø6,5X6MM	Ø 6,5 x 6,0
	018-15-F70-00-19	STOP DRILL Ø6,5X7MM	Ø 6,5 x 7,0
	018-15-F70-00-20	STOP DRILL Ø6,5X8MM	Ø 6,5 x 8,0
	018-15-F70-00-21	STOP DRILL Ø6,5X10MM	Ø 6,5 x 10,0
	018-15-F70-00-22	STOP DRILL Ø6,5X11,5MM	Ø 6,5 x 11,5
	018-15-F70-00-23	STOP DRILL Ø6,5X13MM	Ø 6,5 x 13,0
	018-15-F70-00-24	STOP DRILL Ø6,5X15MM	Ø 6,5 x 15,0
	018-15-F50-01-76	STOP DRILL GUIADA Ø3,0X5	Ø 3,0 x 5,0
	018-15-F50-01-77	STOP DRILL GUIADA Ø3,0X6	Ø 3,0 x 6,0
	018-15-F50-01-78	STOP DRILL GUIADA Ø3,0X7	Ø 3,0 x 7,0
	018-15-F50-01-79	STOP DRILL GUIADA Ø3,0X8	Ø 3,0 x 8,0
018-15-F50-01-80	STOP DRILL GUIADA Ø3,0X10	Ø 3,0 x 10,0	
018-15-F50-01-81	STOP DRILL GUIADA Ø3,0X11,5	Ø 3,0 x 11,5	
018-15-F50-01-82	STOP DRILL GUIADA Ø3,0X13	Ø 3,0 x 13,0	
018-15-F50-01-83	STOP DRILL GUIADA Ø3,0X15	Ø 3,0 x 15,0	
018-15-F50-01-84	STOP DRILL GUIADA Ø5,0X5	Ø 5,0 x 5,0	
018-15-F50-01-85	STOP DRILL GUIADA Ø5,0X6	Ø 5,0 x 6,0	
018-15-F50-01-86	STOP DRILL GUIADA Ø5,0X7	Ø 5,0 x 7,0	
018-15-F50-01-87	STOP DRILL GUIADA Ø5,0X8	Ø 5,0 x 8,0	
018-15-F50-01-88	STOP DRILL GUIADA Ø5,0X10	Ø 5,0 x 10,0	
018-15-F50-01-89	STOP DRILL GUIADA Ø5,0X11,5	Ø 5,0 x 11,5	
018-15-F50-01-90	STOP DRILL GUIADA Ø5,0X13	Ø 5,0 x 13,0	
018-15-F50-01-91	STOP DRILL GUIADA Ø5,0X15	Ø 5,0 x 15,0	

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SURGICAL ACCESSORY

INSTRUCTIONS FOR USE

TECHNICAL NAME: Dental Surgical Instruments

TRADE NAME: Surgical Accessory

MODEL: For a detailed description of the products, see tables 1 to 7.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

Carefully read these Instructions for Use prior to using the product. Follow the recommendations, warnings and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications. Product for Medical and Dental Use. Product Supplied without Sterilization. Sterilize before use. Valid for an Indefinite Period of Time. Reusable and Resterilizable Product.

ANVISA REGISTRATION NO.: 81684340004.

CODE, BATCH NO., AND DATE OF MANUFACTURE: See product labeling.

Rev. 05

SUMMARY

1. DETAILED PRODUCT DESCRIPTION	3
2. PRODUCT COMPOSITION	5
3. PRODUCT PRESENTATION	5
4. PRODUCT RECOMMENDATIONS, PURPOSE OR INTENDED USE	5
5. PRINCIPLES OF INTENDED PRODUCT FUNCTION AND PERFORMANCE	5
6. INSTRUCTIONS FOR USE OF THE PRODUCT	6
7. PRECAUTIONS, RESTRICTIONS AND WARNINGS	6
8. CONTRAINDICATIONS	7
9. POTENTIAL ADVERSE EFFECTS	7
10. FINISHING AND SURFACE TREATMENTS	7
11. DESCRIPTION OF SYMBOLS	7
12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS	7
13. STERILIZATION	8
14. EXPIRATION DATE	8
15. IDENTIFICATION AND TRACEABILITY	8
16. CONDITIONS OF STORAGE, TRANSPORT, CONSERVATION AND HANDLING	8
17. DISPOSING OF THE PRODUCT	9
18. CUSTOMER SERVICE	9
19. LIMITED WARRANTY	9

1. DETAILED PRODUCT DESCRIPTION

The various models of the Surgical Accessory are recommended for use in dental procedures for the purpose of assisting the dental surgeon in accessing and preparing bone tissue, as well as during the support procedures for drilling, preparation of pilot holes, and adjustment procedures used in techniques of installation of Plenum Dental Implants and their respective prosthetic components in order to restore masticatory function and aesthetic appearance of patients' oral cavity.

The Surgical Accessory line consists of a family of metallic, non-implantable, non-cutting surgical accessories developed for application in dental procedures for the purpose of assisting dental surgeons in support and complement processes involved in the technique of installation of Plenum® Dental Implants.

The Surgical Accessory is a reusable product, supplied in non-sterile conditions, and must be sterilized before use (see topic Sterilization). The product is available in the models listed below.

Surgical Accessory - Depth Probe: the probe is a surgical accessory with a cylindrical blunt tip end used to assist in the preparation of bone tissue in relation to the lateral walls and depth. The probe has graduated height markings allowing the body's height index to be observed.

At its other flattened end it is possible to identify two 7 and 5 mm dimensions that are used to verify the space between teeth.

The probes are available in the codes, models and measurements listed in table 1 below.

TABLE 1: Illustrative image of probes, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F10-00-90	Ø2.8 DEPTH PROBE	Ø 2,80
	018-15-F10-00-91	Ø1.6 DEPTH PROBE	Ø 1,60

Surgical Accessory - Alignment Pin: the alignment pin is a cylindrical surgical accessory fitted with 1.6 and 2.1 mm diameters with height markings used to assist in directing drilling during preparation of the bone tissue in relation to the implants and or lateral teeth.

Alignment pins are available in the codes, models and dimensions listed in table 2 below.

TABLE 2: Illustrative images of alignment pins, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F10-00-92	ALIGNMENT PIN	Ø 1,60 x 2,1
	018-15-F10-00-93	25° ALIGNMENT PIN	Ø 1,60 ´ 2,1

Surgical Accessory - Depth gauge: the depth gauge is a cylindrical surgical accessory equipped with height markings of 1 by 1 mm. It is recommended that the depth gauge be fastened to the inside of the Plenum Dental Implant so that professionals are able to verify the height of the transmucosal tissue for improved prosthetic planning. The depth gauge is supplied in 3 types of prosthetic seatings: SL, RE and SH.

Depth gauges are available in the codes, models and dimensions listed in table 3 below.

TABLE 3: Illustrative images of depth gauges, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F10-00-94	DEPTH GAUGE SL	Ø 2.0
	018-15-F10-00-99	DEPTH GAUGE RE	Ø 2.5
	018-15-F10-01-00	DEPTH GAUGE SH	Ø 2.5
	018-15-F10-02-82	DEPTH GAUGE UF	Ø 2.5
	018-15-F10-02-83	DEPTH GAUGE MF	Ø 2.5

Surgical Accessory - Plenum Flag Guide: the Plenum Flag Guide is a surgical accessory with a 1.6 mm rod recommended for use in assisting in the preparation of bone tissue in relation to lateral distances with dimensional reference markings: 5, 7, 9 and 11 mm. The Plenum Flag Guide is available in the code and dimension listed in table 4 below.

TABLE 4: Illustrative image of the Plenum Flag Guide, code and dimension.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F10-00-95	PLENUM FLAG GUIDE	Ø 1.60

Surgical Accessory - Guided Drill Guide-Pilot: the Guided Drill Guide is a cylindrical surgical accessory available in two models with different external diameters and several internal diameters. The Guided Drill Guide is recommended for use in the application of guided surgery techniques in order to properly guide Plenum drills during the preparation and drilling of bone tissue. Guided Drill Guides are available in the codes, models and dimensions listed in table 5 below.

TABLE 5: Illustrative image of pilot drill guides, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F50-00-38	GUIDED DRILL GUIDE-PILOT	Ø 2.0
	018-15-F50-00-39	Ø1.65 GUIDED DRILL GUIDE	Ø 1.65
	018-15-F50-00-40	Ø2.00 GUIDED DRILL GUIDE	Ø 2.0
	018-15-F50-00-41	Ø2.15 GUIDED DRILL GUIDE	Ø 2.15
	018-15-F50-00-42	Ø2.85 GUIDED DRILL GUIDE	Ø 2.85
	018-15-F50-00-43	Ø3.35 GUIDED DRILL GUIDE	Ø 3.35
	018-15-F50-00-44	Ø3.85 GUIDED DRILL GUIDE	Ø 3.85
	018-15-F50-00-45	Ø4.35 GUIDED DRILL GUIDE	Ø 4.35
	018-15-F50-00-46	Ø4.85 GUIDED DRILL GUIDE	Ø 4.85
	018-15-F50-00-47	Ø5.35 GUIDED DRILL GUIDE	Ø 5.35
	018-15-F50-00-48	Ø5.85 GUIDED DRILL GUIDE	Ø 5.85
	018-15-F50-00-49	Ø2.8 IMPLANT GUIDE-SLEEVE 5	Ø 2.8
	018-15-F50-00-50	Ø3.0 IMPLANT GUIDE-SLEEVE 5	Ø 3.0
	018-15-F50-00-51	Ø3.5 IMPLANT GUIDE-SLEEVE 5	Ø 3.5
	018-15-F50-00-52	Ø3.0 IMPLANT GUIDE	Ø 3.0
	018-15-F50-00-53	Ø3.5 IMPLANT GUIDE	Ø 3.5
	018-15-F50-00-54	Ø4.0 IMPLANT GUIDE	Ø 4.0
	018-15-F50-00-55	Ø4.5 IMPLANT GUIDE	Ø 4.5
	018-15-F50-00-56	Ø5.0 IMPLANT GUIDE	Ø 5.0
	018-15-F50-00-57	Ø6.0 IMPLANT GUIDE	Ø 6.0
	018-15-F50-00-58	Ø2.8 IMPLANT GUIDE	Ø 2.8
	018-15-F50-00-59	Ø3.75 IMPLANT GUIDE	Ø 3.75
	018-15-F50-00-60	Ø4.1 IMPLANT GUIDE	Ø 4.1
	018-15-F50-00-61	Ø5.5 IMPLANT GUIDE	Ø 5.5

Surgical Accessory - Stop: the stop is a cylindrical surgical accessory recommended for use in assisting professionals in the preparation and drilling of bone tissue with regards to the depth of the drilling carried out by the Plenum drill. The various diameters and heights are directly related to the various models of drills available and the different needs for depth in the drilling of bone tissue.

Stops are available in the codes, models and dimensions listed in table 6 below.

TABLE 6: Illustrative image of stops, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F40-00-13	Ø0.5X3 DRILL STOP	Ø 0.5 x 3.0
	018-15-F40-00-14	Ø0.5X5 DRILL STOP	Ø 0.5 x 5.0
	018-15-F40-00-15	Ø0.5X7 DRILL STOP	Ø 0.5 x 7.0
	018-15-F70-00-01	Ø3.5X5 DRILL STOP	Ø 3.5 x 5.0
	018-15-F70-00-02	Ø3.5X6 DRILL STOP	Ø 3.5 x 6.0
	018-15-F70-00-03	Ø3.5X7 DRILL STOP	Ø 3.5 x 7.0
	018-15-F70-00-04	Ø3.5X8 DRILL STOP	Ø 3.5 x 8.0
	018-15-F70-00-05	Ø3.5X10 DRILL STOP	Ø 3.5 x 10.0
	018-15-F70-00-06	Ø3.5X11 DRILL STOP	Ø 3.5 x 11.5
	018-15-F70-00-07	Ø3.5X13 DRILL STOP	Ø 3.5 x 13.0
	018-15-F70-00-08	Ø3.5X15 DRILL STOP	Ø 3.5 x 15.0
	018-15-F70-00-09	Ø4.5X5 DRILL STOP	Ø 4.5 x 5.0
	018-15-F70-00-10	Ø4.5X6 DRILL STOP	Ø 4.5 x 6.0
	018-15-F70-00-11	Ø4.5X7 DRILL STOP	Ø 4.5 x 7.0
	018-15-F70-00-12	Ø4.5X8 DRILL STOP	Ø 4.5 x 8.0
	018-15-F70-00-13	Ø4.5X10 DRILL STOP	Ø 4.5 x 10.0
	018-15-F70-00-14	Ø4.5X11 DRILL STOP	Ø 4.5 x 11.5
	018-15-F70-00-15	Ø4.5X13 DRILL STOP	Ø 4.5 x 13.0
	018-15-F70-00-16	Ø4.5X15 DRILL STOP	Ø 4.5 x 15.0
	018-15-F70-00-17	Ø6.5X5 DRILL STOP	Ø 6.5 x 5.0
	018-15-F70-00-18	Ø6.5X6 DRILL STOP	Ø 6.5 x 6.0
	018-15-F70-00-19	Ø6.5X7 DRILL STOP	Ø 6.5 x 7.0
	018-15-F70-00-20	Ø6.5X8 DRILL STOP	Ø 6.5 x 8.0
	018-15-F70-00-21	Ø6.5X10 DRILL STOP	Ø 6.5 x 10.0
	018-15-F70-00-22	Ø6.5X11 DRILL STOP	Ø 6.5 x 11.5
	018-15-F70-00-23	Ø6.5X13 DRILL STOP	Ø 6.5 x 13.0
	018-15-F70-00-24	Ø6.5X15 DRILL STOP	Ø 6.5 x 15.0
	018-15-F50-01-76	Ø3.0X5 GUIDED DRILL STOP	Ø 3.0 X 5.0
	018-15-F50-01-77	Ø3.0X6 GUIDED DRILL STOP	Ø 3.0 X 6.0
	018-15-F50-01-78	Ø3.0X7 GUIDED DRILL STOP	Ø 3.0 X 7.0
	018-15-F50-01-79	Ø3.0X8 GUIDED DRILL STOP	Ø 3.0 X 8.0
018-15-F50-01-80	Ø3.0X10 GUIDED DRILL STOP	Ø 3.0 X 10.0	
018-15-F50-01-81	Ø3.0X11.5 GUIDED DRILL STOP	Ø 3.0 X 11.5	
018-15-F50-01-82	Ø3.0X13 GUIDED DRILL STOP	Ø 3.0 X 13.0	
018-15-F50-01-83	Ø3.0X15 GUIDED DRILL STOP	Ø 3.0 X 15.0	
018-15-F50-01-84	Ø5.0X5 GUIDED DRILL STOP	Ø 5.0 X 5.0	
018-15-F50-01-85	Ø5.0X6 GUIDED DRILL STOP	Ø 5.0 X 6.0	
018-15-F50-01-86	Ø5.0X7 GUIDED DRILL STOP	Ø 5.0 X 7.0	
018-15-F50-01-87	Ø5.0X8 GUIDED DRILL STOP	Ø 5.0 X 8.0	
018-15-F50-01-88	Ø5.0X10 GUIDED DRILL STOP	Ø 5.0 X 10.0	
018-15-F50-01-89	Ø5.0X11.5 GUIDED DRILL STOP	Ø 5.0 X 11.5	
018-15-F50-01-90	Ø5.0X13 GUIDED DRILL STOP	Ø 5.0 X 13.0	
018-15-F50-01-91	Ø5.0X15 GUIDED DRILL STOP	Ø 5.0 X 15.0	

Surgical Accessory – Abutment Remover: the Abutment Remover is a cylindrical surgical accessory with a 1.2 mm hexagonal adaptation fitted at one end and a thread at the opposing end. The Abutment Remover is recommended for use in the removal of prosthetic components that have already been fastened to their respective implant (SL, RE and SH) and which need to be removed for a specific reason.

Abutment Removers are available in the codes, models and dimensions listed in table 7 below.

TABLE 7: Illustrative images of Abutment Removers, codes, models and dimensions.

Illustrative Image	Code	Model	Main Dimensions (mm)
	018-15-F60-00-11	ABUTMENT REMOVER SH	Ø 2.5
	018-15-F60-00-12	ABUTMENT REMOVER RE	Ø 2.5
	018-15-F60-00-13	ABUTMENT REMOVER SL	Ø 2.0

2. PRODUCT COMPOSITION

The Surgical Accessory models described in the tables above are manufactured in titanium according to requirements specified by ASTM F67 and ISO 5832-2 for application of the product.

3. PRODUCT PRESENTATION


The Surgical Accessory line is not provided with accessories.

Surgical accessories in the models described in the tables above are provided in single packaging in non-sterile conditions and presented as follows:

- Polyethylene terephthalate (PET) blister packaging thermo-sealed on rigid cardboard and labeling;
- A printed version of the Instructions for Use can be requested at no additional cost (including the cost of shipping) at the telephone number (11) 3109-9000, or through the electronic address (relacionamento@plenum.bio), or through the website (www.plenum.bio).
- Illustrative images of the primary and external packaging of available commercialized Surgical Accessory product models can be found below.

FRONT VIEW OF SURGICAL ACCESSORY PACKAGING.

VIEW OF REVERSE SIDE OF SURGICAL ACCESSORY PACKAGING.

4. PRODUCT RECOMMENDATIONS, PURPOSE OR INTENDED USE

The various models of the Surgical Accessory are recommended for use in dental procedures for the purpose of assisting the dental surgeon in accessing and preparing bone tissue, as well as during the support procedures used in techniques for installation of Plenum Dental Implants and their respective prosthetic components in order to restore masticatory function and aesthetic appearance of patients' oral cavity.

5. PRINCIPLES OF INTENDED PRODUCT FUNCTION AND PERFORMANCE

The main characteristics and principles of operation in Surgical Accessories are related to their mechanism of action, in which instruments are fitted with a non-cutting surface. These mechanisms are used by the dental surgeon to assist in support and complement processes in the installation technique for Dental Implants and their respective prosthetic components.

6. INSTRUCTIONS FOR USE OF THE PRODUCT

Carefully read the instructions for use, recommendations and guidelines prior to using Surgical Accessories;

Verify product packaging before use. Do not use products if the packaging has been tampered with or damaged;

Surgical Accessories are reusables product provided in a non-sterile condition, and must be sterilized before use. For more details, see the section “Sterilization”;

Before starting procedures, the dentist must verify that the set of surgical accessories selected for use are intact and complete. Surgical procedures must not be initiated if the surgical accessory is incomplete or any item is damaged. Any surgical accessories that show signs of wear, oxidation, or physical damage must be identified and disposed of in a suitable location;

The dental surgeon using the surgical accessory is responsible for final judgment with regards to its integrity;

Surgical Accessories are available in a wide range of models, diameters and sizes, and each of these surgical accessories offer very specific functions. The dental surgeon is responsible for choosing the product model with the characteristics and dimensions that are the most appropriate for the procedure to be performed;

Before the surgical procedure, the dental surgeon must implement rigorous preoperative planning before defining a surgical plan in advance, selecting the models and sizes of the surgical accessories that they intend to use;

The surgical technique to be used must be indicated by the dental surgeon using their clinical and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). It is therefore not recommended that a general surgical technique be applied to all patients. Surgical protocols must be carried out according to references and the previous experience of the dental surgeon, taking the most appropriate choice of model and size of surgical accessory, the sequence of use, drilling speed, and use of a cooling system into consideration;

The use, cleaning, disinfection, and sterilization of the instruments must only be carried out by personnel that are trained and qualified for such operations;

7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

- Dental surgeons must not use the Surgical Accessory product line until they have carefully read the instructions for use and the information provided on the packaging. All information on use of instruments and the warnings, precautions and recommendations included in these instructions for use must be taken into consideration. Failure to observe these points may lead to complications for the dental surgeon and especially for the patient;

- Use of these products is restricted to dental applications. Surgical Accessories must only be used by properly trained and qualified dental surgeons with proven experience and knowledge of implant dentistry techniques;

- The surgical accessory’s product labels include the date of manufacture, code, product lot number, basic dimensions, ANVISA (Brazilian National Healthy Regulatory Agency) registration number, technician responsible, company information, and the respective symbology that allows the surgical accessory to be properly identified and traced;

- The dental surgeon must understand and plan all stages involved in the surgical procedure, and verify whether the model and dimension of the

surgical accessory that they intend to use meets requirements for pre-operative planning. Ensure that there are replacements or alternates available for the accessories that you plan to use;

- Before starting the surgery, the dental surgeon must submit the patient to a rigorous oral antiseptis process to prevent the surgical accessories from coming into contact with contaminating substances;

- Incorrect handling of surgical accessories with small dimensions inside the patient’s mouth may present a risk of aspiration and/or swallowing of these elements;

- Improper selection of the model and dimension of the surgical accessory that will be used, as well as errors in the handling and use of the product, can result in damage to the instrument and the patient’s bone tissue, resulting in failure of the procedure;

- The use of inappropriate techniques and procedures for surgical accessories and inadequate surgical conditions may cause harm to patient and compromise the corresponding implant dentistry procedure;

- The primary and most likely cause of failure of a surgical accessory include: improper selection of surgical accessory model and size, lack of proper maintenance, excessive use, use at a rotation or torque above those established, a lack of cooling during use, maintenance performed by third parties, a lack of professional training, technical errors during use, of the use of instruments belonging to other systems. Dental surgeons must take the entirety of these aspects into consideration during preoperative planning;

- Prior to the surgical procedure, accessories must be subject to a technical wash for decontamination purposes, in order to ensure that they are properly cleaned. Cleaning must be carried out at the health service establishment at which the surgery will take place (see item “Recommended Cleaning and Decontamination Methods”);

- Surgical accessories may be damaged and corrode when used, handled, cleaned and sterilized inappropriately;

- It is recommended to fragments of any surgical accessory that may have come apart during surgery be removed. The patient must be informed in cases in which these elements cannot be removed;

- The manufacturer will not be held responsible for reuse of these instruments in cases in which they are used on patients with Creutzfeldt-Jakob disease or an infection resulting from the HIV virus;

- These surgical accessories are considered invasive and come into contact with body tissues and fluids. As a result, there is a risk of biological contamination and transmission of infectious diseases, such as hepatitis and HIV, among others. These accessories must be therefore be treated as potentially contaminating materials;

- In the event that the product code, diameter, or scale of dimensions are no longer legible due to excessive use, the surgical accessory must be immediately replaced;

- Users must never attempt to sharpen a cutting instrument;

- Before each surgery, check to ensure that attachable parts fit properly;

- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the use of these surgical accessories, as well as from the oral implantology procedure;

- It is of the utmost importance that the surgeon explain the entirety of the above topics to their patient, and that patients in turn have a clear and complete understanding of these procedural aspects.

8. CONTRAINDICATIONS

Surgical Accessory products do not present contraindications and side effects, but it is recommended that they be handled by an instrument specialist, dental surgeon or other specialized professional with knowledge in the handling techniques used for this type of product, exclusively for the specific purposes for which they were designed.

9. POTENTIAL ADVERSE EFFECTS

There were no adverse effects related to the use of Surgical Accessories found. However, the following adverse effects may occur due to trauma during surgical procedures or improper use of surgical accessories by qualified professionals, including: injury to nerves, vessels and tissues; hemorrhaging; hematomas; bone necrosis; bone fracture and sensitivity in soft tissues.

10. FINISHING AND SURFACE TREATMENTS



Surgical Accessories are coated with DLC (Diamond-like Carbon) film. DLC is a carbon film with physical and chemical characteristics similar to that of diamonds. The inherent characteristics of DLC include: a high degree of hardness, low coefficients of friction, resistance to chemical corrosion, biocompatibility, bactericidal and anti-coagulant properties, which provides instruments with a high cutting capacity and resistance to wear and corrosion.

11. DESCRIPTION OF SYMBOLS

Table 8 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 8: Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize

	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS

In order to reduce the risk of infection, it is recommended that all Surgical Accessories be cleaned and decontaminated prior to sterilization. Surgical accessories must be removed from packaging and washed with plenty of water at room temperature. Surgical accessories can be disinfected using water and enzymatic detergent and must be packed into their respective cases and drawers to ensure that parts do not come into contact. A broad-spectrum bactericidal and fungicidal solution may also be used for disinfection. In both cases, the following aspects must be observed:

- Always read the instructions from the manufacturer of the disinfectant or detergent product, adhering to the recommended concentrations and immersion period;
- The use of strong mineral acids (hydrofluoric, nitric, hydrochloric, etc.), or products such as zinc chloride, sodium hypochlorite and compositions containing sodium hydroxide is prohibited;
- Do not use abrasive products or instruments (sponges, metal brushes);
- In cases in which thermal baths are used, the water temperature must never exceed 45 ° C, as temperatures above this level can result in the

coagulation of blood and denaturation of proteins, thereby generating solid residues in surgical accessories;

- Manual cleaning, when performed properly, results in less damage to the surgical accessory's surface. In such cases, the following aspects must be observed during cleaning:
- The water temperature must not exceed the ambient temperature;
- The concentrations of cleaning agents used must be appropriate according to the manufacturer's recommendations;
- Use a brush with soft polymeric bristles that is suitable for this type of cleaning. Never use a steel brush or wool sponges, as this may cause the passive layer of the surgical accessory to rupture, causing corrosion;
- Products must be thoroughly rinsed after cleaning. Distilled water must preferably be used since running water may contain high concentration of metal ions and facilitate the corrosion process;
- Dry the product immediately after cleaning using a dry cloth or clean paper towel with fibers that remain intact;
- After cleaning and decontamination, verify the condition of each surgical accessory;
- After cleaning and decontamination, sterilization packages containing the surgical accessory must be stored in a well-ventilated dry area, away from light and the action of the weather. Necessary care must be taken to ensure that the package remains intact before sterilization

13. STERILIZATION

- Surgical Accessory products are provided in a non-sterile condition and must be sterilized by the hospital/health service establishment before use;
- It is recommended that the moist heat sterilization method involving the use of an autoclave be used, adhering to the parameters and procedures established in Plenum's internal protocols and the standard ISO 17665-1. However, the hospital/health service establishment is solely responsible for the methods for cleaning and microbial load control, sterilization equipment, process controls and sterilization procedures used;
- Prior to sterilization, the products that are part of this family must be packaged in an appropriate location or in packaging suitable for sterilization, with the aid of clean tongs and/or gloves to avoid contamination. Sterilization packaging is subject to the method of sterilization, transportation and storage, which has a considerable influence on sterilization results. Packaging must be selected in a manner that ensures that products fit properly inside the packaging. Use a sterilization indicator and write the date and validity of sterilization on the packaging;
- After sterilization, the sterile packaging must be checked for any damage. Sterilization indicators must also be verified. The validity of sterilization must be verified once products have been sterilized. Do not use the products if the sterilization period has expired;
- During the sterilization process, ensure that the maximum load allowed for the autoclave chamber does not exceed the values recommended in the manufacturer's instructions;
- Recommended parameters for moist heat sterilization: 121°C for 20 min, and 30 min drying time. It is important to note, however, that each user must first read and verify the parameters for their equipment (autoclave) and implement the procedure according to the instruction manual for the respective brand and model of autoclave used by designated health professional;

- Validity of sterilization: 7 to 15 days, provided that the product is packaged in a clean, dry, well-ventilated area away from the action of weather;
- **ATTENTION:** Do not perform autoclaving of this product in its original packaging.
- After sterilization, the products must be opened only at the time of use and handled in sterile surgical fields.
- In cases in which products are reused through the moist heat sterilization process, adhere to the respective sterilization procedure. The cleaning procedure used in parts consists of the manual removal of dirt from the surface of the instrument using enzymatic detergent, a brush and/or sponge and water. This step is extremely important in reducing the amount of biological material and allow for the action of the sterilization process to occur. The enzymatic detergent acts on proteins and fats that deposit and adhere to instruments, thereby facilitating cleaning.

14. EXPIRATION DATE

Surgical Accessories are provided in a non-sterile condition. These products have an indefinite validity in the form in which they were originally commercialized.

15. IDENTIFICATION AND TRACEABILITY

- Surgical Accessory products are identified using an engraving that includes the manufacturer's logo and, when applicable, the instrument's dimensions and scale;
- In addition to the identification above, the packaging of Surgical Accessory products includes labeling providing the information necessary in complete identification according to ABNT NBR ISO 15223-1.

16. CONDITIONS OF STORAGE, TRANSPORT, CONSERVATION AND HANDLING

- The manufacturer's packaging must be intact at the moment at which the product is received. Do not use the product if the packaging is damaged or shows signs of tampering. If the packaging is damaged or has been tampered with, the product must be discarded;
- The product must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The product must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the product near lightbulbs to ensure that the information contained in the manufacturer's packaging is not affected;
- Surgical accessories must only be handled by personnel that are qualified and trained to carry out this procedure.

17. DISPOSING OF THE PRODUCT

- Discard the product if the packaging is damaged or open;
- At the end of the surgical accessory's service life, products that are not subject to analysis or subsequent studies must be disposed of by the health service establishment;
- It is recommended that surgical accessories be completely disassembled to ensure that they cannot be reused. The health service establishment is solely responsible for disassembling the product to ensure that it is not reused;
- After disassembly, the instrument must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards and applicable environmental legislation.

18. CUSTOMER SERVICE

- If the product is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, product code and lot, and patient data; directly notifying Plenum's Client Relations Service;
- If necessary, send the product, which has been cleaned, packaged, sterilized and properly identified, and a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Telephone/Fax: (11) 3109- 9000; Email: customerservice@plenum.bio;
- ANVISA also makes a channel available to report incidents, adverse events and technical complaints related to the use of products and services subject to health regulations, which can be found on the website of the National Health Regulatory Agency - ANVISA at the electronic address www.anvisa.gov.br, using the NOTIVISA link.
- Plenum does not perform any type of maintenance and/ or repair on instruments. If the professional identifies any irregularities such as: problems related to quality, or loss of cutting capacity prior to the period indicated (see chapter 6), the product must be sent to the company for analysis using standard customer service procedures. If a failure in the product is confirmed, it will be replaced by another with the same characteristics. It is not recommended that maintenance and sharpening be performed on cutting instruments.
- FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

19. LIMITED WARRANTY

Products manufactured by Plenum have a factory warranty against manufacturing defect. The occurrence of defects in the product must be reported immediately to the manufacturer, observing the legal deadline, attaching a technical report with the details of the non-compliance with

the product. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the section Customer service above) must be observed in cases involving a risk to the patient's health.

The Warranty of products manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. The lack of such compliance and follow-up exempts the manufacturer and distributor from any type of liability, in particular in the following cases:

- Maintenance of the product performed by the user;
- Improper use or use other than in a surgical setting;
- Product not used according to instructions for use.

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