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GUIDED
ACCESSORY
English

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GUIDED ACCESSORY

INSTRUCTIONS FOR USE

TECHNICAL NAME: Dental Surgical Instruments.

TRADE NAME: Guided Accessory.

MODEL: See tables 1 to 4 for a detailed description of the products.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these Instructions for Use prior to using the product. Follow the recommendations, warnings and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications. Product for Medical and Dental Use. Product Supplied without Sterilization. Sterilize before use. Valid for an Indefinite Period of Time. Reusable and Resterilizable Product.

ANVISA REGISTRATION NO.: 81684340015

CODE, BATCH NO. AND DATE OF MANUFACTURE: See product labeling.

Rev. 03

SUMMARY

| | |
|--------------------------------------------------------------|---|
| 1. DETAILED PRODUCT DESCRIPTION | 3 |
| 2. PRODUCT COMPOSITION | 4 |
| 3. PRODUCT PRESENTATION FORM | 4 |
| 4. PRODUCT RECOMMENDATIONS, PURPOSE OR INTENDED USE | 4 |
| 5. PRINCIPLES OF INTENDED PRODUCT FUNCTION AND PERFORMANCE | 4 |
| 6. INSTRUCTIONS FOR USE OF THE PRODUCT | 4 |
| 7. PRECAUTIONS, RESTRICTIONS AND WARNINGS | 5 |
| 8. CONTRAINDICATIONS | 6 |
| 9. POTENTIAL ADVERSE EFFECTS | 6 |
| 10. FINISHING AND SURFACE TREATMENTS | 6 |
| 11. DESCRIPTION OF SYMBOLS | 6 |
| 12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS | 7 |
| 13. STERILIZATION | 7 |
| 14. EXPIRATION DATE | 8 |
| 15. IDENTIFICATION AND TRACEABILITY | 8 |
| 16. STORAGE CONDITIONS, TRANSPORT, CONSERVATION AND HANDLING | 8 |
| 17. DISPOSING OF THE PRODUCT | 8 |
| 18. CLIENT SERVICES | 8 |
| 19. LIMITED WARRANTY | 8 |

1. DETAILED PRODUCT DESCRIPTION

The Guided Accessory is composed of a family of metallic, non-implantable surgical accessories, of non-current profile, developed for application in dental procedures, with the purpose of assisting the dental surgeon in guided surgery procedures, ensuring the stability and fixation of the surgical guide, assisting the support and complement processes in the technique of installing Plenum Dental Implants.

The Guided Accessory is a reusable product, supplied in non-sterile condition, and must be sterilized before use (see the Sterilization topic), and is available in the models listed below.

Guided Accessory - Guide Fixation Pin: The guide fixation pin is a long cylindrical tip guided accessory indicated to assist in fastening the surgical guide in guided surgery procedures in total or partial edentulous patients.

The guide fixation pins are available in the codes, motors and dimensions listed in table 1 below.

TABLE 1: Illustrative image of guide fixation pins, codes, models and dimensions.

| Illustrative Image | Code | Model | Main Dimensions (mm) |
|-------------------------------------------------------------------------------------|------------------|-------------------------|----------------------|
|  | 018-15-F50-02-25 | GUIDE FIXATION PIN 20mm | Ø1.5 x 20 |
| | 018-15-F50-02-26 | GUIDE FIXATION PIN 15mm | Ø1.5 x 15 |
| | 018-15-F50-02-27 | GUIDE FIXATION PIN 25mm | Ø1.5 x 25 |

Guided Accessory - Guide Fixation Screw: The guide fixation screw is a cylindrical-pointed guided accessory that is indicated to assist in fastening the surgical guide in procedures of guided surgery in total or partial edentulous patients.

Guide fixation screws are available in the codes, models and dimensions listed in table 2 below

TABLE 2: Illustrative image of guide fixation screws, codes, models and dimensions

| Illustrative Image | Code | Model | Main Dimensions (mm) |
|-------------------------------------------------------------------------------------|------------------|---------------------------|----------------------|
|  | 018-15-F50-02-28 | GUIDE FIXATION SCREW 20mm | Ø1.5 x 20 |
| | 018-15-F50-02-29 | GUIDE FIXATION SCREW 15mm | Ø1.5 x 15 |
| | 018-15-F50-02-30 | GUIDE FIXATION SCREW 25mm | Ø1.5 x 25 |

Guided Accessory - GUIDE PIN: The guide pin is a guided accessory composed of a guide body and a threaded tip to assist in the fixation of the surgical guide in implants previously installed in guided surgery procedures in total or partial edentulous patients. Guide pins are available in the codes, models and dimensions listed in table 3 below.

TABLE 3: Illustrative image of guide fasteners, codes, models and dimensions.

| Illustrative Image | Code | Model | Main Dimensions (mm) |
|--------------------------------------------------------------------------------------|------------------|------------------|----------------------|
|  | 018-15-F50-02-31 | GUIDE PIN SL-H8 | Ø3.3 x 8 |
| | 018-15-F50-02-32 | GUIDE PIN SL-H9 | Ø3.3 x 9 |
| | 018-15-F50-02-33 | GUIDE PIN SL-H10 | Ø3.3 x 10 |
| | 018-15-F50-02-34 | GUIDE PIN SL-H11 | Ø3.3 x 11 |
| | 018-15-F50-02-35 | GUIDE PIN SL-H12 | Ø3.3 x 12 |
| | 018-15-F50-02-36 | GUIDE PIN RE-H8 | Ø5.3 x 8 |
| | 018-15-F50-02-37 | GUIDE PIN RE-H9 | Ø5.3 x 9 |
| | 018-15-F50-02-38 | GUIDE PIN RE-H10 | Ø5.3 x 10 |
| | 018-15-F50-02-39 | GUIDE PIN RE-H11 | Ø5.3 x 11 |
| | 018-15-F50-02-40 | GUIDE PIN RE-H12 | Ø5.3 x 12 |
| | 018-15-F50-02-41 | GUIDE PIN SH-H8 | Ø6.3 x 8 |
| | 018-15-F50-02-42 | GUIDE PIN SH-H9 | Ø6.3 x 9 |
| | 018-15-F50-02-43 | GUIDE PIN SH-H10 | Ø6.3 x 10 |
| | 018-15-F50-02-44 | GUIDE PIN SH-H11 | Ø6.3 x 11 |
| | 018-15-F50-02-45 | GUIDE PIN SH-H12 | Ø6.3 x 12 |
| | 018-15-F50-03-28 | GUIDE PIN UF-H10 | Ø5.3 x 10 |
| | 018-15-F50-03-29 | GUIDE PIN MF-H10 | Ø5.3 x 10 |

Guided Accessory - Sleeves: The sleeve is a guide attachment that after installed in the surgical guide is responsible for guiding the drilling and installation of the implant according to the implant position provided for in the planning.

The sleeves are available in the codes, models and dimensions listed in table 4 below.

TABLE 4: Illustrative image of sleeves, codes, models and dimensions

| Illustrative Image | Code | Model | Main Dimensions (mm) |
|---------------------------------------------------------------------------------------|------------------|-------------------------------|----------------------|
|  | 018-15-F50-02-46 | GUIDE FIXATION PIN - SLEEVE | Ø1.5 x 8 |
| | 018-15-F50-02-47 | GUIDE FIXATION SCREW - SLEEVE | Ø1.5 x 8 |
|  | 018-15-F50-02-48 | GUIDED SLEEVE Ø3.0 | Ø3.3 x 4 |
| | 018-15-F50-02-49 | GUIDED SLEEVE Ø3.5 | Ø3.8 x 4 |
| | 018-15-F50-02-50 | GUIDED SLEEVE Ø4.0 | Ø4.3 x 4 |
| | 018-15-F50-02-51 | GUIDED SLEEVE Ø4.5 | Ø4.8 x 4 |
| | 018-15-F50-02-52 | GUIDED SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-02-53 | GUIDED SLEEVE Ø5.5 | Ø5.8 x 4 |
| | 018-15-F50-02-54 | GUIDED SLEEVE Ø6.0 | Ø6.3 x 4 |
| | 018-15-F50-03-50 | UF GUIDED SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-03-51 | UF GUIDED SLEEVE Ø6.0 | Ø6.3 x 4 |
| | 018-15-F50-03-52 | MF GUIDED SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-03-53 | MF GUIDED SLEEVE Ø6.0 | Ø6.3 x 4 |

| | | | |
|-----------------------------------------------------------------------------------|------------------|----------------------------------|----------|
|  | 018-15-F50-03-21 | GUIDED FLANGELESS SLEEVE Ø3.0 | Ø3.3 x 4 |
| | 018-15-F50-03-22 | GUIDED FLANGELESS SLEEVE Ø3.5 | Ø3.8 x 4 |
| | 018-15-F50-03-23 | GUIDED FLANGELESS SLEEVE Ø4.0 | Ø4.3 x 4 |
| | 018-15-F50-03-24 | GUIDED FLANGELESS SLEEVE Ø4.5 | Ø4.8 x 4 |
| | 018-15-F50-03-25 | GUIDED FLANGELESS SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-03-26 | GUIDED FLANGELESS SLEEVE Ø5.5 | Ø5.8 x 4 |
| | 018-15-F50-03-27 | GUIDED FLANGELESS SLEEVE Ø6.0 | Ø6.3 x 4 |
| | 018-15-F50-03-54 | UF GUIDED FLANGELESS SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-03-55 | UF GUIDED FLANGELESS SLEEVE Ø6.0 | Ø6.3 x 4 |
| | 018-15-F50-03-56 | MF GUIDED FLANGELESS SLEEVE Ø5.0 | Ø5.3 x 4 |
| | 018-15-F50-03-57 | MF GUIDED FLANGELESS SLEEVE Ø6.0 | Ø6.3 x 4 |

2. PRODUCT COMPOSITION

The Guided Accessories in the models described above in the tables are manufactured in ELI Titanium alloy according to requirements specified by ASTM-F136.

3. PRODUCT PRESENTATION FORM

The Guided Accessory has no accessories.

The Guided Accessories in the models described in tables 1 to 3 are supplied packaged unitarily, and for the commercial models described in table 4 are provided packaged in a set of 5 units in each package. All are provided in non-sterile condition and presented as follows:

- Polyethylene terephthalate (PET) blister pack has been sealed on rigid cardboard and labeling;
- If you want to obtain the Instructions for Use in print format and at no additional cost (including shipping), request it through the telephone (11) 3109-9000, or through the electronic address (relacionamento@plenum.bio), or through the website (www.plenum.bio);
- See below the illustrative images of the packaging of a Guided Accessory model in the form in which it will be commercialized.

Illustrative images of the front of the Guided Accessory packaging.

Illustrative images of the back of the Guided Accessory packaging.

4. PRODUCT RECOMMENDATIONS, PURPOSE OR INTENDED USE

The Guided Accessory in its various models is indicated for application in dental procedures associated with guided surgeries. It aims to assist the dental surgeon in fixing and maintaining the surgical guides in the correct position, previously virtually planned, in guided insertion processes of Plenum Dental Implants and their respective prosthetic components in guided surgery procedures.

5. PRINCIPLES OF INTENDED PRODUCT FUNCTION AND PERFORMANCE

The main feature and basis of operation of the Guided Accessory are related to its mechanism of action, in the Pin, Screw and Guide pin models, where it has at one end a smooth or self-threading tip to perform its specific function of assisting in the fixation of surgical guides, installation of implants or mini-implants of Plenum, in the insertion of Plenum Dental Implants using the technique of guided surgery, thus ensuring stability and maintenance of the surgical guide in its position throughout the surgical procedure.

To ensure a better performance of use of the Guided Accessories, they must be positioned and inserted in an area with a sufficient bone quality that allows the proper installation of these accessories. Additionally, the number of Pins, Screws and Guided Pins must be adapted to the patient's anatomy, type of surgical guide, number and position of implants to be installed.

The commercial models and Sleeves have the guide function for drilling and insertion of the implant at the time of surgery, and the sleeve for fixation pins, have the guide function for pin and/or fixation screw in cases of protocols.

6. INSTRUCTIONS FOR USE OF THE PRODUCT

- Carefully read all instructions for use, recommendations and guidelines, before using the Guided Accessory;
- Verify product packaging before use. Do not use if the packaging is tampered with or damaged;
- The Guided Accessory is a reusable product, supplied in non-sterile condition, and must be sterilized before use. For more details, see the topic "Sterilization";
- Before starting the procedure, the dental surgeon will verify that the Guided Accessory set selected for use is intact and complete. The guided surgical procedure should not be initiated if the Guided Accessory is incomplete or if any items are damaged. Any Guided Accessory that shows signs of wear, oxygenation, or physical damage must be identified and disposed of in a suitable place (see item Product Disposal);
- The judgment as to the integrity of any Guided Accessory is always from the dental surgeon who uses it;
- The Guided Accessory is available in models, diameters and varied dimensions, and each of these guided accessories have very specific functions, and the dentist is responsible for choosing the product model with the most appropriate characteristics and dimensions for the procedure that will be performed according to the virtuous planning previously approved for guided surgery;

- Before the guided surgical procedure, the dental surgeon must carry out a rigorous preoperative planning, through clinical examinations, images (radiographs and computed tomography) and virtual planning software. They should also consider the need for a detailed history and complementary exams on the general health of the patient (complete blood count, coagulation, calcium dosage, etc.), since this information directly influences the body's biological response to the use of the device, surgical guides, and then define a surgical plan in advance, with the selection of the most appropriate implants and accessories for the patient's condition. In view of this, define a surgical plan in advance, with the approval of the planning and selection of the guide model, the models and sizes of the instruments and guided accessory that it intends to use;

- The surgical technique to be used must be indicated by the dental surgeon using their clinical and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). It is therefore not recommended that a general surgical technique be applied to all patients. Surgical protocols must be carried out according to references and the previous experience of the dental surgeon, taking the most appropriate choice of model and size of guided accessory, the sequence of use, drilling speed, and use of a cooling system into consideration;

- In the procedure for using surgical guides, the surgeon must use only appropriate and sterile surgical instruments;

- The surgeon must use previously approved virtual clinical/ surgical planning and specific guided surgery literatures and techniques for using surgical guides for the immediate and late pre- and postoperative periods;

- After starting the surgical procedure, the surgeon may optionally use the guided accessories to ensure the stability and adaptation of the surgical guide in the region of the surgical bed and perform the guided surgical procedure in a safe and stable manner.

Handling and Manipulation of Guided Accessories

- In the procedure for handling and installing guided accessories, commercial models Pin, Screw and Guide Pin, the surgeon must use only appropriate and sterile surgical instruments;

- The handling of the product must be carried out with care, especially during the process of installing the pin, screw and guide pin to the bone or when screwing fasteners into the Plenum implant. The pins, screws, and guide pin must be installed with care so that their body does not buckle and/or cause bone fractures. As for the fixators, avoid forced fitting, as the use of excessive force or insertion in an oblique manner may cause damage to the implant threads;

- The sleeves must be fixed to the surgical guides as determined in the virtual planning, using the same procedure used in the printing of the guides and performing the cleaning process in isopropyl alcohol and UV chamber curing. Use the same post-processing sequence for Surgical Guides.

Installation

- Prior to installation of the guided accessories, the surgeon must rigorously verify that the accessories are intact, complete and pass freely through the surgical guide sleeve.


Figure 1: Installation and use of guide fixation pins and screws.


Figure 2: Installation and use of guide pins

ATTENTION. The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

- The use, cleaning, disinfection, and sterilization of the instruments must only be done by personnel trained and qualified for these operations.

7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

The dental surgeon should not use the Guide Accessory before carefully reading all the instructions for use and the information contained in the packaging. All information on use of instruments and the warnings, precautions and recommendations included in these instructions for use must be taken into consideration. Failure to observe these points may lead to complications for the dental surgeon and especially for the patient;

- Guided Accessories are provided in the non-sterile condition. Sterilization must be carried out in a health facility (hospitals, offices and clinics) by means of a moist heat method, steam under pressure in an autoclave. Do not use products if the packaging has been tampered with or expired;

- Use of these products is restricted to dental applications. The Guided Accessory should only be used by trained, qualified dentists with experience and knowledge of implant dentistry techniques for guided surgery.

- The instruments' product labels include the date of manufacture, code, product lot number, basic dimensions, ANVISA (Brazilian National Healthy Regulatory Agency) registration number, technician responsible, company information, and the respective symbology that allows the guided accessory to be properly identified and traced;

- The dental surgeon must understand and plan all stages of the surgical procedure, and verify that the model and dimension of the guided accessory they intend to use is attentive to the requirements of their preoperative planning. Always have spare or replacement accessories available for what you want to use;

- Before starting the surgery, the dental surgeon must submit the patient to a rigorous oral antisepsis process to prevent the guided accessories from coming into contact with contaminating substances;

- Incorrect handling of surgical accessories with small dimensions inside the patient's mouth may present a risk of aspiration and/or swallowing of these elements;

- Improper selection of the model and dimension of the guided accessory that will be used, as well as errors in the handling and use of the product, can result in damage to the instrument and the patient's bone tissue, resulting in failure of the procedure;

- The use of surgical accessories with techniques, procedures, and unsuitable surgical conditions may impair the patient, causing unsatisfactory results in the implantation procedure;

- The main and most probable causes of failure of a guided accessory are: wrong choice of model and size of the instrument, lack of maintenance,

excessive number of use, torque above the established, maintenance done by third parties, lack of professional training, technical error in use, instrument belonging to other systems. The dentist should consider all these aspects in their preoperative planning;

- Prior to the surgical procedure, accessories must be subject to a technical wash for decontamination purposes, in order to ensure that they are properly cleaned. Cleaning must be carried out at the health service establishment at which the surgery will take place (see item “Recommended Cleaning and Decontamination Methods”);
- Guided accessories can be damaged and corrode, when used, handled, cleaned and sterilized inappropriately;
- In cases of flapless surgery, no force should be applied to the pins/fasteners to avoid damage to the soft tissues;
- Any product that is damaged during handling must be discarded according to the guidelines described in the Product Disposal item;
- It is recommended to remove fragments of any guided accessory that may have fractured during surgery. The patient must be informed in cases in which these elements cannot be removed;
- The manufacturer will not be held responsible for reuse of the instruments in cases in which they are used on patients with Creutzfeldt-Jakob disease or an infection resulting from the HIV virus;
- These guided accessories are considered invasive and come into contact with body tissues and fluids, therefore, there is a risk of biological contamination and transmission of infectious diseases, such as hepatitis, HIV, among others. These accessories must be therefore be treated as potentially contaminating materials;
- In the event that the product code, diameter, or scale of dimensions are no longer legible due to excessive use, the guided accessory must be immediately replaced;
- Before each surgery, check to ensure that attachable parts fit properly;
- The wrong choice or mistakenly validated planning of the surgical guide by the surgeon, as well as errors in the indication, handling, preparation of bone tissue or soft tissues, and installation of the product, can cause damage to the physical structure and contamination of the product and contribute to its failure;
- The surgical guides and their accessories are auxiliary materials for the installation of implants that require a more accessible anatomical reference. Exclusive use of the Plenum Oss hp alone does not guarantee the success of the treatment. The success of therapy depends on the intervention technique and the patient’s systemic reaction;
- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the use of these guided accessories, as well as from the oral implantology procedure;
- It is of the utmost importance that the surgeon explain the entirety of the above topics to their patient, and that patients in turn have a clear and complete understanding of these procedural aspects.

8. CONTRAINDICATIONS

Guided Accessory products do not present contraindications and side effects, but it is recommended that they be handled by an instrument specialist, dental surgeon or other specialized professional with knowledge in the handling techniques used for this type of product, exclusively for the specific purposes for which they were designed.

9. POTENTIAL ADVERSE EFFECTS

There were no adverse effects directly related to the use of Guided Accessory devices found. However, the following adverse effects may occur due to trauma during surgical assignment or improper use of guided accessories by qualified professionals, including: injury to nerves, vessels and tissues; hemorrhaging; hematomas; bone necrosis; bone fracture and sensitivity in soft tissues.

In addition, possible adverse effects associated with the use of the Guided Accessories may include, but are not limited to:

- Allergic reaction to the raw material of the Guided Accessory (see item Product Composition);
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the dental surgeon of all these complications and the risks arising from surgeries with the use of surgical guides.

10. FINISHING AND SURFACE TREATMENTS

The model of Guided Accessory, Guide Pin, is coated by DLC (Diamond-like Carbon) film deposition. DLC is a thin film of carbon with physical and chemical characteristics close to that of diamond. The inefficient characteristics of DLC include: a high degree of hardness, low coefficients of friction, resistance to chemical corrosion, bio-compatibility, bactericidal and anti-coagulant properties, which provides instruments with a high cutting capacity and resistance to wear and corrosion.

11. DESCRIPTION OF SYMBOLS

Table 5 describes the symbols that may be on the product labeling. Please check the physical labeling for the symbols applied to the product.

TABLE 5. Symbols representing standards that are relevant to the product and compliance in use.

| | |
|---------------------------------------------------------------------------------------|----------------------------------------------|
|  | Manufacturer |
|  | Country and date of manufacture |
|  | Use-by date |
|  | Batch code |
|  | Catalog number |
|  | Sterilized using ethylene oxide |
|  | Sterilized using Irradiation |
|  | Sterilized using vaporized hydrogen peroxide |
|  | Non-sterile |
|  | Do not resterilize |

| | |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Do not use if package is damaged |
|  | Single sterile barrier system |
|  | Double sterile barrier system |
|  | Single sterile barrier system with protective packaging inside |
|  | Single sterile barrier system with protective packaging outside |
|  | Keep away from sunlight |
|  | Keep dry |
|  | Lower limit of temperature |
|  | Upper limit of temperature |
|  | Temperature limit |
|  | Do not re-use |
|  | Consult electronic instructions for use (plenum.bio/ifu) |
|  | Unique device identifier |
|  | Device manufacturing material |
|  | Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market) |
| Proibido reprocessar | Reprocessing prohibited (Brazilian Requirement) |

12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS

To reduce the risk of infection, it is recommended that all Guided Accessories be cleaned and decontaminated before sterilization. The guided accessory must be unpacked and washed with plenty of water at room temperature. Disinfection can be carried out with water and enzymatic detergent, with the guided accessories packed in their own cases and drawers, ensuring that there is no friction between the parts. A broad-spectrum bactericidal and fungicidal solution may also be used for disinfection. In both cases, the following aspects must be observed:

- Always read the instructions from the manufacturer of the disinfectant or detergent product, adhering to the recommended concentrations and immersion period;
- The use of strong mineral acids (hydrofluoric, nitric, hydrochloric, etc.), or products such as zinc chloride, sodium hypochlorite and compositions containing sodium hydroxide is prohibited;
- Do not use abrasive products or instruments (sponges metal brushes);
- In case of use of thermal baths, never exceed the water temperature at 45 ° C, as temperatures above this can induce blood coagulation and protein denaturation generating solid residues in the guided accessories;
- Manual cleaning when performed properly causes less damage to the surface of the guided accessory, and in this case, the following topics must be observed:

- The water temperature must not exceed the ambient temperature;
- The concentrations of cleaning agents must be correct, according to the recommendations of the manufacturer;
- Use a brush with soft polymeric bristles that is suitable for this type of cleaning. Never use a steel brush or wool sponges, as this may cause the passive layer of the surgical accessory to rupture, causing corrosion;
- Products must be thoroughly rinsed after cleaning. Distilled water must preferably be used since running water may contain high concentration of metal ions and facilitate the corrosion process;
- Dry the product immediately after cleaning using a dry cloth or clean paper towel with fibers that remain intact;
- After cleaning and decontamination, verify the condition of each guided accessory;
- After cleaning and decontamination, sterilization packages containing the guided accessory must be stored in a well-ventilated dry area, away from light and the action of the weather. Necessary care must be taken to ensure that the package remains intact before sterilization.

13. STERILIZATION

- The Guided Accessory is provided in non-sterile condition, and must be sterilized before use, by the hospital/health service establishment;
 - It is recommended that the moist heat sterilization method involving the use of an autoclave be used, adhering to the parameters and procedures established in Plenum's internal protocols and the standard ISO 17665-1. However, the cleaning methods, microbial load control, methods, sterilization equipment, process controls and sterilization procedures used are the sole responsibility of the hospital/health service establishment;
 - Prior to sterilization, the products that are part of this family must be packaged in an appropriate location or in packaging suitable for sterilization, with the aid of clean tongs and/or gloves to avoid contamination. Sterilization packaging is subject to the method of sterilization, transportation and storage, which has a considerable influence on sterilization results. Packaging must be selected in a manner that ensures that products fit properly inside the packaging. Use a sterilization indicator and write the date and validity of sterilization on the packaging;
 - After sterilization, the sterile packaging must be checked for any damage. Sterilization indicators must also be verified. The validity of sterilization must be verified once products have been sterilized. Do not use the products if the sterilization period has expired;
 - During the sterilization process, ensure that the maximum load allowed for the autoclave chamber does not exceed the values recommended in the manufacturer's instructions;
 - Recommended parameters for moist heat sterilization: 121°C for 20 min, and 30 min drying time. It is important to note, however, that each user must first read and verify the parameters for their equipment (autoclave) and implement the procedure according to the instruction manual for the respective brand and model of autoclave used by designated health professional;
 - Validity of sterilization: 7 to 15 days, provided that the product is packaged in a clean, dry, well-ventilated area away from the action of weather;
- ATTENTION: Do not perform autoclaving of this product in its original packaging.
- After sterilization, the products must be opened only at the time of use and handled in sterile surgical fields.
 - In cases in which products are reused through the moist heat sterilization

process, adhere to the respective sterilization procedure. The cleaning procedure used in parts consists of the manual removal of dirt from the surface of the instrument using enzymatic detergent, a brush and/or sponge and water. This step is extremely important in reducing the amount of biological material and allow for the action of the sterilization process to occur. The enzymatic detergent acts on proteins and fat that deposit and adhere to instruments, thereby facilitating cleaning.

14. EXPIRATION DATE

- Indicated on product labeling.
- Before use, check the expiration date, do not use the product if it has expired.

15. IDENTIFICATION AND TRACEABILITY

- The packaging of the guided accessory is accompanied by labeling with the necessary information for complete identification according to ISO 15223-1.

16. STORAGE CONDITIONS, TRANSPORT, CONSERVATION AND HANDLING

- The manufacturer's packaging must be intact at the moment at which the product is received. Do not use the product if the packaging is damaged or shows signs of tampering. If the packaging is damaged or has been tampered with, the product must be discarded;
- The product must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The product must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the product near lightbulbs to ensure that the information contained in the manufacturer's packaging is not affected;
- The handling of the guided accessory should only be done by qualified and trained personnel to perform this procedure.

17. DISPOSING OF THE PRODUCT

- Discard the product if the packaging is damaged or open;
- After the end of the useful life of the guided accessory, the product that is not subject to further analysis or study must be discarded by the health service establishment;
- It is recommended that guided accessory products be completely disassembled to ensure that they cannot be reused. The health service establishment is solely responsible for disassembling the product to ensure that it is not reused;
- After disassembly, the instrument must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards and applicable environmental legislation.

18. CLIENT SERVICES

- If the product is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, product code and lot, and patient data; directly notifying Plenum's Client Relations Service;
- If necessary, send the product, which has been cleaned, packaged, sterilized and properly identified, and a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Telephone/Fax: (11) 3109- 9000; Email: relacionamento@plenum.bio;
- ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA.
- Plenum does not perform any type of maintenance and/ or repair on instruments. If the professional identifies any irregularities such as: problems related to quality, or loss of cutting capacity prior to the period indicated (see chapter 6), the product must be sent to the company for analysis using standard client service procedures. If a failure in the product is confirmed, it will be replaced by another with the same characteristics. It is not recommended that maintenance and sharpening be performed on cutting instruments.

19. LIMITED WARRANTY

Products manufactured by Plenum have a factory warranty against manufacturing defect. The occurrence of defects in the product must be reported immediately to the manufacturer, observing the legal deadline, attaching a technical report with the details of the non-compliance with the product. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the section Client Services above) must be observed in cases involving a risk to the patient's health.

The Warranty of products manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. The lack of such compliance and follow-up exempt the manufacturer and distributor from any type of liability, in particular in the following cases:

- Maintenance of the product performed by the user;
- Improper use or use other than in a surgical setting;
- If the product has not been used in full compliance with the instructions for use.

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