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ORGANIZER  
BOX  
**English**

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TECHNICAL NAME: Not applicable

TRADE NAME: Organizer Box

MODEL: Detailed description of the device, see table 1

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Website: [www.plenum.bio](http://www.plenum.bio)

#### ATTENTION:

Carefully read these Instructions for Use prior to using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to complications. Device for Dental Use. Device Supplied without Sterilization. Sterilize before use. Valid for an Indefinite Period of Time. Reusable and Resterilizable Device.

ANVISA REGISTRATION NO.: Not applicable

CODE, LOT NO., AND DATE OF MANUFACTURE AND STERILIZATION: See device labeling.

REV.10

## SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	9
3. DEVICE PRESENTATION FORM	9
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	9
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	9
6. INSTRUCTIONS FOR USE	9
7. PRECAUTIONS, RESTRICTIONS AND WARNINGS	10
8. CONTRAINDICATIONS	10
9. POSSIBLE ADVERSE EFFECTS	10
10. SYMBOLS DESCRIPTION	10
11. RECOMMENDED CLEANING AND DECONTAMINATION METHODS	11
12. STERILIZATION	11
13. EXPIRATION DATE	11
14. IDENTIFICATION AND TRACEABILITY	11
15. STORAGE, TRANSPORT, CONSERVATION AND HANDLING	12
16. DISPOSAL OF DEVICE	12
17. CUSTOMER SERVICES	12
18. LIMITED WARRANTY	12

# 1. DETAILED DEVICE DESCRIPTION











The Organizer Box family is composed of several models of metal boxes developed for the packaging of Plenum products so that these products are protected and organized so the dental surgeon can perform the surgical or prosthetic procedure in the best possible way.
















On the inside of the Organizer Box cover, it is possible to find the sequence of surgical protocols when applicable.









IMPORTANT: The Plenum dental products that will be stored in these box models are not part of these instructions for use. They have a separate ANVISA registration, so they must be purchased separately.

The Organizer Boxes are available in the models described in table 1 below.

**TABLE 1:** Illustrative images, codes, models, and dimensions of the Organizer Boxes.



Illustrative Images		Code	Model	Dimensions: Width x Length x Height (mm)
		018-15-F00-00-44	GRAFT ORANGE	135 x 109 x 51
		018-15-F00-00-45	GRAFT GOLD	135 x 109 x 51
		018-15-F00-00-47	GRAFT BLUE	135 x 109 x 51
		018-15-F00-00-48	GRAFT GREEN	135 x 109 x 51
		018-15-F00-00-49	GRAFT WHITE	135 x 109 x 51
		018-15-F00-00-50	GRAFT PINK	135 x 109 x 51
		018-15-F00-00-07	GRAFT GRAY	135 x 109 x 51
		018-15-F00-00-14	SURGICAL COMPACT - ORANGE	135 x 109 x 51
		018-15-F00-00-17	SURGICAL COMPACT - GOLD	135 x 109 x 51
		018-15-F00-00-23	SURGICAL COMPACT - BLUE	135 x 109 x 51
		018-15-F00-00-26	SURGICAL COMPACT - GREEN	135 x 109 x 51

		018-15-F00-00-29	SURGICAL COMPACT - WHITE	135 x 109 x 51
		018-15-F00-00-32	SURGICAL COMPACT - PINK	135 x 109 x 51
		018-15-F00-00-08	SURGICAL COMPACT - GRAY	135 x 109 x 51
		018-15-F00-01-41	HYBRID GRAY	135 x 109 x 51
		018-15-F00-00-13	PROSTHETIC - ORANGE	135 x 109 x 51
		018-15-F00-00-16	PROSTHETIC - GOLD	135 x 109 x 51
		018-15-F00-00-22	PROSTHETIC - BLUE	135 x 109 x 51
		018-15-F00-00-25	PROSTHETIC - GREEN	135 x 109 x 51
		018-15-F00-00-28	PROSTHETIC - WHITE	135 x 109 x 51
		018-15-F00-00-31	PROSTHETIC - PINK	135 x 109 x 51
		018-15-F00-00-04	PROSTHETIC - GRAY	135 x 109 x 51
		018-15-F00-00-12	RE SURGICAL - ORANGE	225 x 143.5 x 53
		018-15-F00-00-15	RE SURGICAL - GOLD	225 x 143.5 x 53
		018-15-F00-00-21	RE SURGICAL - BLUE	225 x 143.5 x 53
		018-15-F00-00-24	RE SURGICAL - GREEN	225 x 143.5 x 53

		018-15-F00-00-27	RE SURGICAL - WHITE	225 x 143.5 x 53
		018-15-F00-00-30	RE SURGICAL - PINK	225 x 143.5 x 53
		018-15-F00-00-03	RE SURGICAL - GRAY	225 x 143.5 x 53
		018-15-F00-01-00	UF SURGICAL - ORANGE	225 x 143.5 x 53
		018-15-F00-01-01	UF SURGICAL - GOLD	225 x 143.5 x 53
		018-15-F00-01-02	UF SURGICAL - BLUE	225 x 143.5 x 53
		018-15-F00-01-03	UF SURGICAL - GREEN	225 x 143.5 x 53
		018-15-F00-01-04	UF SURGICAL - WHITE	225 x 143.5 x 53
		018-15-F00-01-05	UF SURGICAL - PINK	225 x 143.5 x 53
		018-15-F00-01-07	UF SURGICAL - GRAY	225 x 143.5 x 53
		018-15-F00-01-08	MF SURGICAL - ORANGE	225 x 143.5 x 53
		018-15-F00-01-09	MF SURGICAL - GOLD	225 x 143.5 x 53
		018-15-F00-01-10	MF SURGICAL - BLUE	225 x 143.5 x 53
		018-15-F00-01-11	MF SURGICAL - GREEN	225 x 143.5 x 53
		018-15-F00-01-12	MF SURGICAL - WHITE	225 x 143.5 x 53

		018-15-F00-01-13	MF SURGICAL - PINK	225 x 143.5 x 53
		018-15-F00-01-15	MF SURGICAL - GRAY	225 x 143.5 x 53
		018-15-F00-00-52	GUIDED ORANGE UF/MF	225 x 143.5 x 53
		018-15-F00-00-53	GUIDED GOLD UF/MF	225 x 143.5 x 53
		018-15-F00-00-55	GUIDED BLUE UF/MF	225 x 143.5 x 53
		018-15-F00-00-56	GUIDED GREEN UF/MF	225 x 143.5 x 53
		018-15-F00-00-57	GUIDED WHITE UF/MF	225 x 143.5 x 53
		018-15-F00-00-58	GUIDED PINK UF/MF	225 x 143.5 x 53
		018-15-F00-00-10	GUIDED GRAY UF/MF	225 x 143.5 x 53
		018-15-F00-00-95	SURGICAL GUIDED ORANGE	225 x 143.5 x 53
		018-15-F00-00-96	SURGICAL GUIDED GOLD	225 x 143.5 x 53
		018-15-F00-00-97	SURGICAL GUIDED BLUE	225 x 143.5 x 53
		018-15-F00-00-98	SURGICAL GUIDED GREEN	225 x 143.5 x 53
		018-15-F00-00-99	SURGICAL GUIDED WHITE	225 x 143.5 x 53
		018-15-F00-01-24	SURGICAL GUIDED PINK	225 x 143.5 x 53

		018-15-F00-01-25	SURGICAL GUIDED GRAY	225 x 143.5 x 53
		018-15-F00-01-16	GUIDED SURGICAL - ORANGE	225 x 143.5 x 53
		018-15-F00-01-17	GUIDED SURGICAL - GOLD	225 x 143.5 x 53
		018-15-F00-01-18	GUIDED SURGICAL - BLUE	225 x 143.5 x 53
		018-15-F00-01-19	GUIDED SURGICAL - GREEN	225 x 143.5 x 53
		018-15-F00-01-20	GUIDED SURGICAL - WHITE	225 x 143.5 x 53
		018-15-F00-01-21	GUIDED SURGICAL - PINK	225 x 143.5 x 53
		018-15-F00-01-23	GUIDED SURGICAL - GRAY	225 x 143.5 x 53
		018-15-F00-01-34	SURGICAL ORANGE UF/MF	225 x 143.5 x 53
		018-15-F00-01-35	SURGICAL GREEN UF/MF	225 x 143.5 x 53
		018-15-F00-01-36	SURGICAL PINK UF/MF	225 x 143.5 x 53
		018-15-F00-01-37	SURGICAL BLUE UF/MF	225 x 143.5 x 53
		018-15-F00-01-38	SURGICAL GOLD UF/MF	225 x 143.5 x 53
		018-15-F00-01-39	SURGICAL WHITE UF/MF	225 x 143.5 x 53
		018-15-F00-01-40	SURGICAL GRAY UF/MF	225 x 143.5 x 53

		018-15-F00-01-26	SURGICAL LONG DRILLS ORANGE	135 x 109 x 51
		018-15-F00-01-27	SURGICAL LONG DRILLS GOLD	135 x 109 x 51
		018-15-F00-01-28	SURGICAL LONG DRILLS BLUE	135 x 109 x 51
		018-15-F00-01-29	SURGICAL LONG DRILLS GREEN	135 x 109 x 51
		018-15-F00-01-30	SURGICAL LONG DRILLS WHITE	135 x 109 x 51
		018-15-F00-01-31	SURGICAL LONG DRILLS PINK	135 x 109 x 51
		018-15-F00-01-33	SURGICAL LONG DRILLS GRAY	135 x 109 x 51

## 2. DEVICE COMPOSITION

The Organizer Boxes are manufactured with the following raw materials:

- Structure in Aluminum, which gives resistance and durability to the product, fundamental characteristics for the process of storage, transport and sterilization of Plenum instruments;
- Support for fixation of the instruments in Non-toxic Silicon;
- Stainless steel hinges, according to ASTM F899 standard.

## 3. DEVICE PRESENTATION FORM

The Organizer Box does not have accessories.

In the models described in Table 1, the Organizer Boxes are provided in non-sterile condition and packed unitarily in a carton paper box with its respective labeling.

If you want to obtain the Instructions for Use in print format and at no additional cost (including shipping), request it through the telephone +55 (11) 3109-9000, or the electronic address ([customerservice@plenum.bio](mailto:customerservice@plenum.bio)), or the website ([www.plenum.bio/ifu](http://www.plenum.bio/ifu));

See below the illustrative images of the packaging of an Organizer Box model in the form in which it will be marketed.

### ILLUSTRATIVE IMAGES OF THE PACKAGING OF THE ORGANIZER BOX.

## 4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Organizer Box, in its various models, is indicated for the packaging and organization of Plenum auxiliary products in surgical or prosthetic procedures related to the Plenum Dental Implants and Components product line. The Organizer Box keeps the products protected and organized for the processes of storage, transportation, sterilization, and surgical and prosthetic procedures related to implant dentistry so that the dental surgeon can perform the surgical procedure in the best possible way.

## 5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The main feature and basis of operation of the Organizer Box are related to its ability to package the products in an organized and protected manner, ensuring the performance of the sterilization process of the products in all internal points of the box. In addition, the boxes are manufactured with raw materials that guarantee resistance and durability in the processes of reuse and sterilization.

## 6. INSTRUCTIONS FOR USE

- Carefully read all instructions for use, recommendations, and guidelines before using an Organizer Box model;
- Verify product packaging before use. Do not use products if the packaging has been tampered with or damaged;
- The Organizer Box is reusable and supplied in non-sterile condition and must be sterilized before use. For more details, see the section “Sterilization”;
- The Organizer Box is identified by means of laser engraving containing the manufacturer’s logo and models of the products that will be packaged therein; in addition, laser engraved surgical protocol when applicable;
- The products that will be packed in the Organizer Box are not part of these instructions for use, have separate ANVISA registration, and must be purchased separately. These products must be allocated in an organized manner, in the exact places indicated in the recordings contained in the box, so that the dental surgeon has an adequate sequence of use at the time of the surgery or prosthetic procedure;
- Before starting the procedure, the dental surgeon will verify that the set of products selected for use and allocated inside the box is intact and complete. Surgery should not be initiated if the product set is incomplete or if any product or box is damaged. Any item that shows signs of wear, oxidation, or physical damage must be identified and disposed of in a suitable place;
- The final judgment as to the integrity of any product is always from the dental surgeon who uses it;
- The Organizer Box is available in models with varied sets and dimensions, and each of these boxes has characteristics and compositions of very specific product sets for each type of surgical or prosthetic procedure and the dentist is responsible for choosing the box model with the most appropriate characteristics and compositions for the procedure that will be performed;
- Before the surgical procedure, the dental surgeon must carry out rigorous preoperative planning and previously define a surgical plan, with the selection of the set, models, and sizes of the products he/she intends to use;

- The surgical technique to be used must be indicated by the surgeon using their clinic and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). It is, therefore, not recommended that a general surgical technique be applied to all patients. The surgical protocol must be carried out according to references and previous experiences of the dental surgeon, always considering the most appropriate choice of the set, model, and size of the products and the sequence of use;
- The use, cleaning, disinfection, and sterilization of the box must only be done by trained and qualified personnel for these operations.

## 7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The dental surgeon should not use an Organizer Box model before carefully reading all instructions for use and the information contained in the packaging. All information for use, warnings, precautions, and recommendations mentioned in these instructions for use should be considered. Failure to observe these points can lead to complications for the surgeon and especially for the patient;
- Use of these products is restricted to dental applications. The Organizer Box and the products that are packaged therein must only be used by trained, qualified dental surgeons with experience and knowledge of implant dentistry techniques;
- The box label shows the date of manufacture, code, batch number of the product, ANVISA registration number, technical manager, company data, and the respective symbol, allowing the full identification and traceability of the product;
- The dental surgeon must understand and plan all stages of the surgical procedure and verify that the model and the composition of the box they intend to use meet the requirements of their preoperative planning.
- Improper selection and choice of the model and composition of the box that will be used, as well as errors in the handling and use of the product, can contribute to unsatisfactory results in the implant dentistry procedure;
- The main and most probable causes of failure of an organizer box and set of surgical instruments are wrong choice of model and composition of the box, lack of maintenance, excessive number of uses, maintenance done by third parties, lack of professional training, technical error in use, products belonging to other systems. The dental surgeon should consider all these aspects in the preoperative planning;
- Before the surgical procedure, the box must undergo a technical wash for decontamination in order to ensure the correct cleaning of it. Cleaning must be carried out at the health service establishment at which the surgery will take place (see item “Recommended Cleaning and Decontamination Methods”);
- In the case of reuse, the box can lose its original functional characteristics through continuous use, handling, and constant cleaning and sterilization processes. For the same reasons, hinges, cover latches, and instrument holders may deteriorate, loosen, or break. As a result, every box must be carefully checked before use. Do not use a product that is damaged;
- In the event of illegibility of the recorded information (dimensions, description of the instruments, suggestion of protocol, etc.) due to excessive use, the box must be immediately replaced;
- Maintenance, rewriting, or replacement of parts (hinges, structure, and supports) in an organizer box on behalf of the user is prohibited;
- Before each surgery, check the perfect fit and order of the products packed in the box because products out of order or misplaced can lead the dental

surgeon to errors;

- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the use of these products, as well as from the oral implantology procedure;
- It is of the utmost importance that the surgeon explains the entirety of the above topics to their patient and that patients, in turn, have a clear and complete understanding of these procedural aspects.

## 8. CONTRAINDICATIONS

The Organizer Box models do not present any indications and side effects, but it is recommended that the handling be carried out by an instrument specialist, dental surgeon or other specialized professional and knowledgeable in the techniques of handling this type of product. Only use an organizer box for the specific purposes for which it was designed.

## 9. POSSIBLE ADVERSE EFFECTS

No adverse effects related directly to the use of the Organizer Box were found. However, adverse effects may occur due to trauma of the surgical procedure or improper use of this product, such as a sequence of misuse of instructions due to the poor packaging of the products in the box.

## 10. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be printed on product labeling. Please check the physical labeling for the symbols applied to the product. The symbols express information related to the identification and use of the product.

**TABLE 2** - Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Non-sterile
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not re-use

 plenum.bio/ifu	Consult electronic instructions for use
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist.  (Notification required by FDA for United States market)

**NOTE:** The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

## 11. RECOMMENDED CLEANING AND DECONTAMINATION METHODS

To reduce the risk of infection, it is recommended that every Organizer Box be cleaned and decontaminated before sterilization. The box must be removed from the packaging and washed with plenty of water at room temperature. Disinfection can be carried out with water and enzymatic detergent, ensuring that there is no friction between parts. A broad-spectrum bactericidal and fungicidal solution may also be used for disinfection. In both cases, the following topics should be noted:

- Always read the instructions from the manufacturer of the disinfectant or detergent product, adhering to the recommended concentrations and immersion period
- The use of strong mineral acids (hydrofluoric, nitric, hydrochloric, etc.) or products such as zinc chloride, sodium hypochlorite, and compositions containing sodium hydroxide is prohibited;
- Do not use abrasive products or instruments (sponges, metal brushes);
- In case of use of thermal baths, never exceed the water temperature at 45°C, as temperatures above this can induce blood coagulation and protein denaturation, generating solid residues in the box;

Manual cleaning, when performed properly, causes less damage to the box surface, in which case the following topics should be noted:

- The water temperature must not exceed room temperature;
- The concentrations of cleaning agents used must be appropriate according to the manufacturer's recommendations;
- Use a brush with soft polymeric bristles that is suitable for this type of cleaning. Never use a steel brush or wool sponges, as this may cause the passive layer of the box to rupture, causing corrosion;
- The box should be thoroughly rinsed after cleaning. Distilled water must preferably be used since running water may contain high concentrations of metal ions and facilitate the corrosion process;
- Dry the box immediately after cleaning with a dry cloth or clean paper, which does not loosen fibers;
- After cleaning and decontamination, check the condition and condition of the product;
- After cleaning and decontamination, sterilization packages containing the box must be stored in a well-ventilated dry area, away from light and the action of the weather. Necessary care must be taken to ensure that the package remains intact before sterilization.

## 12. STERILIZATION

- The Organizer Box is provided in non-sterile condition and must be sterilized before use by the hospital/health service establishment;
- It is recommended that the moist heat sterilization method involving the use of an autoclave be used, adhering to the parameters and procedures established in Plenum's internal protocols and the standard ISO 17665-1. However, the hospital/health service establishment is solely responsible for the methods for cleaning and microbial load control, sterilization equipment and methods, process controls, and sterilization procedures used;
- Before sterilization, the organizer box containing Plenum products must be packed in an appropriate place or in packaging suitable for sterilization, with the aid of clean gloves to avoid contamination. The sterilization package is subject to the sterilization, transport and storage method, which has considerable influence on the sterilization result. The packaging must be selected so that the product fits properly inside it. Use a sterilization indicator and write the date and validity of sterilization on the packaging;
- After sterilization, the sterile packaging must be checked for any damage. Sterilization indicators must also be verified. The validity of sterilization must be verified once products have been sterilized. Do not use the product if the sterilization period has expired;
- Recommended parameters for moist heat sterilization: 121°C for 20 min and 30 min drying time. It is important to note, however, that each user must first read and verify the parameters for their equipment (autoclave) and implement the procedure according to the instruction manual for the respective brand and model of autoclave used by the designated health professional;
- **ATTENTION:** Do not perform autoclaving of this product in its original packaging.
- After sterilization, the product must be opened only at the time of use and handled in sterile surgical fields;
- In cases in which boxes are reused through the moist heat sterilization process, adhere to the respective sterilization procedure. The product cleaning procedure consists of manual removal of dirt through physical action applied to the surface of the box using enzymatic detergent, brush and/or sponge and water. This step is extremely important to reduce the amount of biological material and allow the action of the sterilization process. The enzymatic detergent acts on proteins and fats that deposit and adhere to the box, thereby facilitating cleaning.

## 13. EXPIRATION DATE

The Organizer Box is provided in non-sterile condition. This product has an indefinite validity in the form in which it was originally marketed.

## 14. IDENTIFICATION AND TRACEABILITY

The packaging of the box is accompanied by labeling with the necessary information for complete identification according to ABNT NBR ISO 15223-1.

## 15. STORAGE, TRANSPORT, CONSERVATION AND HANDLING

- The manufacturer's packaging must be intact at the time of receipt, do not use the box if the packaging is damaged or tampered with. If the packaging is damaged or has been tampered with, the product must be discarded;
- The box must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The box must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the product near lightbulbs to ensure that the information contained in the manufacturer's packaging is not affected;
- The handling of the product should only be done by qualified and trained personnel to perform this procedure.

## 16. DISPOSAL OF DEVICE

- Discard the product if the packaging is damaged or open;
- At the end of the box's service life, products that are not subject to analysis or subsequent studies must be disposed of by the health service establishment;
- It is recommended that the box be completely disassembled to ensure that they cannot be reused. The health service institution/establishment is solely responsible for disassembling the product to ensure that it is not reused;
- After de-features, the box must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards in accordance with applicable current environmental legislation;

## 17. CUSTOMER SERVICES

- If the product is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, product code and lot, and patient data; directly notifying Plenum's Client Relations Service;
- If necessary, send the clean, decontaminated, sterilized, packaged and properly identified product, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Telephone: (11) 3109-9000; Email: [customerservice@plenum.bio](mailto:customerservice@plenum.bio)
- ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at [www.anvisa.gov.br](http://www.anvisa.gov.br), link NOTIVISA.
- Plenum does not perform any type of maintenance and/or repair on instruments and boxes. If the professional finds any irregularity and quality problems, it must be forwarded to the company for analysis, following the normal customer service procedures. If the failure of the product is confirmed, it will be replaced by another with the same characteristics. The maintenance and laser rewriting of the content described in the respective Organizer Box models is not recommended.

## 18. LIMITED WARRANTY

The products manufactured by Plenum have a factory warranty against manufacturing defects. The occurrence of defects in the product must be reported immediately to the manufacturer, observing the legal deadline and attaching a technical report with the details of the non-compliance with the product. This report must be generated and signed by a licensed and qualified health professional, also observing in case of risk to the patient's health, the standards belonging to ANVISA (see Customer Service).

The Warranty of products manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and follow-up exempts the manufacturer and distributor of the products from any kind of liability, in particular in the following cases:

- Maintenance of the product performed by the user;
- Improper use or use other than in a surgical setting;
- Product not used according to instructions for use.

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