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DRIVER
English

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DRIVER

INSTRUCTIONS FOR USE

TECHNICAL NAME: Instruments for Dental Use

TRADE NAME: Driver

MODEL: Detailed description of the device see table 1

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these Instructions for Use prior to using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to complications. Device for Dental Use. Device Supplied without Sterilization. Sterilize before use. Valid for an Indefinite Period of Time. Reusable and Resterilizable Device.

ANVISA REGISTRATION NO.: 81684340002

CODE, LOT NO., AND DATE OF MANUFACTURE AND STERILIZATION: See device labeling.

Rev. 06

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	4
3. DEVICE PRESENTATION FORM	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	5
6. INSTRUCTIONS FOR USE	5
7. PRECAUTIONS, RESTRICTIONS AND WARNINGS	5
8. CONTRAINDICATIONS	6
9. POSSIBLE ADVERSE EFFECTS	6
10. SURFACE FINISH	6
11. SYMBOLS DESCRIPTION	6
12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS	7
13. STERILIZATION	7
14. EXPIRATION DATE	8
15. IDENTIFICATION AND TRACEABILITY	8
16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING	8
17. DISPOSAL OF DEVICE	8
18. CUSTOMER SERVICES	8
19. LIMITED WARRANTY	8

1. DETAILED DEVICE DESCRIPTION

The Driver family consists of several models, with hexagonal, slit or square profile connection options, with activation by the dental surgeon through the use of a contra-angle or manually. These drivers were developed for application in dental procedures in order to assist in the processes of insertion and installation of dental implants of Plenum and its respective prosthetic components.

The Driver in its various models is a reusable product, supplied in non-sterile condition, and must be sterilized before use (see topic Sterilization).

The Driver is available in the models described in table 1 below.

TABLE 1: Illustrative images, codes, models, and dimensions of the Driver

Illustrative image	Code	Model	Main dimensions (mm)
	018-15-F10-00-65	IMPLANT DRIVER SL-12mm	1.6
	018-15-F10-00-66	IMPLANT DRIVER SL-17mm	1.6
	018-15-F10-00-67	IMPLANT DRIVER SL-10mm	1.6
	018-15-F10-00-68	IMPLANT DRIVER RE-12mm	2.0
	018-15-F10-00-69	IMPLANT DRIVER RE-17mm	2.0
	018-15-F10-00-70	IMPLANT DRIVER SH-10mm	2.0
	018-15-F10-00-71	IMPLANT DRIVER SH-11mm	2.0
	018-15-F10-00-72	IMPLANT DRIVER SH-16mm	2.0
	018-15-F10-00-73	IMPLANT DRIVER SH-9mm	2.0
	018-15-F10-02-70	IMPLANT DRIVER UF-10mm	2.0
	018-15-F10-02-71	IMPLANT DRIVER UF-12mm	2.0
	018-15-F10-02-72	IMPLANT DRIVER UF-17mm	2.0
	018-15-F10-02-73	IMPLANT DRIVER MF-9mm	2.0
	018-15-F10-02-74	IMPLANT DRIVER MF-11mm	2.0
	018-15-F10-02-75	IMPLANT DRIVER MF-16mm	2.0
	018-15-F50-00-25	GUIDED IMPLANT DRIVER SL-8	1.6
	018-15-F50-00-26	GUIDED IMPLANT DRIVER RE-8	2.0
	018-15-F50-00-27	GUIDED IMPLANT DRIVER SH-8	2.0
	018-15-F10-00-75	HEXALOBULAR DRIVER 1.2-31mm	1.2
	018-15-F20-00-01	HEXALOBULAR DRIVER 1.2-24mm	1.2
	018-15-F20-00-02	HEXALOBULAR DRIVER 1.2-36mm	1.2
	018-15-F10-00-76	HEXAGONAL DRIVER 1.2-31mm	1.2
	018-15-F20-00-03	HEXAGONAL DRIVER 1.2-24mm	1.2
	018-15-F20-00-04	HEXAGONAL DRIVER 1.2-36mm	1.2
	018-15-F10-00-77	HEXAGONAL DRIVER 2.0-17mm	2.0
	018-15-F20-00-05	HEXAGONAL DRIVER 2.0-7mm	2.0
	018-15-F20-00-06	HEXAGONAL DRIVER 2.0-23mm	2.0
	018-15-F40-00-10	PLENUM GRAFT DRIVER - 5mm	2.0
	018-15-F40-00-11	PLENUM GRAFT DRIVER - 10mm	2.0
	018-15-F40-00-12	PLENUM GRAFT DRIVER - 20mm	2.0
	018-15-F40-00-18	PLENUM GRAFT DRIVER - 50mm	2.0

	018-15-F10-00-78	IMPLANT RATCHET DRIVER SL 12mm	1.6
	018-15-F10-00-79	IMPLANT RATCHET DRIVER SL 17mm	1.6
	018-15-F10-00-80	IMPLANT RATCHET DRIVER SL 10mm	1.6
	018-15-F10-00-81	IMPLANT RATCHET DRIVER RE 12mm	2.0
	018-15-F10-00-82	IMPLANT RATCHET DRIVER RE 17mm	2.0
	018-15-F10-00-83	IMPLANT RATCHET DRIVER RE 10mm	2.0
	018-15-F10-00-84	IMPLANT RATCHET DRIVER SH 11mm	2.0
	018-15-F10-00-85	IMPLANT RATCHET DRIVER SH 16mm	2.0
	018-15-F10-00-86	IMPLANT RATCHET DRIVER SH 9mm	2.0
	018-15-F10-02-76	IMPLANT RATCHET DRIVER UF 10mm	2.0
	018-15-F10-02-77	IMPLANT RATCHET DRIVER UF 12mm	2.0
	018-15-F10-02-78	IMPLANT RATCHET DRIVER UF 17mm	2.0
	018-15-F10-02-79	IMPLANT RATCHET DRIVER MF 9mm	2.0
	018-15-F10-02-80	IMPLANT RATCHET DRIVER MF 11mm	2.0
	018-15-F10-02-81	IMPLANT RATCHET DRIVER MF 16mm	2.0
	018-15-F50-00-28	GUIDED IMP MANUAL DRIVER SL-8	1.6
018-15-F50-00-29	GUIDED IMP MANUAL DRIVER RE-8	2.0	
018-15-F50-00-30	GUIDED IMP MANUAL DRIVER SH-8	2.0	
	018-15-F10-00-87	MANUAL HEXAGONAL 1.2-14mm	1.2
	018-15-F10-00-88	MANUAL HEXAGONAL 1.2-9mm	1.2
	018-15-F10-00-89	MANUAL HEXAGONAL 1.2-19mm	1.2
	018-15-F20-00-07	HEXAGONAL 0.9-16mm	0.9
	018-15-F20-00-28	HEXAGONAL 0.9-6mm	0.9
	018-15-F20-00-08	HEXAGONAL 0.9-21mm	0.9
	018-15-F20-00-09	HEXAGONAL 1.2-16mm	1.2
	018-15-F20-00-10	HEXAGONAL 1.2-6mm	1.2
	018-15-F20-00-11	HEXAGONAL 1.2-21mm	1.2
	018-15-F20-00-12	HEXALOBULAR 1.2-16mm	1.2
	018-15-F20-00-13	HEXALOBULAR 1.2-6mm	1.2
	018-15-F20-00-14	HEXALOBULAR 1.2-21mm	1.2
	018-15-F20-00-29	HEXAGONAL 1.6-16mm	1.6
	018-15-F20-00-15	HEXAGONAL 1.6-6mm	1.6
	018-15-F20-00-16	HEXAGONAL 1.6-21mm	1.6
	018-15-F20-00-23	HEXAGONAL 2.0-17mm	2.0
	018-15-F20-00-24	HEXAGONAL 2.0-7mm	2.0
	018-15-F20-00-25	HEXAGONAL 2.0-23mm	2.0
	018-15-F20-00-17	ABUTMENT SQUARE DRIVER-16mm	1.3
	018-15-F20-00-18	ABUTMENT SQUARE DRIVER-6mm	1.3
	018-15-F20-00-19	ABUTMENT SQUARE DRIVER-21mm	1.3
	018-15-F20-00-20	ABUTMENT SLIT DRIVER-16mm	-
	018-15-F20-00-21	ABUTMENT SLIT DRIVER-6mm	-
	018-15-F20-00-22	ABUTMENT SLIT DRIVER-21mm	-
	018-15-F10-00-98	MANUAL IMPLANT DRIVER	4.0
	018-15-F40-00-17	PLENUM GRAFT MANUAL DRIVER	2.0
	018-15-F50-01-92	GUIDED IMPLANT DRIVER SL-9	1.7
	018-15-F50-01-93	GUIDED IMPLANT DRIVER SL-10	1.7
	018-15-F50-01-94	GUIDED IMPLANT DRIVER SL-11	1.7
	018-15-F50-01-95	GUIDED IMPLANT DRIVER SL-12	1.7

	018-15-F50-01-96	GUIDED IMPLANT DRIVER SL	1.7
	018-15-F50-01-97	GUIDED IMPLANT DRIVER RE-9	2.0
	018-15-F50-01-98	GUIDED IMPLANT DRIVER RE-10	2.0
	018-15-F50-01-99	GUIDED IMPLANT DRIVER RE-11	2.0
	018-15-F50-02-00	GUIDED IMPLANT DRIVER RE-12	2.0
	018-15-F50-02-01	GUIDED IMPLANT DRIVER RE	2.0
	018-15-F50-02-02	GUIDED IMPLANT DRIVER SH-9	2.5
	018-15-F50-02-03	GUIDED IMPLANT DRIVER SH-10	2.5
	018-15-F50-02-04	GUIDED IMPLANT DRIVER SH-11	2.5
	018-15-F50-02-05	GUIDED IMPLANT DRIVER SH-12	2.5
	018-15-F50-02-06	GUIDED IMPLANT DRIVER SH	2.5
	018-15-F50-03-30	GUIDED IMPLANT DRIVER RI5 UF	2.0
	018-15-F50-03-31	GUIDED IMPLANT DRIVER RI5 MF	2.5
	018-15-F50-03-34	GUIDED IMPLANT DRIVER RI5 UF-10	2.0
	018-15-F50-03-35	GUIDED IMPLANT DRIVER RI5 MF-10	2.5
	018-15-F50-03-38	GUIDED IMPLANT DRIVER RI6 UF	2.0
	018-15-F50-03-39	GUIDED IMPLANT DRIVER RI6 UF-10	2.0
	018-15-F50-03-40	GUIDED IMPLANT DRIVER RI6 MF	2.5
	018-15-F50-03-41	GUIDED IMPLANT DRIVER RI6 MF-10	2.5
	018-15-F50-03-42	GUIDED IMPLANT DRIVER RI5 SX	2.5
	018-15-F50-03-43	GUIDED IMPLANT DRIVER RI5 SX-10	2.5
	018-15-F50-02-07	GUIDED IMP MANUAL DRIVER SL-9	1.7
	018-15-F50-02-08	GUIDED IMP MANUAL DRIVER SL-10	1.7
	018-15-F50-02-09	GUIDED IMP MANUAL DRIVER SL-11	1.7
	018-15-F50-02-10	GUIDED IMP MANUAL DRIVER SL-12	1.7
	018-15-F50-02-11	GUIDED IMP MANUAL DRIVER SL	1.7
	018-15-F50-02-12	GUIDED IMP MANUAL DRIVER RE-9	2.0
	018-15-F50-02-13	GUIDED IMP MANUAL DRIVER RE-10	2.0
	018-15-F50-02-14	GUIDED IMP MANUAL DRIVER RE-11	2.0
	018-15-F50-02-15	GUIDED IMP MANUAL DRIVER RE-12	2.0
	018-15-F50-02-16	GUIDED IMP MANUAL DRIVER RE	2.0
	018-15-F50-02-17	GUIDED IMP MANUAL DRIVER SH-9	2.5
	018-15-F50-02-18	GUIDED IMP MANUAL DRIVER SH-10	2.5
	018-15-F50-02-19	GUIDED IMP MANUAL DRIVER SH-11	2.5
	018-15-F50-02-20	GUIDED IMP MANUAL DRIVER SH-12	2.5
	018-15-F50-02-21	GUIDED IMP MANUAL DRIVER SH	2.5
	018-15-F50-03-32	GUIDED IMP MANUAL DRIVER RI5 UF	2.0
	018-15-F50-03-33	GUIDED IMP MANUAL DRIVER RI5 MF	2.5
	018-15-F50-03-36	GUIDED IMP MANUAL DRIVER RI5 UF-10	2.0
	018-15-F50-03-37	GUIDED IMP MANUAL DRIVER RI5 MF-10	2.5
	018-15-F50-03-44	GUIDED IMP MANUAL DRIVER RI6 UF	2.0
	018-15-F50-03-45	GUIDED IMP MANUAL DRIVER RI6 UF-10	2.0
	018-15-F50-03-46	GUIDED IMP MANUAL DRIVER RI6 MF	2.5
	018-15-F50-03-47	GUIDED IMP MANUAL DRIVER RI6 MF-10	2.5
	018-15-F50-03-48	GUIDED IMP MANUAL DRIVER RI5 SX	2.5
	018-15-F50-03-49	GUIDED IMP MANUAL DRIVER RI5 SX-10	2.5
	018-15-F10-01-26	MANUAL HEXALOBULAR-14mm	1.4
	018-15-F10-01-27	MANUAL HEXALOBULAR-9mm	1.4
	018-15-F10-01-28	MANUAL HEXALOBULAR-19mm	1.4
	018-15-F50-02-22	GUIDED FIX SCREW-MED.	1.4
	018-15-F50-02-23	GUIDED FIX SCREW-SHORT	1.4
	018-15-F50-02-24	GUIDED FIX SCREW-LONG	1.4
	018-15-F20-00-44	RESCUE HEXAG DRIVER 1.2-31mm	1.2
	018-15-F20-00-45	RESCUE HEXAG DRIVER 1.2-24mm	1.2
	018-15-F20-00-46	RESCUE HEXAG DRIVER 1.2-36mm	1.2
	018-15-F50-04-81	GUIDED IMPLANT DRIVER RI3.5 MF	2.5
	018-15-F50-04-79	GUIDED IMPLANT DRIVER RI3.5 RE	2.0
	018-15-F50-04-80	GUIDED IMPLANT DRIVER RI3.5 UF	2.0
	018-15-F50-04-85	GUIDED IMP MANUAL DRIVER RI3.5 MF	2.5

	018-15-F50-04-83	GUIDED IMP MANUAL DRIVER RI3.5 RE	2.0
	018-15-F50-04-84	GUIDED IMP MANUAL DRIVER RI3.5 UF	2.0
	018-15-F50-04-82	GUIDED IMP MANUAL DRIVER RI3 SL	1.7
	018-15-F50-04-90	GUIDED IMP MANUAL DRIVER RI6 SH	2.5
	018-15-F50-04-86	GUIDED IMP MANUAL DRIVER RI5 RE	2.0
	018-15-F50-04-87	GUIDED IMP MANUAL DRIVER RI5 UF	2.0
	018-15-F50-04-88	GUIDED IMP MANUAL DRIVER RI5 MF	2.5
	018-15-F50-04-89	GUIDED IMP MANUAL DRIVER RI5 SX	2.5
	018-15-F50-04-90	GUIDED IMP MANUAL DRIVER RI5 SH	2.5

2. DEVICE COMPOSITION

The Drivers in the models described above in the tables are made of stainless steel, whose raw material provides the products with high strength and durability, according to the requirements specified by ASTM F899, according to the application of the product.

3. DEVICE PRESENTATION FORM

The Driver has no accessories.

The Driver in the models described in table 1 is provided in unitary packaging, in non-sterile condition, and presented as follows:

- Thermo-sealed polyethylene terephthalate (PET) blister packaging on rigid cardboard and labeling;
- A printed version of the Instructions for Use including shipping), can be requested at no additional cost at the telephone number (11) 3109-9000, or through the electronic address (customerservice@plenum.bio), or through the website (plenum.bio/ifu);
- See below the illustrative images of the packaging of a model of the Driver in the form in which it will be marketed.

Front illustrative image of Driver packaging

Back illustrative image of Driver packaging

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Driver in its various models is indicated for application in dental procedures, with activation through the use of a contra-angle or manually, with the purpose of assisting the dental surgeon in the processes of insertion and installation of Plenum Dental Implants and their respective prosthetic components, with the final purpose of restoring the masticatory function and aesthetics of the patient.


5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The main characteristic and basis of operation of the Driver are related to its mechanism of action, where it has at one end a hexagonal, slit, or square profile connection, which is triggered through the use of a contra-angle or manually by the dental surgeon, to perform its specific function of assisting him in the processes of insertion and installation of Plenum Dental Implants and their respective prosthetic components.

6. INSTRUCTIONS FOR USE

Carefully read all instructions for use, recommendations, and guidelines, before using the Driver;

- Verify product packaging before use. Do not use products if the packaging has been tampered with or damaged;
- The Driver is provided in non-sterile condition, and must be sterilized before use. For more details, see the section “Sterilization”;
- Before starting the procedure, the dentist will check that the drivers selected for use are clean and functional. Special attention must be paid to the hexagonal, slit or square connection of the driver, as wear or damage in this region of the instrument may cause damage to the connection of the implant and/or components, proving the integrity and proper fixation of the products;
- Any driver that shows signs of wear, oxidation, or physical damage must be identified and discarded in an appropriate location;
- The final judgment as to the integrity of any instrument is always from the dental surgeon who uses it;
- The Driver is available in models and varied dimensions, and each of these drivers have has very specific functions., The and the dentist is responsible for choosing the product model with the most appropriate characteristics and dimensions for the procedure that will be performed;
- The sequential and correct use of the drivers reduces the possibility of injuries to the products and allows the correct technique of insertion and installation of the implant and components. An inappropriate use may compromise the primary fixation of the implant next to the bone tissue and its subsequent osseointegration;
- Before the surgical procedure, the dental surgeon must see a rigorous pre-operative planning carried out, and then previously define a surgical plan, with the selection of models and sizes of drivers that they want to use;
- The surgical technique to be used must be indicated by the dental surgeon using their clinical and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). It is, therefore, not recommended that a general surgical technique be applied to all patients. The surgical protocol must be performed according to the references and previous experiences of the dental surgeon, always considering the most appropriate choice of the model and size of the drivers, and the sequence of use;
- The use, cleaning, disinfection, and sterilization of the driver must only be done by well-trained and capable personnel for these operations.

ATTENTION. Figures are merely illustrative. They do not represent the actual dimensions and characteristics of the product.

7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The dental surgeon should not use the Driver before carefully reading all the instructions for use and the information contained in the packaging. All information on use of the ratchet and the warnings, precautions, and recommendations included in these instructions for use must be taken into consideration. Failure to observe these points may lead to complications for the dental surgeon and especially for the patient;
- Use of these products is restricted to dental applications. The Driver must only be used by properly trained and qualified dental surgeons with proven experience and knowledge of implant dentistry techniques;
- The instruments' product labels include the date of manufacture, code, product lot number, basic dimensions, ANVISA (Brazilian National Healthy Regulatory Agency) registration number, technician responsible, company information, and the respective symbology that allows the driver to be properly identified and traced;
- The dental surgeon must understand and plan all stages of the surgical procedure, and always verify that the model and dimension of the driver they intend to use meets the requirements of his preoperative planning. Always have spare or replacement drivers available that you want to use;
- Before starting the surgery, the dental surgeon must submit the patient to a rigorous oral antisepsis to prevent the driver from coming into contact with contaminating substances;
- Incorrect handling of the driver with small dimensions inside the patient's mouth may present a risk of aspiration and/or swallowing of these elements;
- Improper selection and choice of the model and size of the driver to be used, as well as errors in the handling and use of the product, can cause damage to the instruments, implant and/or component, contributing to the failure of the procedure;
- The use of drivers with inadequate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results in the implantology procedure;
- The main and most likely causes of failure of a surgical instrument,

including the drivers, are: choosing the wrong model and size of the instrument, lack of maintenance, excessive number of times it is used, maintenance performed by third parties, lack of professional training, technical error in its use, use with instruments belonging to other systems. The dental surgeon should consider all of these aspects in their preoperative planning;

- Before the surgical procedure, the driver must go through a technical wash for decontamination, in order to ensure the correct cleaning. Cleaning must be carried out at the health service establishment at which the surgery will take place (see item “Recommended Cleaning and Decontamination Methods”);
- The driver can be damaged and suffer corrosion, when used, handled, cleaned and sterilized improperly;
- It is recommended to fragments of any instrument that may have come apart during surgery be removed. The patient must be informed in cases in which these elements cannot be removed;
- In the event of possible use of this instrument in patients with Creutzfeldt-Jakob disease or an infection caused by the HIV or hepatitis virus, the manufacturer disclaims any responsibility for reuse;
- These instruments are considered invasive and come into contact with body tissues and fluids. As a result, there is a risk of biological contamination and transmission of infectious diseases, such as hepatitis and HIV, among others. Therefore, this instrument must be treated as potentially contaminating material;
- In case of illegibility of the code, size or scale of dimensions due to excessive use, the driver must be immediately replaced;
- Before each surgery, check to ensure that attachable parts fit properly;
- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the use of these surgical instruments, as well as from the oral implantology procedure;
- It is of the utmost importance that the surgeon explain the entirety of the above topics to their patient, and that patients in turn have a clear and complete understanding of these procedural aspects.

8. CONTRAINDICATIONS

The Driver does not present contraindications and side effects, but it is recommended that the handling of this instrument be performed by an instrument specialist, dental surgeon or other specialized professional and knowledgeable in the techniques of handling this type of product. Only use Plenum Drivers for the specific purposes for which they were designed.

9. POSSIBLE ADVERSE EFFECTS

No adverse effects related directly to the use of the Driver were found. However, the following adverse effects may occur due to the trauma of the surgical procedure or inadequate use of instruments: nerve, vessel, and tissue injuries; hemorrhage; hematomas; superficial or deep infection; thrombosis; pulmonary embolism; bone necrosis; bone fracture and soft tissue sensitivity.






10. SURFACE FINISH

The Driver is coated by DLC film deposition (Diamond-like Carbon). DLC is a thin film of carbon with physical and chemical characteristics close to that of a diamond. The inherent characteristics of DLC include: a high degree of hardness, low coefficients of friction, resistance to chemical corrosion, biocompatibility, bactericidal and anti-coagulant properties, which provides provide instruments with a high cutting capacity and resistance to wear and corrosion.

11. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 2: Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature

	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
Proibido reprocessor	Reprocessing prohibited (Brazilian Requirement)

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

12. RECOMMENDED CLEANING AND DECONTAMINATION METHODS

To reduce the risk of infection, it is recommended that the entire Driver be cleaned and decontaminated prior to sterilization. The product must be unpacked and washed with plenty of water at room temperature. Disinfection can be carried out with water and enzymatic detergent, with the driver packed in its own cases and drawers, ensuring that there is no friction between the parts. A broad-spectrum bactericidal and fungicidal solution may also be used for disinfection. In both cases, the following aspects must be observed:

- Always read the instructions from the manufacturer of the disinfectant or detergent product, adhering to the recommended concentrations and immersion period;
- The use of strong mineral acids (hydrofluoric, nitric, hydrochloric, etc.), or products such as zinc chloride, sodium hypochlorite, and compositions containing sodium hydroxide is prohibited;
- Do not use abrasive products or instruments (sponges, metal brushes);
- In cases in which thermal baths are used, the water temperature must never exceed 45 ° C, as temperatures above this level can result in the coagulation of blood and denaturation of proteins, thereby generating solid residues in instruments;
- Manual cleaning, when performed properly, results in less damage to the instrument's surface. In such cases, the following aspects must be observed during cleaning:
 - The water temperature must not exceed room temperature;
 - The concentrations of cleaning agents used must be appropriate according to the manufacturer's recommendations;
 - Use a brush with soft polymeric bristles that is suitable for this type of cleaning. Never use a steel brush or wool sponges, as this may cause the passive layer of the instrument to rupture, causing corrosion;
 - The driver should be thoroughly rinsed after cleaning. Distilled water must preferably be used since running water may contain high concentrations of metal ions and facilitate the corrosion process;
 - Dry the product immediately after cleaning using a dry cloth or clean paper towel with fibers that remain intact;

- After cleaning and decontamination, verify the condition of each driver;
- After cleaning and decontamination, sterilization packages containing the product must be stored in an a well-ventilated dry area, away from light and the action of the weather. Necessary care must be taken to ensure that the package remains intact before sterilization.
- The driver should be thoroughly rinsed after cleaning. Distilled water must preferably be used since running water may contain high concentrations of metal ions and facilitate the corrosion process;
- Dry the product immediately after cleaning using a dry cloth or clean paper towel with fibers that remain intact;
- After cleaning and decontamination, verify the condition of each driver;
- After cleaning and decontamination, sterilization packages containing the product must be stored in an well-ventilated dry area, away from light and the action of the weather. Necessary care must be taken to ensure that the package remains intact before sterilization.

13. STERILIZATION

- The Driver is provided in non-sterile condition, and must be sterilized before use by the hospital/health service establishment;
- It is recommended that the moist heat sterilization method involving the use of an autoclave be used, adhering to the parameters and procedures established in Plenum's internal protocol and the standard ISO 17665-1. However, the methods of cleaning and controlling the microbial load, methods and equipment of sterilization procedures, process controls and sterilization procedures used are the sole responsibility of the hospital/health service establishment;
- Before sterilization, the driver must be packed in an appropriate place or in packaging suitable for sterilization, with the aid of clean tweezers and/or gloves to avoid contamination. Sterilization packaging is subject to the method of sterilization, transportation and storage, which has a considerable influence on sterilization results. Packaging must be selected in a manner that ensures that products fit properly inside the packaging. Use a sterilization indicator and write the date and validity of sterilization on the packaging;
- After sterilization, the sterile packaging must be checked for any damage. Sterilization indicators must also be verified. After the driver is sterile, check the validity of the sterilization. Do not use the product if the sterilization period has expired;
- During the sterilization process, ensure that the maximum load allowed for the autoclave chamber does not exceed the values recommended by the manufacturer's instructions;
- Recommended parameters for moist heat sterilization: 121°C for 20 min and 30 min drying time. It is important to note, however, that each user must first read and verify the parameters for their equipment (autoclave) and implement the procedure according to the instruction manual for the respective brand and model of autoclave used by the designated health professional;
- Validity of sterilization: 7 to 15 days, provided that the product is packaged in a clean, dry, well-ventilated area away from the action of weather;
ATTENTION: Do not perform autoclaving of this product in its original packaging.
- After sterilization, the driver must be opened only at the time of use and handled in sterile surgical fields;
- In cases in which products are reused through the moist heat sterilization process, adhere to the respective sterilization procedure. The cleaning procedure

used in parts consists of the manual removal of dirt from the surface of the driver using enzymatic detergent, a brush and/or sponge, and water. This step is extremely important in reducing the amount of biological material and allow allowing for the action of the sterilization process to occur. The enzymatic detergent acts on proteins and fats that deposit and adhere to the product, thereby facilitating cleaning.

14. EXPIRATION DATE

The Driver is provided in non-sterile condition. This product has an indefinite validity in the form in which it was originally marketed.

15. IDENTIFICATION AND TRACEABILITY

- The Driver is identified by means of laser engraving, having the manufacturer's type, and when applicable, the size and dimensional scale;
- In addition to the identification above, the packaging of Instrumental products includes labeling providing the information necessary in for complete identification according to NBR ISO 15223-1.

16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING

- The manufacturer's packaging must be intact at the moment at which the product is received. Do not use the product if the packaging is damaged or shows signs of tampering. If the packaging is damaged or has been tampered with, the product must be discarded;
- The product must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight, and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The product must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the product near lightbulbs to ensure that the information contained in the manufacturer's packaging is not affected;
- The handling of the product should only be done by qualified and trained personnel to perform this procedure.

17. DISPOSAL OF DEVICE

- Discard the product if the packaging is damaged or open;
- At the end of the ratchet's service life, devices that are not subject to analysis or subsequent studies must be disposed of by the health service establishment;
- It is recommended that Ratchets be completely disassembled to ensure that they cannot be reused. The health service establishment is solely responsibly responsible for disassembling the device to ensure that it is not reused;
- After disassembly, the instrument must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards and applicable environmental legislation.

18. CUSTOMER SERVICES

If the product is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, product code and lot, and patient data; directly notifying Plenum's Customer Services;

- If necessary, send the product, which has been cleaned, packaged, sterilized, and properly identified, and a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Aina Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Telephone/Fax: (11) 3109- 9000; E-mail: customerservice@plenum.bio
- ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA.
- Plenum does not perform any type of maintenance and/or repair on the drivers. If the professional finds any irregularity and or quality problems in the product, the same must be sent to the company for analysis, following the normal customer service procedures. If a failure in the product is confirmed, it will be replaced by another with the same characteristics. The maintenance and laser rewriting of the substance described in the respective Drivers models is are not recommended.

19. LIMITED WARRANTY

Products manufactured by Plenum have a factory warranty against manufacturing defect. Product defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-conformity identified in the product. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the section Customer Services above) must be observed in cases involving a risk to the patient's health.

The Warranty of products manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. The lack of such compliance and follow-up exempt the manufacturer and distributor from any type of liability, in particular in the following cases:

- Maintenance of the product performed by the user;
- Improper use or use other than in a surgical setting;
- Product not used according to instructions for use.

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