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DIGITAL  
WORKFLOW  
**English**

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### INSTRUCTIONS FOR USE

---

TECHNICAL NAME: Dental Implant Components

TRADE NAME: Digital Workflow.

MODEL: For a detailed description of the product, see Table 1.

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E-mail: [customerservice@plenum.bio](mailto:customerservice@plenum.bio)

Website: [www.plenum.bio](http://www.plenum.bio)

Carefully read these Instructions for Use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile Device: Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only, do not resterilize or reuse.

ANVISA REGISTRATION NO.: 81684340006.

DEVICE CODE, BATCH CODE, AND MANUFACTURING DATE: See device labeling.

## SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	3
3. DEVICE PRESENTATION FORM	3
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLES OF AND INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	4
7. PRECAUTIONS	4
8. WARNINGS	5
9. CONTRAINDICATIONS	5
10. POTENTIAL ADVERSE EFFECTS	5
11. SYMBOLS DESCRIPTION	5
12. AUXILIARY INSTRUMENTS	6
13. STERILIZATION	6
14. EXPIRATION DATE	6
15. IDENTIFICATION AND TRACEABILITY	6
16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING CONDITIONS	7
17. DISPOSAL OF DEVICE	7
18. CUSTOMER SERVICES	7
19. LIMITED WARRANTY	7

## 1. DETAILED DEVICE DESCRIPTION

The Digital Workflow is a dental device, polymeric, available in the Scanbody model, as listed in Table 1.

The Digital Workflow- Scanbody Model attaches to the implant and/or intermediary component for digital reading during the scanning procedure. Scanning the Scanbody allows you to detect the spatial position of the implants later to transfer the information to a virtual planning software.

Table 1 below shows an illustrative image, codes, models, and main dimensions of the Scanbody.

**TABLE 1:** Illustrative images of Scanbody - codes, models, and dimensions

Illustrative image	Code	Model	Dimensions (mm)
	018-21-C06-20-00	Implant scanbody SL	Ø3x12.7
	018-22-C06-20-00	implant scanbody RE	Ø3x12.6
	018-23-C06-20-00	implant scanbody SH	Ø4x12
	018-20-C06-05-33	Ti-base scanbody Ø3.9	Ø3.9x6.4
	018-20-C06-05-45	Ti-base scanbody Ø4.5	Ø4.5x6.4
	018-20-C06-20-00	Multi abutment scanbody	Ø4.8x10
	018-20-C06-50-33	Plenum abt. Scanbody Ø3.5x4	Ø3.5x10
	018-20-C06-51-33	Plenum abt. Scanbody Ø3.5x6	Ø3.5x10
	018-20-C06-50-45	Plenum abt. Scanbody Ø4.5x4	Ø4.5x10
	018-20-C06-51-45	Plenum abt. Scanbody Ø4.5x6	Ø4.5x10
	018-20-C06-20-36	Multi abutment scanbody Ø3.6	Ø3.6x10
	018-20-C06-90-35	Univ. Abt. Scanbody Ø3.5 H4	Ø3.5x10
	018-20-C06-91-35	Univ. Abt. Scanbody Ø3.5 H6	Ø3.5x10
	018-20-C06-90-45	Univ. Abt. Scanbody Ø4.5 H4	Ø4.5x10
	018-20-C06-91-45	Univ. Abt. Scanbody Ø4.5 H6	Ø4.5x10

## 2. DEVICE COMPOSITION

The Digital Workflow in the models described in Table 1 is manufactured in PEEK polymer (polyether ether ketone), according to the specifications of ASTM F2026 and ABNT NBR 15723-8.

The fixation screw is made of titanium alloy (Ti-6Al- 4V), according to ASTM F136.

## 3. DEVICE PRESENTATION FORM

The Digital Workflow has no accessories.

The Digital Workflow is supplied packaged unitarily, in the sterile condition, and presented as follows:

- Primary packaging is a polyethylene terephthalate (PET) blister and heat-sealed with surgical grade paper;
- External packaging is a rigid cardboard box containing the primary packaging;
- A printed version of the Instructions for Use (including shipping) can be requested at no additional cost at the telephone number +55 11 3109-9000 or through the electronic address (customerservice@plenum.bio), or the website (www.plenum.bio/ifu);
- See below the illustrative images of the primary and external packaging of a Digital Workflow model in the form in which it is commercialized.

**FIGURE 1: PRIMARY PACKAGING OF THE DIGITAL WORKFLOW + FIXATION SCREW (EXCEPT TI-BASE SCANBODIES).**

**FIGURE 2: EXTERNAL PACKAGING OF THE DIGITAL WORKFLOW.**

## 4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Digital Workflow – Scanbody models are dental devices that aim to transfer the positioning of implants or intermediate components through a digital reading performed through a scanning procedure, which can be performed using an intraoral or benchtop scanner.

Digital Components are for transient use; where they are used in implant and component positioning transfer procedures that last less than sixty (60) minutes.

## 5. PRINCIPLES OF AND INTENDED DEVICE FUNCTION AND PERFORMANCE

The Digital Workflow is an auxiliary element in the virtual planning of a dental prosthesis, which will be installed on the Plenum Dental Implant.

With the use of Digital Workflow, it is possible to transfer the positioning of the implant or intermediate components using an intraoral or bench scanner.

Digital Workflow can be used either directly in the patient's mouth or the plaster model, commonly used in prosthesis laboratories.

In virtual planning, it is possible to select the models and dimensions of the prosthetic components to be used later, allowing greater predictability in the prosthetic results.

## 6. INSTRUCTIONS FOR USE

- Carefully read all instructions for use, recommendations, and guidelines, before using the Digital Workflow;
- Before use, check the device packaging. Do not use the component if the packaging is tampered with or damaged;
- Device supplied in sterile condition using Gamma Irradiation;
- The Digital Workflow family comprises Scanbody models, whose primary function is to be used as auxiliary elements in the virtual planning of a dental prosthesis on implantation (see specific functions in the topic "Device Description"), being up to the dental surgeon the choice of device models with the most appropriate characteristics and dimensions, after the individual clinical analysis of each clinical case;
- The techniques for the use and application of the component models that are part of this family vary according to the preference of the dental and prosthetic surgeon. The dental surgeon is responsible for selecting the best clinical approach;
- The Digital Workflow is always used in conjunction with the Plenum Dental Implant and other Plenum prosthetic components, which are provided separately;
- Carefully observe the correct combination and fit between the Digital Workflow and the respective implant or intermediate component;
- For the proper installation of the prosthetic component, it is essential to use Plenum instruments, as recommended by Plenum (see information on the website: [www.plenum.bio](http://www.plenum.bio));

- Before starting the procedure, make sure that the necessary components and instruments are intact and complete;
- The technique, scanning equipment, and software to be used must be chosen by the dental surgeon, who must consider, among other things, their previous clinical experience and the techniques established by the scientific literature;
- Special care must be taken in the hexagonal connection existing in the fixation screws of some models of Digital Workflow and in the driver in its introduction, because an unsuitable fit between the two elements can cause damage to the connection of the device, compromising its integrity and proper fixation;
- Care must also be taken with the existing contact surface between the components. Both surfaces must be free of particles and solid materials that may prevent the perfect fitting and fixation of the device;
- To protect the patient and the surgeon, it is important to ensure the correct and complete traceability of the component used through the annotation in the patient's protocol of the device code and batch.

### HOW TO USE THE DIGITAL WORKFLOW

- 1) The dentist or prosthetic surgeon must select the Digital Workflow – Scanbody Model according to the installed implant or intermediate component if used;
- 2) After selecting the Digital Workflow, installation must be carried out with the adaptation to the implant or intermediate component installed in the mouth or plaster model; For Scanbody models installed directly on Implants or Mini Abutments, the fixation must be carried out with the use of a 1.2 mm hexagonal driver and manual torque (not exceeding 10 Ncm). Scanbody models installed on Ti-Bases must be adapted and fitted on them. The verification of the fitting must be carried out rigorously by the dentist or prosthetic surgeon;  
  
NOTE: It is recommended that the veneering surface be positioned towards the vestibular surface and not adjacent to the proximal teeth.
- 3) Once the Digital Workflow is installed in the mouth or plaster model, virtual scanning is performed using an intraoral scanner or bench scanner to capture the positioning of the implant or intermediate component;  
  
NOTE: The Digital Workflow is indicated for single use, and reuse is not recommended.
- 4) The files obtained from the scan performed are validated and used in software selected for the planning of the prostheses to be made. The planning can be carried out either by the clinician or by the laboratory. The Digital Workflows can be selected in the Virtual Library installed in the planning software, at this stage the Ti-Bases and Multi Abutments can be chosen according to the needs of the clinical case;
- 5) After the planning, a digital file of the prosthesis to be made will be obtained for sending to manufacturing;
- 6) The prostheses will be obtained by machining, using specific milling machines according to the selected material;
- 7) Resin models can be made using 3D printers so that the finishing of the prostheses can be performed before installation on the patient;

## 7. PRECAUTIONS

- The dentist and/or prosthetic surgeon must not use the Digital Workflow before carefully reading all the instructions for use and the

information contained in the packaging. All usage information, warnings, precautions, and recommendations mentioned in these instructions for use should be considered. A failure to observe these points may result in complications during procedures;

- This device is restricted to dental use and should only be used by trained dental surgeons with experience and knowledge of implant dentistry techniques and dental prostheses on implants;
- Date of manufacture, device code, lot number, and other information regarding identification and traceability are provided on the device labeling;
- The opening of the packaging for the use of the device must only be done by qualified and trained personnel to perform this procedure.
- Before the procedure, the dental surgeon must carry out rigorous clinical planning, selecting the model and size of the prosthetic component most appropriate to the patient's condition;
- Surgeons must be aware of the patient's clinical conditions. In the preoperative examination, the site where the component will be used must be rigorously verified by means of clinical examinations and imaging. The dental surgeon must consider the general clinical status of the patient and analyze local or systemic factors that may interfere with the success of the procedure;
- Patients with inadequate oral hygiene, clinical conditions with hematological disorders, and other uncontrolled systemic changes present restrictions on the use of this device;
- For safety and effectiveness and not to compromise the intended use of components, only prosthetic components and implants manufactured by Plenum should be used. Components and/ or implants manufactured by other companies can cause dimensional incompatibility, improper fitting, micro-movements, increase the risk of corrosion (electrolytic or galvanic) and fracture of these devices, compromising the result of the surgical and prosthetic procedure;
- Before use, special attention should be paid to the insertion drivers. Damage to the tip of the drivers can prevent the correct installation of the component, compromising its intended performance. Therefore, this device must not be used if there is any damaged instrument for specific use;
- Improper selection of the component that will be used, as well as errors in the indicated use, handling, and installation of the device, may result in damage to the structure of intermediary components, which may contribute to failure due to fatigue, fracturing, or loosening of the component;
- Improper selection of the component that will be used, as well as errors in the indicated use, handling, and installation of the device, may interfere with virtual planning and may result in a poorly adapted prosthesis;
- Do not force introducing a prosthetic component beyond the established depth. This action may cause deformation in the hexagonal connection of the component and/or the introduction driver and fracturing of the device;
- Incorrect handling of components of small sizes inside the patient's mouth may incur a risk of aspiration and/or swallowing of these elements;
- The success of a procedure for placing a dental implant and prosthetic component is related to the correct indication and selection of devices, proper surgical technique, as well as the time and quality of osseointegration, which can be affected by numerous biological, clinical, and biomechanical factors, which are extrinsic to the quality of the devices. These factors can cause failure and early loss of devices and compromise the surgical procedure. Therefore, faithful obedience to all indications for use, contraindications, precautions, warnings, and other observations contained in these instructions for use are essential for the completion of the procedure;

- Elderly patients, patients with mental health problems, individuals with chemical dependency, or any other difficulty in adhering to postoperative recommendations may represent a higher risk of failure of this device since these patients may ignore instructions for this procedure;

- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the use of prosthetic components, as well as from the oral implantology procedure;

- Inform the need for the patient to maintain follow-up with periodic returns to the dental surgeon so that the proper position and integration of the implant and components can be evaluated;

- It is fundamental that the dental surgeon expresses all the above topics to his patient so that the patient has a complete and clear understanding of all these aspects.

## 8. WARNINGS

- Single-use only, do not resterilize or reuse.
- Do not use if package is opened or damaged or if expiration date has been exceeded.

## 9. CONTRAINDICATIONS

The Digital Workflow is always used with the Plenum Dental Implant and other Plenum prosthetic components. Therefore, as it is a device for use with transient contact to the tissues surrounding the implant and prosthetic components, the Digital Workflow does not present contraindications and side effects to the patient. However, it is recommended that Laboratory Components be used by trained professionals that have specialized in the area and are knowledgeable in the techniques of handling this type of component. These components must only be used for the specific purpose they were designed.

- Contraindication only for patients with allergy or sensitivity to the raw material used in the manufacture of the Digital Workflow (see "Device Composition").

## 10. POTENTIAL ADVERSE EFFECTS

As it is a device for transient contact with tissues surrounding the implant and prosthetic components, it does not present adverse effects to the patient directly related to the use of Digital Workflow. However, wrong choice of model or component size or the use of inappropriate techniques may make the preparation of a suitable dental prosthesis unfeasible

## 11. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

**TABLE 2** – Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not re-sterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist.  (Notification required by FDA for United States market)
Proibido reprocessar	Reprocessing prohibited (Brazilian Requirement)

**NOTE:** The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

## 12. AUXILIARY INSTRUMENTS

- All surgical instruments developed exclusively for the installation of the Digital Workflows are considered auxiliary instruments;
- The list of surgical instruments can be found in the respective printed or electronic communication materials or on the Plenum website (www.plenum.bio). The use of Plenum instruments must follow the manufacturer's instructions;
- Surgical instruments are not the main purpose of these instructions for use, that is, they have separate registrations with ANVISA, therefore, they must be purchased separately. Use only instruments manufactured by Plenum, so there is no incongruity between the fittings and irreversible damage to the components and instruments.

## 13. STERILIZATION

- The Digital Workflow is supplied, in a sterile condition, sterilized using Gamma Irradiation;
- Single-use device, do not reprocess, re-sterilize, or reuse;
- Before use, check the integrity of the packaging and the expiration date. Do not use the device if the package is violated or with the expiration date expired;
- Date of manufacture, sterilization expiration date, code, batch number, and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments;
- The packaging opening for the device's use must only be done by qualified and trained personnel to perform this procedure.

## 14. EXPIRATION DATE

- Indicated on device labeling.
- Before use, check the expiration date, do not use the device if it has expired.

## 15. IDENTIFICATION AND TRACEABILITY

The packaging of the Digital Workflow is accompanied by labeling with the necessary information for complete identification according to NBR ISO 15223-1.

## 16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING CONDITIONS

- The device packaging must be intact at the time of receipt, do not use the device if the packaging is damaged. If the packaging is damaged, the device must be discarded;
- The data for identification is on the device labeling;
- The device must be stored and transported in a clean environment, protected from heat or humidity sources, away from the incidence of sunlight, and free from the action of weathering, to prevent any damage or change in the packaging or its physical form that makes it impossible to use;
- The storage of the device on the shelves or drawers/cabinets must be in such a way that facilitates the visualization and manipulation of the operator. Do not store the device near lamps so that the information contained in the packaging is not affected;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

## 17. DISPOSAL OF DEVICE

- After use, the prosthetic component that is not the subject of further analysis or study must be discarded by the hospital/health service establishment;
- It is recommended that the device be completely disassembled to ensure that it cannot be reused. The de-featuring of the device is the sole responsibility of the health service institution /establishment;
- After de-featuring, the device must be identified as being unfit for use and disposed of accordingly to the waste management standards of health services (hospitals, offices, and clinics) in accordance with applicable current environmental legislation;
- The device that is unused during surgery or laboratory processing must also be discarded following the guidelines described above.

## 18. CUSTOMER SERVICES

- If the device is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, device code, batch code and patient data, directly notifying Plenum's Customer Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and properly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Phone: +55 11 3109-9000; Email: [customerservice@plenum.bio](mailto:customerservice@plenum.bio);
- ANVISA also provides a channel to receive notifications of incidents, adverse events, and technical complaints related to the use of devices and services under health surveillance, which can be found on the website of the

National Health Surveillance Agency - ANVISA, at [www.anvisa.gov.br](http://www.anvisa.gov.br), link NOTIVISA.

- FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

## 19. LIMITED WARRANTY

Devices manufactured by Plenum have a factory warranty against manufacturing defects. The occurrence of defects in the device must be reported immediately to the manufacturer, observing the legal deadline and attaching a technical report with the details of the non-compliance with the device. This report must be generated and signed by a licensed and qualified health professional, also observing, in case of risk to the patient's health, the standards belonging to ANVISA (see Customer Service).

The Warranty of the devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications of use, contraindications, precautions, warnings and other information contained in the instructions for use. The lack of such compliance and follow-up exempts the manufacturer and distributor of the devices from any kind of liability, in particular in the following cases:

- Reprocessing, resterilization or reuse of the device by the user;
- Improper use or use other than in a surgical setting;
- Device not used according to instructions for use.

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