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HEALING  
ABUTMENT  
**English**

COMPONENTE  
PEEK  
**Español**

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### INSTRUCTIONS FOR USE

---

TECHNICAL NAME: Dental Implant Components

TRADE NAME: Healing Abutment

MODEL: For a detailed description of the device, see Table 1.

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E-mail: [customerservice@plenum.bio](mailto:customerservice@plenum.bio)

Website: [www.plenum.bio](http://www.plenum.bio)

ATTENTION: Carefully read these Instructions for Use prior to using the device. Follow the recommendations, warnings and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile Device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited to Reprocess, Resterilize, and Reuse.

ANVISA DEVICE REGISTRATION NO.: 81684340009.

DEVICE CODE, BATCH CODE, AND DATE OF MANUFACTURE AND STERILIZATION: See device labeling.

## SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	5
3. DEVICE PRESENTATION FORM	5
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	5
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	6
6. INSTRUCTIONS FOR USE	6
7. PRECAUTIONS	6
8. WARNING	7
9. CONTRAINDICATIONS	7
10. POSSIBLE ADVERSE EFFECTS	8
11. SYMBOLS DESCRIPTION	8
12. AUXILIARY INSTRUMENTS	9
13. STERILIZATION	9
14. EXPIRATION DATE	9
15. IDENTIFICATION AND TRACEABILITY	9
16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING	9
17. DISPOSAL OF DEVICE	10
18. CUSTOMER SERVICES	10
19. LIMITED WARRANTY	10

# 1. DETAILED DEVICE DESCRIPTION

The PEEK Components family includes prosthetic abutments that are manufactured in two parts, a transmucosal part made from PEEK polymer (polyether ether ketone) and a fastening (Healing Abutment) or embedded (Healing Cap) screw in titanium alloy (Ti-6Al-4V) manufactured according to ASTM F136. For Healing Abutments models, the fastening screw has a fitting for the 1.2mm hexagonal driver on one end, and a connection for installation, and its other end is fitted with an interface compatible with the prosthetic interface used in Plenum Implants/Components. The Healing Cap is fastened onto the Abutment through a thread inside the Cap and an embedded fastening screw. This screw is fitted with a hexagonal socket and must be used with a 1.2 mm hexagonal driver.

PEEK Components are intended for transitory use and must be installed exclusively in Plenum Implants or Components in order to guide the proper

healing or remodeling of the gingival tissue for the installation of a dental prosthesis. Healing Abutment devices can also be used to digitally read implants during intraoral scanning.

The devices that are part of the PEEK Components family – Healing Cap – are available for use in prosthetic adaptations of Plenum Abutments and Multi Abutments.

Each commercial code for the devices that make up the PEEK Component family consists of the following parts represented in Figures 1A and 1B according to Tables 1A and 1B. The Components made available in the models described in the table below must only be used with Plenum Dental Implants, Prosthetic Components, and Instruments (these Implants, components, and Instruments are not described in these instructions for use and must be purchased separately).

Tables 1A and 1B present illustrative images, codes, models, and main dimensions of the PEEK Components.

**TABLE 1A:** Illustrative images, codes, models, prosthetic platforms, and dimensions of PEEK Components, Healing Abutment models.

Illustrative image	Codes*	Models	Model composition*	Diameter (mm)	Length (mm)	Prosthetic interface
	018-21-C01-35-01	HEALING ABT SL Ø3.5X1mm	HEALING ABT SL Ø3.5X1mm+HEALING ABT SCREW SL 1mm	3.5	1.0	SL
	018-21-C01-35-02	HEALING ABT SL Ø3.5X2mm	HEALING ABT SL Ø3.5X2mm+HEALING ABT SCREW SL 2mm		2.0	
	018-21-C01-35-03	HEALING ABT SL Ø3.5X3mm	HEALING ABT SL Ø3.5X3mm+HEALING ABT SCREW SL 3mm		3.0	
	018-21-C01-35-04	HEALING ABT SL Ø3.5X4mm	HEALING ABT SL Ø3.5X4mm+HEALING ABT SCREW SL 4mm		4.0	
	018-21-C01-35-05	HEALING ABT SL Ø3.5X5mm	HEALING ABT SL Ø3.5X5mm+HEALING ABT SCREW SL 5mm		5.0	
	018-21-C01-35-06	HEALING ABT SL Ø3.5X6mm	HEALING ABT SL Ø3.5X6mm+HEALING ABT SCREW SL 6mm		6.0	
	018-22-C01-35-01	HEALING ABT RE Ø3.5X1mm	HEALING ABT RE Ø3.5X1mm+HEALING ABT SCREW RE 1mm	3.5	1.0	RE
	018-22-C01-35-02	HEALING ABT RE Ø3.5X2mm	HEALING ABT RE 3.5X2mm+HEALING ABT SCREW RE 2mm		2.0	
	018-22-C01-35-03	HEALING ABT RE Ø3.5X3mm	HEALING ABT RE 3.5X3mm+HEALING ABT SCREW RE 3mm		3.0	
	018-22-C01-35-04	HEALING ABT RE Ø3.5X4mm	HEALING ABT RE Ø3.5X4mm+HEALING ABT SCREW RE 4mm		4.0	
	018-22-C01-35-05	HEALING ABT RE Ø3.5X5mm	HEALING ABT RE 3.5X5mm+HEALING ABT SCREW RE 5mm		5.0	
	018-22-C01-35-06	HEALING ABT RE Ø3.5X6mm	HEALING ABT RE 3.5X6mm+HEALING ABT SCREW RE 6 mm		6.0	
	018-22-C01-45-01	HEALING ABT RE Ø4.5X1mm	HEALING ABT RE Ø4.5X1mm+HEALING ABT SCREW RE 1mm	4.5	1.0	
	018-22-C01-45-02	HEALING ABT RE Ø4.5X2mm	HEALING ABT RE Ø4.5X2mm+HEALING ABT SCREW RE 2mm		2.0	
	018-22-C01-45-03	HEALING ABT RE Ø4.5X3mm	HEALING ABT RE Ø4.5X3mm+HEALING ABT SCREW RE 3mm		3.0	
	018-22-C01-45-04	HEALING ABT RE Ø4.5X4mm	HEALING ABT RE 4.5X4mm+HEALING ABT SCREW RE 4mm		4.0	
	018-22-C01-45-05	HEALING ABT RE Ø4.5X5mm	HEALING ABT RE 4.5X5mm+HEALING ABT SCREW RE 5mm		5.0	
	018-22-C01-45-06	HEALING ABT RE Ø4.5X6mm	HEALING ABT RE 4.5X6mm+HEALING ABT SCREW RE 6 mm		6.0	
	018-22-C01-55-01	HEALING ABT RE Ø5.5X1mm	HEALING ABT RE Ø5.5X1mm+HEALING ABT SCREW RE 1mm	5.5	1.0	
	018-22-C01-55-02	HEALING ABT RE Ø5.5X2mm	HEALING ABT RE 5.5X2mm+HEALING ABT SCREW RE 2mm		2.0	
	018-22-C01-55-03	HEALING ABT RE Ø5.5X3mm	HEALING ABT RE 5.5X3mm+HEALING ABT SCREW RE 3mm		3.0	
	018-22-C01-55-04	HEALING ABT RE Ø5.5X4mm	HEALING ABT RE 5.5X4mm+HEALING ABT SCREW RE 4mm		4.0	

	018-22-C01-55-05	HEALING ABT RE Ø5.5X5mm	HEALING ABT RE 5.5X5mm+HEALING ABT SCREW RE 5mm	5.5	5.0	RE
	018-22-C01-55-06	HEALING ABT RE Ø5.5X6mm	HEALING ABT RE 5.5X6mm+HEALING ABT SCREW RE 6 mm		6.0	
	018-23-C01-45-01	HEALING ABT SH Ø4.5X1mm	HEALING ABT SH Ø4.5X1mm+HEALING ABT SCREW SH 1mm	4.5	1.0	SH
	018-23-C01-45-02	HEALING ABT SH Ø4.5X2mm	HEALING ABT SH 4.5X2mm+HEALING ABT SCREW SH 2mm		2.0	
	018-23-C01-45-03	HEALING ABT SH Ø4.5X3mm	HEALING ABT SH 4.5X3mm+HEALING ABT SCREW SH 3mm		3.0	
	018-23-C01-45-04	HEALING ABT SH Ø4.5X4mm	HEALING ABT SH 4.5X4mm+HEALING ABT SCREW SH 4mm		4.0	
	018-23-C01-45-05	HEALING ABT SH Ø4.5X5mm	HEALING ABT SH 4.5X5mm+HEALING ABT SCREW SH 5mm	5.5	5.0	
	018-23-C01-45-06	HEALING ABT SH Ø4.5X6mm	HEALING ABT SH 4.5X6mm+HEALING ABT SCREW RE 6 mm		6.0	
	018-23-C01-55-01	HEALING ABT SH Ø5.5X1mm	HEALING ABT SH Ø5.5X1mm+HEALING ABT SCREW SH 1mm		1.0	
	018-23-C01-55-02	HEALING ABT SH Ø5.5X2mm	HEALING ABT SH 5.5X2mm+HEALING ABT SCREW SH 2mm		2.0	
	018-23-C01-55-03	HEALING ABT SH Ø5.5X3mm	HEALING ABT SH 5.5X3mm+HEALING ABT SCREW SH 3mm		3.0	
	018-23-C01-55-04	HEALING ABT SH Ø5.5X4mm	HEALING ABT SH 5.5X4mm+HEALING ABT SCREW SH 4mm		4.0	
	018-23-C01-55-05	HEALING ABT SH Ø5.5X5mm	HEALING ABT SH 5.5X5mm+HEALING ABT SCREW SH 5mm	5.0		
	018-23-C01-55-06	HEALING ABT SH Ø5.5X6mm	HEALING ABT SH 5.5X6mm+HEALING ABT SCREW RE 6 mm	6.0		

**IMPORTANT:**

\*Each commercial device code consists of the Healing Abutment model (PEEK part) and its respective screw for the exclusive use of the family (according to the column composition of the model). It is important to note that the screw is an integral and exclusive part of the PEEK Healing Abutment Components model, and without such an item the device will not fulfill its intended purpose of use.

The colors of the fastening screw represent the prosthetic interface: PINK interface with SLIM prosthetic components; YELLOW interface with REGULAR prosthetic components; GREEN interface with SHORT prosthetic components.

Healing Abutment Components are compatible with all Plenum Dental Implants, provided that the prosthetic interface corresponds to the implant's internal connection (SL, RE, and SH).

Technical characteristics of screws: Each prosthetic interface model (SL, RE, and SH) has a specific screw type for each component height. Torque is applied to the screw through means of the hexagonal connection.

**TABLE 1B:** Illustrative images, codes, models, prosthetic platforms and main dimensions of PEEK Components, Cap models.

Illustrative image	Codes*	Models	Model composition*	Diameter (mm)	Length (mm)	Prosthetic interface
	018-20-C70-20-00	HEALING CAP MULTI ABT.	HEALING CAP MULTI ABT.	5	5.4	MULTI ABT.
	018-20-C70-50-33	HEALING CAP PLENUM ABT. Ø3.5x4	HEALING CAP PLENUM ABT. Ø3.5x4	3.5	6.5	PLENUM ABT.
	018-20-C70-51-33	HEALING CAP PLENUM ABT. Ø4.5X4	HEALING CAP PLENUM ABT. Ø4.5X4	3.5	8.5	
	018-20-C70-50-45	HEALING CAP PLENUM ABT. Ø3.5X6	HEALING CAP PLENUM ABT. Ø3.5X6	4.5	6.5	
	018-20-C70-51-45	HEALING CAP PLENUM ABT. Ø4.5X6	HEALING CAP PLENUM ABT. Ø4.5X6	4.5	8.5	

**IMPORTANT:**

\*Each commercial device code consists of the PEEK Component Cap model (according to the column model content). The screw positioned at the top of the Cap is fixed.

\*Each Cap interfaces with its respective abutment, as described in the commercial model. Ex.: Healing Cap Multi Abutment for interfacing with Multi Abutments, Healing Cap Plenum Abt. Ø3.5X6 with interface with Ø3.5 Plenum Abutment at a height of 6mm.

## 2. DEVICE COMPOSITION

PEEK Components are composed of two (2) parts: 1. component part, manufactured in PEEK polymer (polyether-ester-ketone), according to the specifications of ASTM F2026 and ABNT NBR 15723-8; 2. Fitted (Healing Abutment) or recessed (Healing Cap) fastening screw made of titanium alloy (Ti-6Al-4V), according to ASTM F136.

**FIGURE 1A.** Schematic model of the parts that make up the PEEK Component, Healing Abutment.

**FIGURE 1B:** Schematic model of the parts that make up PEEK Component models, Healing Caps.

## 3. DEVICE PRESENTATION FORM

PEEK Components in the models described in the table above are supplied in a single package in sterile conditions, and presented in the following manner:

- Primary packaging is a polyethylene terephthalate (PET) blister and thermo-sealed with surgical grade paper; containing 1 unit of the device, sterilized using Gamma Irradiation;
- External packaging is a cardboard box containing the primary packaging and 3 identification and traceability labels;
- It is the responsibility of the hospital/health care facility to ensure complete traceability of the device, through the annotation in the patient's chart of the code and lot of the device used, or identification through one of the 3 traceability labels that come with the package;
- If you want to obtain the Instructions for Use in printed format and at no additional cost (including shipping), request by telephone + 55 11 3109-9000, email ([customerservice@plenum.bio](mailto:customerservice@plenum.bio)), or by the website ([www.plenum.bio/ifu](http://www.plenum.bio/ifu));

- See the illustrative images of the primary and external packaging of the PEEK Components, in the form in which it will be commercialized below.

### ILLUSTRATIVE IMAGE OF THE PRIMARY PACKAGING OF PEEK COMPONENTS (EX.: HEALING ABUTMENT).

### ILLUSTRATIVE IMAGE OF THE EXTERNAL PACKAGING OF PEEK COMPONENTS.

## 4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Healing Abutment is recommended and used during gingival tissue healing after insertion of the Plenum Dental Implant or at the moment of reopening to remove the cover screw during the preparation and installation of a dental prosthesis. These devices in their various models must be installed in exclusive connection with Plenum Dental Implants aiming at the appropriate procedure for placing a dental prosthesis. These components can also be used to digitally read implants during intraoral scanning.

The Healing Cap is a medical device used to protect its respective abutment during dental prosthesis (crown) preparation. Its function is to preserve the adaptation of the abutment (Plenum Abutment or Multi Abutment), preventing contact between the gingival tissue and the prosthetic component.

- **NOTE:** The PEEK Component must only be used with Dental Implant or Prosthetic Components from Plenum (devices that are not part of these instructions for use and must be purchased separately, subject to a specific ANVISA Device Registration Number). Implants from other brands may result in dimensional incompatibility, improper fitting, and micro-movements, and increase the risk of these devices' corrosion (electrolytic or galvanic) and fracturings, compromising the surgical and prosthetic procedure results.

## 5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The PEEK Component in its various models functions to guide, maintain, and facilitate the proper healing of peri-implant gingival tissue after the implant is installed, in reopening surgery (second surgical phase), as well as after installation of the prosthetic component (Abutment or Multi Abutment), shaping the gingival tissue to create space for subsequent dental prosthesis placement procedures. The height of the component body is indicated according to the patient's interocclusal space.

## 6. INSTRUCTIONS FOR USE

- Carefully read all instructions for use, recommendations, and guidelines before using PEEK Components;
- Before rehabilitating the patient, the dental surgeon must carry out rigorous planning through the use of clinical and imaging examinations (radiographic and tomographic) to determine whether the implants installed are in a position to receive prosthetic components;
- Before starting clinical and surgical procedures, the dental surgeon must subject the patient to a rigorous oral antisepsis to prevent the component from coming into contact with contaminating substances;
- Before the use of PEEK Components, the dental surgeon must strictly verify that the device uses the same prosthetic interface for the implant/prosthetic component that is installed, according to the table;
- Before starting the clinical procedure, the dental surgeon must make sure that the prosthetic components and auxiliary instruments are intact and complete;
- The clinical/surgical procedure, the model and size of the component, the technique, and the sequence of installation of PEEK Components must comply with established methods for implant dentistry;
- To use PEEK Components in a single-stage procedure (esthetic and/or immediate loading), the dental surgeon must respect the occlusion principles established in implant dentistry and use compatible restorative materials to ensure that the procedure has satisfactory results;
- To use PEEK Components in two-stage procedures, it is recommended that peri-implant tissues be prepared in advance with the use of this component;
- Instruments manufactured by Plenum must be exclusively used with PEEK Components, which are considered auxiliary instruments and are not part of these instructions for use and must be purchased separately (see Auxiliary Instruments);
- The handling of the device must be carried out with care, especially during the process of fastening the component onto the dental implant. The Component must be passively seated onto the implant platform, thereby avoiding a forced fit, the use of excessive force, or oblique insertion;
- It is highly recommended that imaging tests be used to certify the proper seating of prosthetic components onto the implant platform;
- Special care must be taken with the component's hexagonal connection, and the driver used to introduce the part, as an unsuitable fit between the two elements may result in damage to the device's connection, compromising its integrity and proper fastening to the implant;
- Care must be taken with the surface of the screw thread and the

Morse Taper (MT) connection in the prosthetic component. Both surfaces must be free of particles and solid materials that can prevent perfect fitting and fastening to the implant and provisional dental prosthesis;

- The dental surgeon must carefully assess the repeated introduction and removal of a prosthetic component in the same patient, as this may result in deterioration of the threads, the Morse Taper, and the device's connection;
- Check the passivity and perform the occlusal and interproximal adjustment after the prosthesis's installation, avoiding compromising the implant/prosthesis set.

**ATTENTION:** The figures are merely illustrative. They do not represent the actual dimensions and characteristics of the device.

## 7. PRECAUTIONS

The dental surgeon should not use this device before carefully reading all the Instructions for Use and information on the packaging. Failure to observe these points may lead to complications for the dental surgeon and especially for the patient;

- The PEEK Component must be discarded after the period of healing of the gingival tissue and or digital transfer (scanning) of the implant at which it was adapted;
- The PEEK Component family is primarily used in tissue preparation and installation of a dental prosthesis, and dentists are responsible for selecting the device models with the most appropriate anatomical and dimensional characteristics, after an individual clinical analysis of each case;
- Postoperative, immediate, and subsequent care and proper oral hygiene techniques are extremely important. The patient must be given instructions regarding the existing risks for maintenance of implants and components; The patient's ability and willingness to follow instructions from the dental surgeon is one of the most important aspects of the dental implant placement procedure;
- Patient must be warned that a failure to comply with instructions for the post-operative period may lead to failure and early loss of the component and implant, requiring additional surgery to remove the implant or replace the component;
- For the safety and efficacy of the treatment, only prosthetic components and implants manufactured by Plenum should be used;
- Components and/or implants from other brands may result in dimensional incompatibility, improper fitting, micro-movements, and increase the risk of corrosion (electrolytic or galvanic) and fracturing of these devices,

compromising the results of the surgical and prosthetic procedure.

- Before use, special attention should be paid to PEEK Component's insertion drivers. Damage to the tip of the drivers can prevent the correct component installation, compromising the good prosthetic result. Therefore, this device must not be used if there is any damaged instrument for specific use;
- Improper selection of the component that will be used, as well as errors in the indicated use, handling and installation of the device may result in damage to the device's structure, which may contribute to failure due to fatigue, fracturing, or loosening of the component;
- The use of excessive torque during the insertion of a prosthetic component may lead to device failure. Do not exceed 10 Ncm of insertion torque for the fastening screw;
- If a prosthetic component and/or the implant is overloaded beyond its functional capacity, excessive bone loss and breakdown of devices and other restoration elements may occur;
- Do not force the introduction of a prosthetic component beyond the depth established by the thread limit, as this may cause deformation in the hexagonal connection of the component and/or the insertion driver and fracture of the device;
- Incorrect handling of components of small sizes inside the patient's mouth may incur a risk of aspiration and/or swallowing of these elements;
- The dental surgeon must rigorously evaluate the initial mechanical stability of the prosthetic component. If it has not obtained ideal and safe mechanical stability, the dental surgeon must evaluate and consider a possible removal and/or replacement thereof;
- The patient must be informed of the risks involved in the application of excessive mechanical load or trauma that affects the operated region, which may imply failure, loss of mechanical stability, loosening, breakage, fracture or loss of the dental implant and its components;
- Elderly patients, patients with mental health problems, individuals with chemical dependency, or any other difficulty in adhering to postoperative recommendations may represent a higher risk of failure of this device since these patients may ignore instructions for this procedure;
- Smokers, users of tobacco and drugs, alcoholics, diabetics, or malnourished patients must be warned about the increased incidence of lack of osseointegration, with consequent failure of implants and components. Additionally, for this group of patients, the rate of insufficient healing of the gingival tissue around the body of the component is considerable and may also be inefficient, compromising the installation of the dental prosthesis;
- If there are visible deformations, scratches, cracks, or dents in the component, it must be discarded;
- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the installation procedure for the dental implant set and prosthetic component;
- The dental surgeon must guide the patient, who must understand the importance of follow-up during the postoperative period, in which the stability, integrity, and positioning of the prosthetic component and dental implant used must be verified, as well as the osseointegration and biological fastening of the implant at the appropriate time;
- To protect the patient and their surgeon, it is important that the correct and complete traceability of the component used be guaranteed by noting the device code and lot in the patient's records or by attaching to the identification and traceability labels that come with the device to these records;
- The success of a procedure for placing a prosthetic component is related

to the correct indication, selection, resistance of the implant and component, correct surgical technique, as well as the quality of osseointegration, which can be affected by numerous biological, clinical, and biomechanical factors, which are external to the quality of the implants and components. These factors may limit the useful life of the devices and consequently cause their failure and early loss. Strict adherence to recommendations, contraindications, precautions, warnings, and other observations contained in these instructions for use are, therefore, essential for the procedure's success.

## 8. WARNING

- Single-use only do not resterilize or reuse.
- Do not use if the package is opened or damaged or if the expiration date has been exceeded.

## 9. CONTRAINDICATIONS

Peek Components are always used with Plenum Implants (provided separately). Dentists must therefore assess the patient's condition prudently before surgery and use of the components and check for any contraindications reported below. Contraindications may include, but are not limited to:

- An implant that is poorly positioned and/or presents insufficient osseointegration to support or anchor a dental prosthesis and its components, or is subject to strong restrictions on the use of this device;
- Patients with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the dental implant and component due to lack of osseointegration;
- Patients that present signs of allergy or hypersensitivity to the composition of the raw materials used in the manufacturing of PEEK Components (see "DEVICE COMPOSITION");
- Insufficient interocclusal or mesiodistal space for insertion of any PEEK Component model;
- Inflammation and/or periodontal infection;
- Recent history of systemic or localized infection;
- Parafunctional habits;
- Patient with insufficient intermandibular space;
- Non-treatable occlusion and joint problems;
- Pregnancy and lactation period;
- Use in patients with immunological changes and systemic diseases that result in poor healing or a lack of integrity in existing tissue at the implant site and gingival tissue;
- Alcoholic patient or drug user who cannot correctly follow postoperative instructions;
- Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;
- Patient unable to maintain good oral hygiene;
- All contraindications included in known dental/scientific literature for using dental implants and their components.

## 10. POSSIBLE ADVERSE EFFECTS

Adverse effects are not expected if the device is used by its instructions for use.

However, all surgical and prosthetic procedures may not achieve positive results. Every procedure presents risks and the possibility of complications, such as infections, allergic reactions to drugs, among others. Possible adverse effects associated with the use of prosthetic components may include, but are not limited to:


- Pain or discomfort due to the presence of the prosthetic component and/or dental implant or due to the surgical procedure;
- Lack of osseointegration, reabsorption of bone tissue, or even delay in bone healing, which can result in the failure and breakage of the prosthetic component and/or dental implant;
- Fracture, displacement, or early or late failure of the component and/or implant due to improper installation technique;
- Failure, displacement, or fracture of the component or implant due to excessive loads applied to the devices and/or inadequate bone fastening;
- Loss of the component or implant due to parafunctional habits or poor oral hygiene;
- Allergic reaction to the raw materials used in PEEK Components (see item 2 - "DEVICE COMPOSITION");
- Inefficient healing process in gingival tissue due to the imbalance of normal forces/stress on the implant/component interface or the emergence of systemic disease(s) after installation of the component, which the patient did not report to the dental surgeon;
- Bone necrosis, insufficient vascularization, bone reabsorption, and insufficient bone formation, which can cause the loss of fastening of the dental implant and consequently the component; these effects are correlated to the patient's systemic facts or surgical procedure that is inappropriate without complying with the instructions for use described for this family of components;
- The patient must be informed and understand that complications in all surgical procedures are possible. The dental surgeon must inform the patient of these complications and the risks of implant surgery.

## 11. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

**TABLE 2** - Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code

	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not re-sterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)

**NOTE:** The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

## 12. AUXILIARY INSTRUMENTS

- All surgical instruments developed exclusively for the installation of PEEK Components are considered to be auxiliary instruments;
- Plenum Instruments are manufactured in stainless steel (AISI 420-10), the raw materials of which provide devices with high strength and durability, according to requirements specified by NBR 13911; ASTM F899, and ISO 5832-1;
- The use of instruments manufactured by Plenum is essential to ensuring safety and effectiveness in the installation of implants according to the surgical protocol recommended by the manufacturer;
- It is essential that the surgeon exclusively use the introducer drivers manufactured by Plenum to avoid damage to the connection of the components and also to introduction drivers;
- For the proper installation of the implants, use only the surgical instruments manufactured by Plenum, so that there is no incongruence between the fittings and irreversible damage to the components.
- A list of surgical instruments can be found in printed or electronic communication materials or on Plenum's website ([www.plenum.bio](http://www.plenum.bio)). Plenum instruments must be used in accordance with the manufacturer's instructions;
- Surgical instruments are not the primary purpose of these instructions for use, that is, they have separate registrations with ANVISA. Therefore, they must be purchased separately.

## 13. STERILIZATION

- PEEK Components are supplied in sterile condition and sterilized using Gamma Irradiation;
- Single-use only; do not reprocess, resterilize or reuse;
- Before use, check the integrity of the packaging and the expiration date. Do not use the device if the package is opened or damaged or if the expiration date has been exceeded;
- Date of manufacture, sterilization expiration date, device code, batch code, and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments;

The opening of the packaging for the use of the device must only be done by qualified and trained personnel to perform this procedure.

## 14. EXPIRATION DATE

- Indicated on device labeling.
- Before use, check the expiration date, do not use the device if it has been exceeded.

## 15. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of PEEK Components and comply with health surveillance requirements, the dental surgeon or someone on their team must keep information about the component used in the patient's medical record. Additionally, this information must also be passed on to the device distributor and patient to complete the traceability cycle. The information necessary for traceability is related to the device used and surgical procedure as follows:

- Name of the dental surgeon;
- Name of hospital/clinic/office;
- Name of manufacturer and distributor Date of surgery;
- Device code;
- Batch code;
- Quantity used;
- ANVISA device registration number.

The responsible dental surgeon and their team must use the information in the 3 traceability labels supplied with PEEK Components, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that they have information about the component used in their procedure.

The labels contain the following information necessary for the traceability of the PEEK Component used.

- Manufacturer's identification;
- Device name or commercial model Device code;
- Batch code;
- ANVISA device registration number.

Traceability information is necessary for notification by the health service and/or by the patient to the National Health Surveillance Agency – ANVISA and to the manufacturer, when complaints, technical complaints and/or adverse events occur, for the conduct of the appropriate investments.

## 16. STORAGE, TRANSPORT, CONSERVATION AND HANDLING

- The packaging must be intact at the time of receipt, do not use the device if the packaging is damaged or if the expiration date has been exceeded. If the packaging is damaged, the device must be considered non-sterile and discarded;
- The device must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight, and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The device must be stored on shelving or drawers/cabinets to allow it to be easily viewed and handled by the operator. Do not store the device near lamps, so that the information contained in the packaging is not affected;
- The device should only be opened at the time of use in a sterile surgical field;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

## 17. DISPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After use of the device or after it is removed from the patient, PEEK Components that are not subject to analysis or subsequent studies must be disposed of by the health service establishment (hospitals, offices, and clinics);
- It is recommended that Provisional devices must be completely disassembled to ensure they cannot be reused. The health service institution/ establishment is solely responsible for disassembling the device;
- After disassembly, the instrument must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- PEEK Components that are not used during surgery must also be discarded following the guidelines above.

## 18. CUSTOMER SERVICES

- If the device is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, device code, and lot and patient data, directly notifying Plenum's Customer Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and correctly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Phone: + 55 11 3109-9000; Email: [customerservice@plenum.bio](mailto:customerservice@plenum.bio).
- ANVISA also makes a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found on the website of the National Health Regulatory Agency - ANVISA at the electronic address [www.anvisa.gov.br](http://www.anvisa.gov.br), using the NOTIVISA link.
- FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

## 19. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. The occurrence of defects in the device must be reported immediately to the manufacturer, observing the legal deadline and attaching a technical report with the details of the non-compliance with the device. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see item 17) must be observed in cases involving a risk to the patient's health.

The Warranty of devices manufactured by Plenum is conditioned to the

strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and follow-up exempts the manufacturer and distributor of the devices from any kind of liability, in particular in the following cases:

- Reprocessing, resterilization, or reuse of the device by the user;
- Improper use or use other than in a surgical setting;
- Device not used according to instructions for use.

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