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DENTAL IMPLANT
English

IMPLANTE DENTAL
Español

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Comprimento do implante:	5mm
Diâmetro do implante	Protocolo cirúrgico
Ø3,0	X
Ø3,5	X
Ø4,0	X
Ø5,0	Ø2,0x5 - Ø3,5x5 - Ø5,0x5 - Macho 5,0*
Ø6,0	Ø2,0x5 - Ø3,5x5 - Ø5,0x5 - Ø6,0x5 - Macho 6,0*

Comprimento do implante:	8mm
Diâmetro do implante	Protocolo cirúrgico
Ø3,0	X
Ø3,5	Ø2,0x5 - Ø3,5x8 - Macho 3,5**
Ø4,0	Ø2,0x5 - Ø3,5x8 - Macho 3,5**
Ø5,0	Ø2,0x5 - Ø3,5x8 - Ø4,0x8 - Ø5,0x8 - Macho 5,0**
Ø6,0	X

Comprimento do implante:	10mm
Diâmetro do implante	Protocolo cirúrgico
Ø3,0	Ø2,0x5 - Ø3,0x10 - Macho 3,0**
Ø3,5	Ø2,0x5 - Ø3,5x10 - Macho 3,5**
Ø4,0	Ø2,0x5 - Ø3,5x10 - Ø4,0x10 - Macho 4,0**
Ø5,0	Ø2,0x5 - Ø3,5x8 - Ø4,0x8 - Ø5,0x10 - Macho 5,0**
Ø6,0	X

Comprimento do implante:	11,5mm*
Diâmetro do implante	Protocolo cirúrgico
Ø3,0	Ø2,0x5 - Ø3,0x13 - Macho 3,0**
Ø3,5	Ø2,0x5 - Ø3,5x13 - Macho 3,5**
Ø4,0	Ø2,0x5 - Ø3,5x10 - Ø4,0x13 - Macho 4,0**
Ø5,0	Ø2,0x5 - Ø3,5x8 - Ø4,0x8 - Ø5,0x13 - Macho 5,0**
Ø6,0	X

Comprimento do implante:	13mm
Diâmetro do implante	Protocolo cirúrgico
Ø3,0	Ø2,0x5 - Ø3,0x13 - Macho 3,0**
Ø3,5	Ø2,0x5 - Ø3,5x13 - Macho 3,5**
Ø4,0	Ø2,0x5 - Ø3,5x10 - Ø4,0x13 - Macho 4,0**
Ø5,0	Ø2,0x5 - Ø3,5x8 - Ø4,0x8 - Ø5,0x13 - Macho 5,0**
Ø6,0	X

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TECHNICAL NAME: Dental Implants (Osseointegrable)

TRADE NAME: Dental Implant

MODEL: For a detailed description of the device, see Table 1.

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Website: www.plenum.bio

ATTENTION:

Carefully read these Instructions for Use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile Device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited to Reprocess, resterilize, and reuse.

ANVISA DEVICE REGISTRATION NO.: 81684340014.

Device code, batch code, and date of manufacture and sterilization: See device labeling.

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	6
3. DEVICE PRESENTATION FORM	6
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	8
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	8
6. INSTRUCTIONS FOR USE	8
7. PRECAUTIONS	13
8. WARNINGS	14
9. CONTRAINDICATIONS	14
10. POSSIBLE ADVERSE EFFECTS	14
11. SYMBOL DESCRIPTIONS	14
12. AUXILIARY INSTRUMENTS	15
13. STERILIZATION	15
14. EXPIRATION DATE	15
15. IDENTIFICATION AND TRACEABILITY	16
16. STORAGE, TRANSPORT, AND HANDLING CONDITIONS	16
17. DISPOSAL OF DEVICE	16
18. CUSTOMER SERVICES	16
19. LIMITED WARRANTY	16

1. DETAILED DEVICE DESCRIPTION

The Dental Implant is an implantable medical device manufactured in titanium alloy (Ti-6Al-4V ELI). The main functions of this device are to restore the masticatory function and aesthetics of the patient by replacing the dental root and supporting a dental prosthesis in cases of total, partial, and unitary restorations, eliminating the need for the use of removable prosthesis (dentures).

The Dental Implant has a cylindrical-conical external body with a Morse Taper connection, a double-inlet trapezoidal thread profile, and a self-threading tip.

Macro and micro-structured surface topography enables bone cell adhesion, allowing effective osseointegration.

The Dental Implant must be used in connection with other Plenum abutments and dental prosthesis.

IMPORTANT: To use the Dental Implant, it is necessary to use Components manufactured by Plenum, which are not part of this family of devices and must be purchased separately.

The Cover Screw is part of this line and must be used temporarily with the Dental Implant to close the implant platform until the dental prosthesis is placed.

Table 1 shows illustrative images of the devices, codes, models, and main dimensions of Dental Implants.

TABLE 1. Illustrative images of Dental Implants, codes, models, and dimensions.

Prosthetic Interface	SL		
Model Composition*	SLIM Implant+ SLCover Screw 1mm		
Illustrative Image			
*Codes	Models	Diameter (mm)	Length (mm)
018-10-110-30-08	SLIM Ø3.0 X 8MM	3.0	8.0
018-10-110-30-10	SLIM Ø3.0 X 10MM		10.0
018-10-110-30-11	SLIM Ø3.0 X 11.5MM		11.5
018-10-110-30-13	SLIM Ø3.0 X 13MM		13.0
018-10-110-30-15	SLIM Ø3.0 X 15MM		15.0

IMPORTANT:

- The interior colors represent the prosthetic interface: PINK for interfacing with SLIM Prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT prosthetic components.

- Dental Implants belonging to this family are compatible with all Prosthetic Components, as long as the prosthetic interface column is matched.

*Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

Prosthetic Interface	RE		
Model Composition*	REGULAR Implant + RE Cover Screw 1mm		
Illustrative Image			
*Codes	Models	Diameter (mm)	Length (mm)
018-10-110-35-08	REGULAR Ø3.5 X 8MM	3.5	8
018-10-110-35-10	REGULAR Ø3.5 X 10MM		10
018-10-110-35-11	REGULAR Ø3.5 X 11.5MM		11.5
018-10-110-35-13	REGULAR Ø3.5 X 13MM		13
018-10-110-35-15	REGULAR Ø3.5 X 15MM		15
018-10-110-35-17	REGULAR Ø3.5X17MM		17
018-10-110-35-19	REGULAR Ø3.5X19MM		19
018-10-110-35-21	REGULAR Ø3.5X21MM	21	
018-10-110-35-23	REGULAR Ø3.5X23MM	23	
018-10-110-37-08	REGULAR Ø3.75 X 8MM	3.75	8
018-10-110-37-10	REGULAR Ø3.75 X 10MM		10
018-10-110-37-11	REGULAR Ø3.75 X 11.5MM		11.5
018-10-110-37-13	REGULAR Ø3.75 X 13MM		13
018-10-110-37-15	REGULAR Ø3.75 X 15MM		15
018-10-110-37-17	REGULAR Ø3.75X17MM		17
018-10-110-37-19	REGULAR Ø3.75X19MM		19
018-10-110-37-21	REGULAR Ø3.75X21MM	21	
018-10-110-37-23	REGULAR Ø3.75X23MM	23	
018-10-110-40-08	REGULAR Ø4.0 X 8MM	4.0	8
018-10-110-40-10	REGULAR Ø4.0 X 10MM		10
018-10-110-40-11	REGULAR Ø4.0 X 11.5MM		11.5
018-10-110-40-13	REGULAR Ø4.0 X 13MM		13
018-10-110-40-15	REGULAR Ø4.0 X 15MM		15
018-10-110-40-17	REGULAR Ø4.0X17MM		17
018-10-110-40-19	REGULAR Ø4.0X19MM		19
018-10-110-40-21	REGULAR Ø4.0X21MM	21	
018-10-110-40-23	REGULAR Ø4.0X23MM	23	
018-10-110-41-08	REGULAR Ø4.1 X 8MM	4.1	8
018-10-110-41-10	REGULAR Ø4.1 X 10MM		10
018-10-110-41-11	REGULAR Ø4.1 X 11.5MM		11.5
018-10-110-41-13	REGULAR Ø4.1 X 13MM		13
018-10-110-41-15	REGULAR Ø4.1 X 15MM		15
018-10-110-45-08	REGULAR Ø4.5 X 8MM	4.5	8
018-10-110-45-10	REGULAR Ø4.5 X 10MM		10
018-10-110-45-11	REGULAR Ø4.5 X 11.5MM		11.5
018-10-110-45-13	REGULAR Ø4.5 X 13MM		13
018-10-110-45-15	REGULAR Ø4.5 X 15MM		15
018-10-110-50-08	REGULAR Ø5.0 X 8MM	5.0	8
018-10-110-50-10	REGULAR Ø5.0 X 10MM		10
018-10-110-50-11	REGULAR Ø5.0 X 11.5MM		11.5
018-10-110-50-13	REGULAR Ø5.0 X 13MM		13
018-10-110-55-08	REGULAR Ø5.5 X 8MM	5.5	8
018-10-110-55-10	REGULAR Ø5.5 X 10MM		10
018-10-110-55-11	REGULAR Ø5.5 X 11.5MM		11.5
018-10-110-55-13	REGULAR Ø5.5 X 13MM		13
018-10-110-60-08	REGULAR Ø6.0 X 8MM	6.0	8
018-10-110-60-10	REGULAR Ø6.0 X 10MM		10
018-10-110-60-11	REGULAR Ø6.0 X 11.5MM		11.5
018-10-110-60-13	REGULAR Ø6.0 X 13MM		13

IMPORTANT:

- The interior colors represent the prosthetic interface: PINK for interfacing with SLIM prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT prosthetic components.

- Dental Implants belonging to this family are compatible with all Prosthetic Components, as long as the prosthetic interface column is matched.

* Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

Prosthetic Interface	SH		
Model Composition*	SHORT Implant + SH Cover Screw 1mm		
Illustrative Image			
*Codes	Models	Diameter (mm)	Length (mm)
018-10-I10-37-06	SHORT Ø3.75 X 6MM	3.75	6
018-10-I10-37-07	SHORT Ø3.75 X 7MM		7
018-10-I10-40-06	SHORT Ø4.0 X 6MM	4.0	6
018-10-I10-40-07	SHORT Ø4.0 X 7MM		7
018-10-I10-41-06	SHORT Ø4.1 X 6MM	4.1	6
018-10-I10-41-07	SHORT Ø4.1 X 7MM		7
018-10-I10-45-05	SHORT Ø4.5 X 5MM	4,5	5
018-10-I10-45-06	SHORT Ø4.5 X 6MM		6
018-10-I10-45-07	SHORT Ø4.5 X 7MM		7
018-10-I10-50-05	SHORT Ø5.0 X 5MM	5.0	5
018-10-I10-50-06	SHORT Ø5.0 X 6MM		6
018-10-I10-50-07	SHORT Ø5.0 X 7MM		7
018-10-I10-55-05	SHORT Ø5.5 X 5MM	5.5	5
018-10-I10-55-06	SHORT Ø5.5 X 6MM		6
018-10-I10-55-07	SHORT Ø5.5 X 7MM		7
018-10-I10-60-05	SHORT Ø6.0 X 5MM	6.0	5
018-10-I10-60-06	SHORT Ø6.0 X 6MM		6
018-10-I10-60-07	SHORT Ø6.0 X 7MM		7

IMPORTANT:

The interior colors represent the prosthetic interface: PINK for interfacing with SLIM prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT prosthetic components.

Dental Implants belonging to this family are compatible with all Prosthetic Components, provided that the prosthetic interface column is matched.

*Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

Prosthetic Interface	SH		
Model Composition*	SX Implant + SH Cover Screw 1mm		
Illustrative Image			
*Codes	Models	Diameter (mm)	Length (mm)
018-10-I11-41-05	SX Ø4.1X5MM	4.1	5.0
018-10-I11-41-06	SX Ø4.1X6MM		6.0
018-10-I11-41-07	SX Ø4.1X7MM		7.0
018-10-I11-41-08	SX Ø4.1X8MM		8.0
018-10-I11-41-10	SX Ø4.1X10MM		10.0
018-10-I11-41-11	SX Ø4.1X11.5MM		11.5
018-10-I11-41-13	SX Ø4.1X13MM		13.0
018-10-I11-41-15	SX Ø4.1X15MM		15.0
018-10-I11-41-17	SX Ø4.1X17MM		17.0
018-10-I11-41-19	SX Ø4.1X19MM		19.0
018-10-I11-41-21	SX Ø4.1X21MM		21.0
018-10-I11-41-23	SX Ø4.1X23MM		23.0
018-10-I11-45-05	SX Ø4.5X5MM		4.5
018-10-I11-45-06	SX Ø4.5X6MM	6.0	
018-10-I11-45-07	SX Ø4.5X7MM	7.0	
018-10-I11-45-08	SX Ø4.5X8MM	8.0	
018-10-I11-45-10	SX Ø4.5X10MM	10.0	
018-10-I11-45-11	SX Ø4.5X11.5MM	11.5	
018-10-I11-45-13	SX Ø4.5X13MM	13.0	
018-10-I11-45-15	SX Ø4.5X15MM	15.0	
018-10-I11-45-17	SX Ø4.5X17MM	17.0	
018-10-I11-45-19	SX Ø4.5X19MM	19.0	
018-10-I11-45-21	SX Ø4.5X21MM	21.0	
018-10-I11-45-23	SX Ø4.5X23MM	23.0	
018-10-I11-50-05	SX Ø5.0X5MM	5.0	
018-10-I11-50-06	SX Ø5.0X6MM		6.0
018-10-I11-50-07	SX Ø5.0X7MM		7.0
018-10-I11-50-08	SX Ø5.0X8MM		8.0
018-10-I11-50-10	SX Ø5.0X10MM		10.0
018-10-I11-50-11	SX Ø5.0X11.5MM		11.5
018-10-I11-50-13	SX Ø5.0X13MM		13.0
018-10-I11-50-15	SX Ø5.0X15MM		15.0
018-10-I11-50-17	SX Ø5.0X17MM		17.0
018-10-I11-50-19	SX Ø5.0X19MM		19.0
018-10-I11-50-21	SX Ø5.0X21MM		21.0
018-10-I11-50-23	SX Ø5.0X23MM		23.0
018-10-I11-55-05	SX Ø5.5X5MM		5.5
018-10-I11-55-06	SX Ø5.5X6MM	6.0	
018-10-I11-55-07	SX Ø5.5X7MM	7.0	
018-10-I11-55-08	SX Ø5.5X8MM	8.0	
018-10-I11-55-10	SX Ø5.5X10MM	10.0	
018-10-I11-55-11	SX Ø5.5X11.5MM	11.5	

018-10-I11-55-13	SX Ø5.5X13MM	5.5	13.0
018-10-I11-55-15	SX Ø5.5X15MM		15.0
018-10-I11-55-17	SX Ø5.5X17MM		17.0
018-10-I11-55-19	SX Ø5.5X19MM		19.0
018-10-I11-55-21	SX Ø5.5X21MM		21.0
018-10-I11-55-23	SX Ø5.5X23MM		23.0
018-10-I11-60-05	SX Ø6.0X5MM	6.0	5.0
018-10-I11-60-06	SX Ø6.0X6MM		6.0
018-10-I11-60-07	SX Ø6.0X7MM		7.0
018-10-I11-60-08	SX Ø6.0X8MM		8.0
018-10-I11-60-10	SX Ø6.0X10MM		10.0
018-10-I11-60-11	SX Ø6.0X11.5MM		11.5
018-10-I11-60-13	SX Ø6.0X13MM		13.0
018-10-I11-60-15	SX Ø6.0X15MM		15.0
018-10-I11-60-17	SX Ø6.0X17MM		17.0
018-10-I11-60-19	SX Ø6.0X19MM		19.0
018-10-I11-60-21	SX Ø6.0X21MM		21.0
018-10-I11-60-23	SX Ø6.0X23MM		23.0

IMPORTANT:

The interior colors represent the prosthetic interface: PINK for interfacing with SLIM Prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT and SX prosthetic components.

Dental Implants belonging to this family are compatible with all Prosthetic Components, as long as the prosthetic interface column is matched.

*Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

Prosthetic Interface	UF		
Model Composition*	UNIFIT Implant		
Illustrative image			
*Codes	Models	Diameter (mm)	Length (mm)
018-10-I18-35-08	UNIFIT Ø3.5X8MM	3.5	8.0
018-10-I18-35-10	UNIFIT Ø3.5X10MM		10.0
018-10-I18-35-11	UNIFIT Ø3.5X11.5MM		11.5
018-10-I18-35-13	UNIFIT Ø3.5X13MM		13.0
018-10-I18-35-15	UNIFIT Ø3.5X15MM		15.0
018-10-I18-40-08	UNIFIT Ø4.0X8MM		4.0
018-10-I18-40-10	UNIFIT Ø4.0X10MM	10.0	
018-10-I18-40-11	UNIFIT Ø4.0X11.5MM	11.5	
018-10-I18-40-13	UNIFIT Ø4.0X13MM	13.0	
018-10-I18-40-15	UNIFIT Ø4.0X15MM	15.0	

018-10-I18-45-08	UNIFIT Ø4.5X8MM	4.5	8.0
018-10-I18-45-10	UNIFIT Ø4.5X10MM		10.0
018-10-I18-45-11	UNIFIT Ø4.5X11.5MM		11.5
018-10-I18-45-13	UNIFIT Ø4.5X13MM		13.0
018-10-I18-45-15	UNIFIT Ø4.5X15MM		15.0
018-10-I18-50-08	UNIFIT Ø5.0X8MM	5.0	8.0
018-10-I18-50-10	UNIFIT Ø5.0X10MM		10.0
018-10-I18-50-11	UNIFIT Ø5.0X11.5MM		11.5
018-10-I18-50-13	UNIFIT Ø5.0X13MM		13.0
018-10-I18-60-08	UNIFIT Ø6.0X8MM	6.0	8.0
018-10-I18-60-10	UNIFIT Ø6.0X10MM		10.0

IMPORTANT:

The interior colors represent the prosthetic interface: PINK for interfacing with SLIM prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT and SX prosthetic components; GRAY for interfacing with UNIFIT prosthetic components and BLUE for interfacing with MORSEFIT prosthetic components

Dental Implants belonging to this family are compatible with all Prosthetic Components, provided that the prosthetic interface column is matched.

*Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

Prosthetic Interface	MF			
Model Composition*	MF Implant			
Illustrative image				
*Codes	Models	Diameter (mm)	Length (mm)	
018-10-I19-35-08	MORSEFIT Ø3.5X8MM	3.5	8.0	
018-10-I19-35-10	MORSEFIT Ø3.5X10MM		10.0	
018-10-I19-35-11	MORSEFIT Ø3.5X11.5MM		11.5	
018-10-I19-35-13	MORSEFIT Ø3.5X13MM		13.0	
018-10-I19-35-15	MORSEFIT Ø3.5X15MM		15.0	
018-10-I19-35-17	MORSEFIT Ø3.5X17MM		17.0	
018-10-I19-35-19	MORSEFIT Ø3.5X19MM		19.0	
018-10-I19-35-21	MORSEFIT Ø3.5X21MM		21.0	
018-10-I19-35-23	MORSEFIT Ø3.5X23MM		23.0	
018-10-I19-38-08	MORSEFIT Ø3.8X8MM		3.8	8.0
018-10-I19-38-10	MORSEFIT Ø3.8X10MM			10.0
018-10-I19-38-11	MORSEFIT Ø3.8X11.5MM			11.5
018-10-I19-38-13	MORSEFIT Ø3.8X13MM			13.0
018-10-I19-38-15	MORSEFIT Ø3.8X15MM	15.0		
018-10-I19-38-17	MORSEFIT Ø3.8X17MM	17.0		
018-10-I19-38-19	MORSEFIT Ø3.8X19MM	19.0		
018-10-I19-38-21	MORSEFIT Ø3.8X21MM	21.0		
018-10-I19-38-23	MORSEFIT Ø3.8X23MM	23.0		

018-10-119-41-08	MORSEFIT Ø4.1X8MM		8.0
018-10-119-41-10	MORSEFIT Ø4.1X10MM		10.0
018-10-119-41-11	MORSEFIT Ø4.1X11.5MM		11.5
018-10-119-41-13	MORSEFIT Ø4.1X13MM		13.0
018-10-119-41-15	MORSEFIT Ø4.1X15MM	4.1	15.0
018-10-119-41-17	MORSEFIT Ø4.1X17MM		17.0
018-10-119-41-19	MORSEFIT Ø4.1X19MM		19.0
018-10-119-41-21	MORSEFIT Ø4.1X21MM		21.0
018-10-119-41-23	MORSEFIT Ø4.1X23MM		23.0
018-10-119-45-08	MORSEFIT Ø4.5X8MM		8.0
018-10-119-45-10	MORSEFIT Ø4.5X10MM		10.0
018-10-119-45-11	MORSEFIT Ø4.5X11.5MM		11.5
018-10-119-45-13	MORSEFIT Ø4.5X13MM		13.0
018-10-119-45-15	MORSEFIT Ø4.5X15MM	4.5	15.0
018-10-119-45-17	MORSEFIT Ø4.5X17MM		17.0
018-10-119-45-19	MORSEFIT Ø4.5X19MM		19.0
018-10-119-45-21	MORSEFIT Ø4.5X21MM		21.0
018-10-119-45-23	MORSEFIT Ø4.5X23MM		23.0
018-10-119-50-08	MORSEFIT Ø5.0X8MM		8.0
018-10-119-50-10	MORSEFIT Ø5.0X10MM		10.0
018-10-119-50-11	MORSEFIT Ø5.0X11.5MM		11.5
018-10-119-50-13	MORSEFIT Ø5.0X13MM		13.0
018-10-119-50-15	MORSEFIT Ø5.0X15MM	5.0	15.0
018-10-119-50-17	MORSEFIT Ø5.0X17MM		17.0
018-10-119-50-19	MORSEFIT Ø5.0X19MM		19.0
018-10-119-50-21	MORSEFIT Ø5.0X21MM		21.0
018-10-119-50-23	MORSEFIT Ø5.0X23MM		23.0
018-10-119-55-08	MORSEFIT Ø5.5X8MM		8.0
018-10-119-55-10	MORSEFIT Ø5.5X10MM		10.0
018-10-119-55-11	MORSEFIT Ø5.5X11.5MM		11.5
018-10-119-55-13	MORSEFIT Ø5.5X13MM		13.0
018-10-119-55-15	MORSEFIT Ø5.5X15MM	5.5	15.0
018-10-119-55-17	MORSEFIT Ø5.5X17MM		17.0
018-10-119-55-19	MORSEFIT Ø5.5X19MM		19.0
018-10-119-55-21	MORSEFIT Ø5.5X21MM		21.0
018-10-119-55-23	MORSEFIT Ø5.5X23MM		23.0
018-10-119-60-08	MORSEFIT Ø6.0X8MM		8.0
018-10-119-60-10	MORSEFIT Ø6.0X10MM		10.0
018-10-119-60-11	MORSEFIT Ø6.0X11.5MM		11.5
018-10-119-60-13	MORSEFIT Ø6.0X13MM		13.0
018-10-119-60-15	MORSEFIT Ø6.0X15MM	6.0	15.0
018-10-119-60-17	MORSEFIT Ø6.0X17MM		17.0
018-10-119-60-19	MORSEFIT Ø6.0X19MM		19.0
018-10-119-60-21	MORSEFIT Ø6.0X21MM		21.0
018-10-119-60-23	MORSEFIT Ø6.0X23MM		23.0

IMPORTANT:

The interior colors represent the prosthetic interface: PINK for interfacing with SLIM prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT and SX prosthetic components; GRAY for interfacing with UNIFIT prosthetic components and BLUE for interfacing with MORSEFIT prosthetic components

Dental Implants belonging to this family are compatible with all Prosthetic Components, provided that the prosthetic interface column is matched.

*Each commercial device code is composed of the Implant model + respective 1mm cover screw model (according to the model composition column).

2. DEVICE COMPOSITION

The Dental Implant is manufactured in Titanium Alloy Ti-6Al-4V ELI, according to the specifications of ASTM F3001 “Standard Specification for Additive Manufacturing Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) with Powder Bed Fusion”.

3. DEVICE PRESENTATION FORM

The Dental Implant (models described in Table 1) is supplied individually packaged, in a sterile condition, and presented as follows:

- Primary packaging is a case with a crystal polystyrene lid (GPPs) and white polypropylene (PP) base;
- Secondary packaging is a polyethylene terephthalate (PET) blister heat-sealed with surgical grade paper, containing 01 unit of the medical device, sterilized using Gamma Irradiation.
- External packaging is a rigid cardboard box containing the primary, secondary, and three (03) identification and traceability labels;
- It is the responsibility of the hospital/health care facility to ensure complete traceability of the device, by keeping notes in the patient’s chart of the device code and batch code, or else the identification through one of the three (03) traceability labels accompanying the package;
- If you want to obtain the Instructions for Use in printed format and at no additional cost (including shipping), request by telephone + 55 11 3109-9000, email (customerservice@plenum.bio) or by the website (www.plenum.bio/ifu);
- See below the illustrative images of the primary (Figure 1), secondary (Figure 2), and external (Figure 3) packaging of the Dental Implant in the form it will be commercialized.

FIGURE 1: ILLUSTRATIVE IMAGE OF THE PRIMARY PACKAGING.


FIGURE 2: ILLUSTRATIVE IMAGE OF THE SECONDARY PACKAGING.

FIGURE 3: ILLUSTRATIVE IMAGES OF THE EXTERNAL PACKAGING.

ACCESSORIES:

The Dental Implant has accessory models called Cover Screws, which must be used temporarily to close the platform of the respective dental implant until the prosthetic component or dental prosthesis is installed. The Cover Screw is produced from Grade 4 Titanium according to the specifications of ASTM F67. Table 2 presents illustrative images, codes, models, and main dimensions of the Cover Screw device.

TABLE 2. Illustrative images, codes, models, and dimensions of the Cover Screw.

Prosthetic Interface	SL			
Model composition*	COVER SCREW SL 1mm / COVER SCREW SL			
Illustrative Image				
Codes*	Models	Diameter (mm)	Length (mm)	
018-10-110-01-01	COVER SCREW SL 1mm	4.0	1	
018-10-110-00-01	COVER SCREW SL		0	
Prosthetic Interface	RE			
Model composition*	COVER SCREW RE 1mm / COVER SCREW RE			
Illustrative Image				
Codes*	Models	Diameter (mm)	Length (mm)	
018-10-110-01-02	COVER SCREW RE 1MM	4.1	1	
018-10-110-00-02	COVER SCREW RE		0	

Prosthetic Interface	SH			
Model composition*	COVER SCREW SH 1mm / COVER SCREW SH			
Illustrative Image				
Codes*	Models	Diameter (mm)	Length (mm)	
018-10-110-01-03	COVER SCREW SH 1MM	4.5	1	
018-10-110-00-03	COVER SCREW SH		0	

Prosthetic Interface	UF			
Model composition*	COVER SCREW UF 1mm / COVER SCREW UF			
Illustrative Image				
Codes*	Models	Diameter (mm)	Length (mm)	
018-10-110-01-04	COVER SCREW UF 1MM	4.1	1.0	
018-10-110-00-04	COVER SCREW UF		0.0	

Prosthetic Interface	MF			
Model composition*	COVER SCREW MF 1mm / COVER SCREW MF			
Illustrative Image				
Codes*	Models	Diameter (mm)	Length (mm)	
018-10-110-01-05	COVER SCREW MF 1MM	4.5	1.0	
018-10-110-00-05	COVER SCREW MF		0.0	

IMPORTANT:

The colors represent the Cover screw / Implant interface: PINK for interfacing with SLIM prosthetic components; YELLOW for interfacing with REGULAR prosthetic components; GREEN for interfacing with SHORT and SX prosthetic components; GRAY for interfacing with UNIFIT prosthetic components and BLUE for interfacing with MORSEFIT prosthetic components

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Dental Implant is indicated for dental replacement, to restore the masticatory function and aesthetics of the patient in cases of total, partial, and unitary restorations, including in grafted areas, with the use of a conventional technique from one or two surgical stages, with immediate or conventionally (delayed) loading. The morphological characteristics of the device, such as cylindrical-conical external shape, trapezoidal thread profile, and macro and micro structured surface topography, associated with the correct indication and installation technique, provide immediate mechanical fixation (implant/bone tissue interface), which is subsequently complemented by bone cell adhesion, allowing effective osseointegration and consequently adequate support for exclusive placement in prosthetic components over the dental implant.

Note: Plenum Dental Implants should only be used with the Components manufactured by Plenum (devices not part of these instructions for use and must be purchased separately, being objects of specific ANVISA registers). Prosthetic components of other brands may cause dimensional incompatibility, improper fitting, micro-movements, and increase the risk of corrosion (electrolytic or galvanic) and fracture of these devices, compromising the result of the surgical and prosthetic procedure.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The performance of the Dental Implant is due to its ability to replace the missing dental element. After being screwed into a bone cavity, either milled or anatomical (tooth socket), in the maxilla and/or mandible of the patient, it can anchor different prosthetic components (not part of these instructions for use and must be purchased separately) and a dental prosthesis that can resist masticatory efforts, thus restoring the function and aesthetics of the patient, in cases of total, partial and single restorations, eliminating the need for denture and/or mucosupported removable prostheses (dentures, removable partial dentures, and fixed partial dentures).

6. INSTRUCTIONS FOR USE

This device is restricted to dental use and should only be used by dentists previously trained and qualified for this procedure, with experience and knowledge of implant dentistry techniques, as well as knowledge of anatomy and anesthetic techniques;

- Before the surgical procedure, the dental surgeon must carry out rigorous preoperative planning:

- Detailed anamnesis - to diagnose possible systemic and psychological health conditions that may interfere with the process of bone healing, soft tissues, and osseointegration or that contraindicate the use of dental implants;

- Clinical examinations - to evaluate the visible amount of ridge in the edentulous area (no teeth), quality of the gingival tissue, parafunctional habits, mesiodistal clinical distance, interocclusal space, and local factors that may interfere with the bone and soft tissue healing process and implant osseointegration;

- Laboratories - it is essential to request tests such as: blood count, coagulogram, calcium dosage, fasting blood glucose, creatinine and/or any other test that is necessary for the safety of the patient during the surgical procedure;

- Imaging - radiographs and/or computed tomography to assist in the three-dimensional planning of the implant installation, in the selection of the most appropriate implant, and for the identification of anatomical structures (e.g., inferior alveolar canal, mental foramen, sinus floor);

- The Dental Implant family is composed of implants with various models, and it is up to the dental surgeon to choose the most suitable size for the anatomical and dimensional characteristics that are most appropriate for the individual needs of each patient;

- In case of implant installation in proximity to the sinus floor, the professional must critically evaluate the amount and density of the bone in the area to select the model and the most appropriate installation torque for the implantation area;

- The clinical/surgical planning and surgical technique to be used must be selected by the dental surgeon, based on the available scientific literature, on their clinical and surgical experience and their technical predicates;

- Before starting the surgery, the dental surgeon must submit the patient to a rigorous oral antisepsis to prevent the implant from coming into contact with contaminating substances;

HANDLING: in the procedure of handling and implantation of Dental Implants, the dental surgeon must use only appropriate and sterile surgical instruments;


THE SEQUENCE OF CASE OPENING AND IMPLANT HANDLING.

Primary packaging containing the implant and the implant cover.


Turn the cover counterclockwise as indicated in the OPEN direction on the cover.

Fit the driver connected to the contra-angle to the internal connection of the implant.

Using your thumb and index finger, press the flaps to release the grip that holds the implant in place.

Fit the driver to the cover screw located at the base of the case for safe removal and without contaminating the product.

ATTENTION: The figures are merely illustrative and do not represent the actual dimensions and characteristics of the device.

- To reduce the potential for bone fractures and/or thermal necrosis, it is recommended that all perforations and/or cavity preparation for the installation of the Dental Implant should be done with the Instrumental Kit drills in a range of 800 - 1500 rpm at low rotation and with irrigation constant of sterile physiological solution (these surgical instruments are considered auxiliary instruments, they have separate registrations with ANVISA. Therefore, they must be purchased separately) according to surgical protocols established in Tables 3 and 4.

TABLE 3. Surgical Protocol for the use of Plenum Cutting Instruments according to bone density.

Region	ANTERIOR JAW (BONE TYPE I)
Implant Diameter	Surgical Protocol
Ø3.0	Pilot - 2.15 - 2.85 - TAP (3.0)
Ø3.5	Pilot - 2.15 - 2.85 - 3.35 TAP (3.5)
Ø4.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 TAP (4.0)
Ø4.5	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 TAP (4.5)
Ø5.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 - 4.85 TAP (5.0)
Ø6.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 - 4.85 - 5.35 - 5.85 TAP (6.0)

Region	POSTERIOR JAW (BONE TYPE II)
Implant Diameter	Surgical Protocol
Ø3.0	Pilot - 2.15 - 2.85
Ø3.5	Pilot - 2.15 - 2.85 - 3.35
Ø4.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85
Ø4.5	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35
Ø5.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 - 4.85
Ø6.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 - 4.85 - 5.35 - 5.85

Region	MAXILLA (BONE TYPE III E IV)
Implant Diameter	Surgical Protocol
Ø3.0	Pilot - 1.65*
Ø3.5	Pilot - 2.15
Ø4.0	Pilot - 2.15 - 2.85
Ø4.5	Pilot - 2.15 - 2.85 - 3.35 - 3.85
Ø5.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35
Ø6.0	Pilot - 2.15 - 2.85 - 3.35 - 3.85 - 4.35 - 4.85

*The cutter with diameter Ø1.65 is indicated that it is used, with the appropriate restrictions, only in surgical procedures for installation of the Plenum implant of 3.0 mm diameter in the upper lateral incisor region.

Note: TAP (screw male) is optional in the surgical protocol sequence.


The complete sequence of Plenum drills for a surgical protocol for placement of Dental Implants - RE Ø6.0 x 13 mm for type I bone.

The complete sequence of Plenum drills for a surgical protocol for placement of Dental Implant - SL Ø3.0 x 13 mm in the upper lateral incisor region. Only use the Ø1.65 mm drill for this indication of the surgical protocol.

INSTALLATION

ATTENTION: The Figures are merely illustrative, and they do not represent the actual dimensions and characteristics of the device.

- Rotation for implant insertion from 25 to 50 rpm;
- Maximum torque limit for insertion 60 Ncm, the torque can be measured with the manual ratchet;
- For 3.0 implants, do not exceed 45 Ncm of torque insertion;
- For implant insertion, a manual ratchet or counter-angle may be used; the choice of insertion instrument will depend on the clinical situation and accessibility.

Manual insertion: for manual insertion of the implant, connect the Ratchet to the Plenum Driver (surgical instruments are considered auxiliary instruments and have separate records at ANVISA. Therefore, they must be purchased separately) and place the implant up to the end of the drilling.

Insertion with contra-angle: connect the hexagonal driver to the contra-angle, and the driver to the implant, with a perfect fit between the connection and the implant. On the surgical motor, use a maximum torque of 35 Ncm and a rotation speed from 25 to 50 rpm; complete the implant insertion with the ratchet with torque reference and ratchet connections, short or long;

- Do not exceed 10 Ncm of torque to install the Cover Screw.

- For immediate loading, the Dental Implant must withstand a torque of at least 35 Ncm at the end of the installation.

INSTALLATION OF IMPLANTS USING THE GUIDED SURGERY TECHNIQUE

TABLE 4 - Surgical Protocol for the use of Plenum Guided Instruments.

Implant Length	5mm
Implant Diameter	Surgical Protocol
Ø3.0	X
Ø3.5	X
Ø4.0	X
Ø5.0	Ø2.0x5 - Ø3.5 x5 - Ø5.0 x5 - TAP 5.0**
Ø6.0	Ø2.0x5 - Ø3.5 x5 - Ø5.0 x5 - Ø6.0x5 - TAP 6.0**

Implant Length	8mm
Implant Diameter	Surgical Protocol
Ø3.0	X
Ø3.5	Ø2.0x5 - Ø3.5 x8 - TAP 3.5**
Ø4.0	Ø2.0x5 - Ø3.5 x8 - Ø4.0x8 - TAP 4.0**
Ø5.0	Ø2.0x5 - Ø3.5 x8 - Ø4.0x8 - Ø5.0x8 - TAP 5.0**
Ø6.0	X

Implant Length	10mm
Implant Diameter	Surgical Protocol
Ø3.0	Ø2.0x5 - Ø3.0x10 - TAP 3.0**
Ø3.5	Ø2.0x5 - Ø3.5 x10 - TAP 3.5**
Ø4.0	Ø2.0x5 - Ø3.5 x10 - Ø4.0x10 - TAP 4.0**
Ø5.0	Ø2.0x5 - Ø3.5 x8 - Ø4.0x8 - Ø5.0 x10 - TAP 5.0**
Ø6.0	X

Implant Length	11.5mm*
Implant Diameter	Surgical Protocol
Ø3.0	Ø2.0x5 - Ø3.0x13 - TAP 3.0**
Ø3.5	Ø2.0x5 - Ø3.5 x13 - TAP 3.5**
Ø4.0	Ø2.0x5 - Ø3.5 x10 - Ø4.0x13 - TAP 4.0**
Ø5.0	Ø2.0x5 - Ø3.5 x8 - Ø4.0x8 - Ø5.0 x13 - TAP 5.0**
Ø6.0	X

Implant Length	13mm
Implant Diameter	Surgical Protocol
Ø3.0	Ø2.0x5 - Ø3.0x13 - TAP 3.0**
Ø3.5	Ø2.0x5 - Ø3.5 x13 - TAP 3.5**
Ø4.0	Ø2.0x5 - Ø3.5 x10 - Ø4.0x13 - TAP 4.0**
Ø5.0	Ø2.0x5 - Ø3.5 x8 - Ø4.0x8 - Ø5.0 x13 - TAP 5.0**
Ø6.0	X


*For 11.5mm length implants, use the 13mm drill and drill to the 11.5mm indication on your guide.

**TAP and punch are optional following the surgical protocol, and the TAP is recommended only for type I bone.

The complete sequence of drills for a guided surgical protocol for placement of Dental Implants - RE Ø5.0 × 13mm for type I bone.

INSTALLATION


ATTENTION: The figures are merely illustrative and do not represent the actual dimensions and characteristics of the device.

- Rotation for implant insertion from 25 to 50 rpm;
- Maximum torque limit for insertion 60 Ncm; the torque can be measured with the manual ratchet;
- For 3.0 implants, do not exceed 45 Ncm of torque insertion;
- For implant insertion, a manual ratchet or counter-angle may be used; the choice of insertion instrument will depend on the clinical situation and accessibility.

INSERTION: The implant insertion can be done manually or with a contra-angle. For manual insertion of the implant, connect the Ratchet to the Plenum Driver (surgical instruments are considered auxiliary instruments, they are registered separately at ANVISA. Therefore, they must be purchased separately) and place the implant until the end of the drilling, and for insertion with the contra-angle connect the hexagonal driver to the contra-angle, the driver to the implant, with perfect fitting between the connection and the implant. In the surgical motor, use a maximum torque of 35 Ncm and rotation from 25 to 50 rpm;

- Do not exceed 10 Ncm of torque for the installation of the Cover Screw;
- For immediate loading, the Dental Implant must withstand a torque of at least 35 Ncm at the end of the installation.


In cases of high resistance to implant insertion:

- Excessive torque during implant installation should be avoided (may cause deformation damage to the implant or internal connection and also cause excessive pressure to the bone tissue).
- Torque exceeding 60 Ncm (RE and SH implants) or 45 Ncm (for SL implants), during implant insertion at any point, stop and rotate the implant counterclockwise one and a half times to allow removal of the implant and return to the initial implant insertion procedure.
- If you still find resistance in the implant insertion, remove the implant and place it in the original packaging, attaching it again to the grips. Extend the bone drilling of the surgical socket preparation by following the correct surgical protocols described in tables 3 and 4 or use the TAP for the implant diameter. Do not exceed the marking of the maximum limit indicated in the marking of the respective TAP.

ATTENTION: drillings performed outside the conditions mentioned above and the strict non-compliance with the surgical procedures described in tables 3 and 4 may exceed the values established and indicate for insertion torque and consequently may directly compromise the osseointegration of the implant;

- The surgeon must always consider that from a biomechanical point of view, the classic technique of placing a dental implant recommends that, whenever possible, the implant installation must respect a correct three-month positioning and must be installed at least 1mm below, respecting the minimum medium-distal tolerance of 1.5 mm between adjacent teeth and implants and 3 mm between implant and implant;

- Carefully observe the correct dimensional combination between the selected implant and its respective prosthetic components so that the masticatory efforts are correctly distributed throughout the implant-component set, thus maintaining the mechanical properties guaranteed by the manufacturer (Plenum components are not part of these instructions for use, and must be purchased separately);

- Special care must be taken between the existing connection on the inside of the Dental Implant and the accessories/ components that fit temporarily or permanently (surgical instrument provided separately), as an unsuitable fit between the two elements or the excessive use of force can cause damage to the implant connection, compromising the integrity of the device and the proper fixation of the components;

- The dental surgeon should carefully evaluate inserting and removing a dental implant more than once in the same patient. This can cause deterioration of the threads, the self-cutting tip, and the implant connection. This operation should only be performed in exceptional cases if the insertion torque is too high. In this case, it is recommended to remove the implant, followed by adjustment and relief of the bone site using a new drill or sequence of drills;

- The dental surgeon should rigorously evaluate the initial mechanical stability of the implant. If optimal and safe mechanical stability has not been achieved, the dental surgeon should assess and consider possible removal and/or replacement;

- Postoperative care and proper hygiene/cleanliness of the operated site are extremely important. The patient must be instructed repeatedly about the hygiene techniques and the existing risks until there is complete osseointegration of the implant;

- The dental surgeon must guide the patient, who must understand the importance of follow-up in the postoperative period. In the postoperative period, the stability, integrity, and positioning of the dental implant used must be checked, as well as the osseointegration and biological completion of the implant at the appropriate time;

- To protect the patient and the dental surgeon, it is important to ensure the correct and complete traceability of the implant used by annotating the patient's protocol with the device code and batch code, or by fixing in the patient's protocol the identification and traceability labels that accompany the device.

7. PRECAUTIONS

- The success of a dental implant placement procedure is related to the correct indication and selection of the implant, the proper surgical technique, the stability and durability of the implants and components, as well as the timing and quality of osseointegration, which can be affected by numerous biological, clinical and biomechanical factors, which are extrinsic to the quality of the implants and components. These factors may limit the useful life of the devices and consequently cause their failure and early loss. Therefore, faithful obedience to all indications for use, precautions, restrictions, warnings, contraindications, and other observations contained in these instructions for use are essential for the success of the procedure;

- The dental surgeon should not use this device before carefully reading all the instructions for use and the information on the packaging. Failure to observe these points may lead to complications for the dental surgeon and especially for the patient;

- The device must be kept in the original packaging until the use time and handled only in sterile surgical fields. If one of the external packaging that ensures sterility is damaged or open, the device should be discarded;

- The device label presents the expiration date, device code, batch code, the ANVISA device registration number, the technical product manager and company information, as well as the respective symbology, allowing full traceability of the implant;

- Before using the Dental Implant, the dental surgeon must perform a detailed analysis of the surgical instruments (Auxiliary Instruments, provided separately) used to place this device. These instruments must be intact and complete. Particular attention should be paid to the drivers used in inserting the implant, the drills' cutting profile, and the ratchet's integrity. Damage to the tip of the driver may prevent the proper placement of the implant, damage to the drills may cause inadequate bone drilling for the placement of the implant, and damaged ratchets may cause an excessive or insufficient force in the placement of the implant, compromising the result of the treatment. Therefore, do not use this device if there is any damaged auxiliary instrument for specific use;

- The patient should be warned that failure to comply with the postoperative period instructions may lead to failure and early loss of the implant, requiring new surgery for review or removal;

- For the safety and effectiveness of the surgery, only Prosthetic Components and Auxiliary Surgical Instruments manufactured by Plenum should be used. Components and surgical instruments from other manufacturers can cause dimensional incompatibility between parts, inadequate fit, micro-movements, increased risk of corrosion (electrolytic or galvanic), fracture of the implants, and early failure of the implants and components, compromising the treatment result;

- Incorrect selection of the drill or implant may cause temporary or permanent paresthesia of the nerves of the mandibular region and drilling of the sinus membrane;

- The wrong choice of implant and/or components to be used, as well as errors in the indication, manipulation, preparation of bone tissue, and insertion of implants, whose procedures can cause damage to the devices, contributing to failure by fatigue, fracture, or loosening, and lack of osseointegration of the implant;

- Use of excessive torque during insertion of an implant may lead to the device failure;

- The patient must be informed of the risks involved in a dental implant that will be overloaded beyond its functional capacity, a fact that may cause bone loss, failure of the implant and its components, and other elements of prosthetic restoration;


- Do not force the insertion of an implant beyond the depth established by the milling, as it may cause deformation in the hexagonal connection of the implant and/or the insertion driver, in addition to fracturing the implant;

- If there is an impact, deformation, scratches, cracks, or dents visible to the implant, which may impair the correct and proper functioning of the device, it must be discarded;

- The main and most probable causes of the lack of bone integration with the consequent loss of an implant are: infections, localized or systemic diseases, low bone quality, insufficient bone quantity, inadequate osteotomy, inadequate oral hygiene, occlusal trauma, inadequate occlusal conditions, mechanical overload of the implant, poorly sized implant, initial mechanical instability of the implant, failure in the bone tissue freshening process, bone necrosis, failure in the irrigation of bone drilling, surgical instructions belonging to other systems, wear of the surgical instrument, lack of previous professional training, contamination of the implant, technical error in the placement of the components and patient who does not comply with the postoperative instructions. The dental surgeon should consider all these aspects in preoperative planning;

- Incorrect handling of implants and/or components of small size inside the patient's mouth can result in the risk of aspiration and/or swallowing of these elements by the patient;

- Because the Ø1.65mm drill has a reduced diameter and consequently lower mechanical resistance to bending, compared to other drills in the Plenum portfolio, we advise caution not to exert tilting movements, as illustrated in the figures below, to avoid fracture of the drill and compromise the surgical procedure;

- Smoking patients, tobacco users or drugs, alcoholics, malnourished, or with any other relevant systemic condition should be alerted to the increase in the incidence of non-osseointegration of the implant;

- Patients with an atrophic mandible or maxilla, inadequate oral hygiene, clinical pictures with hematological disorders, hormonal diseases, and other uncontrolled systemic alterations have restrictions to this device;

- For elderly patients with mental problems, chemical impairments, or with difficulty following postoperative recommendations, this device may present a greater risk of failure;

- The safety and compatibility of Plenum Dental Implants in patients undergoing magnetic resonance were not evaluated. Heating has not been tested, or distortion suffered from dental implants. Doing an MRI on a patient with this device may result in harm to the patient. It is always important to inform the medical team about any metallic medical device present in the body so that they can safely decide which situations are at risk and which are safe for the MRI examination;

- The patient must be informed of all limitations, risks, complications, warnings, and possible adverse effects arising from the procedure for installing a dental implant;
- The patient must be informed of the risks involved in the application of excessive mechanical load or trauma that affects the operated region, which may imply failure, loss of mechanical stability, loosening, breakage, fracture, or loss of the dental implant and its components;
- It is important and fundamental that the dental surgeon expresses all the above topics to his patient and that the same has a complete and clear understanding of all these aspects.

8. WARNINGS

- Single-use only, do not resterilize or reuse.
- Do not use if package is opened or damaged or if expiration date has been exceeded.

9. CONTRAINDICATIONS

The dental surgeon shall prudently evaluate the patient's condition before surgery and use of a Dental Implant and check for any contraindications reported below. Contraindications may include, but are not limited to:

- Patient who has insufficient quantity and bone quality to provide primary mechanical stability to the implant;
- Patient with atrophic alveolar ridge with insufficient quantity and remaining bone quality to provide height and thickness of bone tissue around the implant;
- Patient with poor bone quality due to osteoporosis or osteopenia, which can cause early dental implant failure due to lack of osseointegration;
- Inflammation and/or periodontal infection;
- Recent history of systemic or localized infection;
- Fever;
- Leukocytosis;
- Patient with allergy or sensitivity to the raw material used in the manufacture of Plenum Dental Implant (see "Device Composition");
- Patient in a growth phase, with incomplete skeletal system development;
- Patient with developmental disorders of the skeletal system;
- Parafunctional habits;
- Patient with insufficient intermandibular space;
- Non-treatable occlusion and joint problems;
- Pregnancy and lactation period;
- Use in patients with immunological changes and systemic diseases that result in poor healing or the absence of integrity of existing tissue at the implant site;
- Patient with vascular disorders or coagulation disorders;
- Alcoholic patient or drug user who cannot correctly follow postoperative instructions;
- Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;

- Patient unable to maintain good oral hygiene;
- Immunosuppressive pathologies;
- Patients undergoing radiotherapy or use of medications that compromise the vascularization of bone tissue;
- All contraindications included by the dental-scientific literature known for using dental implants and their components.

10. POSSIBLE ADVERSE EFFECTS


Positive results may not be achieved during all surgical procedures. Every surgical procedure presents risks and possibilities of complications, such as infections, hemorrhages, allergic drug reactions, and anesthetic risks. Possible adverse effects associated with the use of Dental Implant may include, but are not limited to:

- The surgical procedure to install a dental implant can bring complications in the trans- and post-operative period, such as: pain, edema, local bleeding, abscess, fistula, dehiscence, difficulty speaking, nerve injury or compression (sensory or motor paresthesia, partial or total, transient or permanent), soft tissue injury, local or systemic infection, and other complications typical of a dental implant placement procedure;
- Pain or discomfort due to the presence of the dental implant or due to the surgical procedure;
- Sinusitis due to placement of the dental implant in maxillary sinuses;
- Accidental drilling of the cortical bone of the facial sinuses, sinus membrane, bone plate, inferior alveolar canal, or soft tissues;
- Allergic reaction to the raw material of the implant;
- Intense fibrotic tissue reaction around the surgical area;
- Alteration of bone structure due to the imbalance of normal forces/tensions on the bone/implant interface;
- Bone necrosis, insufficient vascularization, bone reabsorption, and insufficient bone formation, which can cause the loss of fixation of the dental implant and consequently the component; these effects are correlated to the patient's systemic facts or surgical procedure that is inappropriate without complying with the instructions for use described for this family of components;
- Premature suture dehiscence, the tension in the coaptation of the surgical edges, or premature exposure of the recipient bed can alter the healing processes and result in implant failure;
- The patient must be informed and understand that complications in all surgical procedures are possible. The dental surgeon must inform the patient of all these complications and risks resulting from implant surgery.

11. SYMBOL DESCRIPTIONS

Table 5 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 5 - Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
Proibido reprocessar	Reprocessing prohibited (Brazilian Requirement)

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

12. AUXILIARY INSTRUMENTS

- All surgical instruments developed exclusively for the installation of the Dental Implant are considered auxiliary instruments;
- Plenum Instruments are manufactured in stainless steel, whose raw material provides devices with high strength and durability, according to requirements specified by NBR 13911; ASTM F899, and ISO 5832-1;
- It is essential to use instruments manufactured by Plenum to ensure safety and effectiveness in the installation of implants according to the surgical protocol recommended by the manufacturer;
- Adequate irrigation, rotation and torque are recommended during the installation of the Plenum Dental Implant; the sequence of drills recommended by the manufacturer is used. Do not exceed the technical specifications of the drills indicated by the manufacturer;
- The Plenum Dental Implant has a hexagonal connection in its internal region. The dental surgeon must use only Plenum drivers to avoid damaging the implant connection and also the insertion drivers;
- For properly installing the implants, use only the surgical instruments manufactured by Plenum. The list of surgical instruments can be found in the respective printed or electronic communication materials or on the Plenum website (www.plenum.bio). The use of Plenum instruments must follow the instructions of the manufacturer;
- Surgical instruments are not the primary purpose of these instructions for use; they have separate registrations with ANVISA. Therefore, they must be purchased separately.
- Use only instruments manufactured or validated by Plenum, so there is no incongruence between the fittings and irreversible damage to the implants and instruments.

13. STERILIZATION

- The Dental Implant is provided in a sterile condition, sterilized using Gamma Irradiation;
- The Dental Implant is a single-use device and cannot be reprocessed, resterilized, or reused;
- Before use, check the integrity of the packaging and the expiration date. Do not use the device if the package is violated or with the sterilization expiration date expired;
- Date of manufacture, sterilization expiration date, device code, batch code, and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments;
- The packaging opening for the device's use must only be done by qualified and trained personnel to perform this procedure.

14. EXPIRATION DATE

- Indicated on device labeling.
- Before use, check the expiration date, do not use the device if it has expired.

15. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Dental Implant and comply with health surveillance requirements, the surgeon or other staff must keep the information regarding the device used in the patient's medical records. Additionally, this information must also be passed on to the device distributor and patient to complete the traceability cycle. The information needed for traceability is related to the device used, surgical procedure, and patient, as follows:

- Patient's name;
- Surgeon's name;
- hospital/clinic/office identification;
- manufacturer and distributor identification;
- Date of surgery;
- Device code used;
- Batch code;
- Quantity used;
- ANVISA device registration number.

The responsible dental surgeon and their team should use the information contained in the three (03) traceability labels supplied with the Dental Implant, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so they have information about the implant used in their procedure.

The labels contain the following information necessary for the traceability of the Dental Implant used.

- Manufacturer identification;
- Device name or commercial model;
- Device code;
- Batch code;
- ANVISA device registration number .

Traceability information is necessary for notification by the health service and/ or the patient to the National Health Surveillance Agency - ANVISA and to the manufacturer when complaints, technical complaints, and/or adverse events occur for the conduct of the appropriate investments.

16. STORAGE, TRANSPORT, AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the device if the packaging is damaged or if the expiration date has expired. If the packaging is damaged, the device must be considered non-sterile and discarded;
- The device must be stored and transported in a clean environment away from sources of heat or humidity, sunlight, and the action of the weather, to prevent any damage or change in the packaging or its physical form that makes it impossible to use;
- The device must be stored on shelving or drawers/cabinets to allow it to be easily viewed and handled by the operator. Do not store the device near lamps so that the information contained in the packaging is not affected;

- The device shall only be opened at the time of use in a sterile surgical field;
- The device handling should only be done by qualified and trained personnel to perform this procedure.

17. DISPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After using the device or removing it from the patient, the Dental Implant that is not subject to analysis or subsequent studies, must be disposed of by the health service establishment (hospitals, offices, and clinics);
- It is recommended that the Dental Implant be disassembled entirely to ensure that it cannot be reused. The health service institution/establishment is solely responsible for disassembling the device;
- After de-featuring, the device must be identified as being unfit for use and disposed of by health service waste management standards (hospitals, offices, and clinics), in accordance with applicable current environmental laws;
- The Dental Implant not used during surgery must also be discarded following the guidelines above.

18. CUSTOMER SERVICES

- If the device presents any adverse event with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is causing any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, batch code, and patient data; directly notifying the Plenum Customer Relations Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and properly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - CEP 13212-213 - Jundiaí - SP - Brazil. Phone: +55 11 3109-9000; E-mail: customerservice@plenum.bio
- ANVISA also provides a channel to receive notifications of incidents, adverse events, and technical complaints related to the use of devices and services under health surveillance, which can be found on the National Health Surveillance Agency website - ANVISA, at www.anvisa.gov.br, link NOTIVISA.
- FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

19. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health

professional and addressed to regulatory authority standards (see the Customer Services section) when it is observed that issues involve a risk to the patient's health.

The warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing, resterilization, or reuse of the device by the user;
- Improper use or use other than in a surgical environment;
- Device not used according to instructions for use.

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