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PLENUM® OSS 3D β
English

PLENUM® OSS 3D β
Español

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PLENUM[®] OSS 3D_β

INSTRUCTIONS FOR USE

TECHNICAL NAME: Bone graft and associated devices

TRADE NAME: Plenum[®] Oss 3D_β

MODEL: For a detailed description of the device, see Table 1

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA DEVICE REGISTRATION N°.: 81684340021.

Device Code, batch code, and date of manufacture and sterilization: See device labeling.

REV.01

SUMMARY




1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	4
3. DEVICE PRESENTATION FORM	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTION FOR USE	5
7. PRECAUTIONS	6
8. WARNINGS	7
9. CONTRAINDICATIONS	7
10. POSSIBLE ADVERSE EFFECTS	8
11. SYMBOLS	8
12. STERILIZATION	8
13. EXPIRATION DATE	9
14. IDENTIFICATION AND TRACEABILITY	9
15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	9
16. DISPOSAL OF DEVICE	9
17. CUSTOMER SERVICE	9
18. LIMITED WARRANTY	9

1. DETAILED DEVICE DESCRIPTION

Plenum® Oss 3D_β is a resorbable bone graft based on bioceramic synthetic (β-Phosphate tricalcium) produced from the additive manufacturing process (3D printing), which allows a standardized

macrostructure with pores of 400 μm, favoring a predictable replacement of the graft with newly formed bone tissue and concomitant slow resorption during the process of tissue repair or regeneration.

TABLE 1: Illustrative images, codes and dimensions of the Plenum® Oss 3D_β commercial models.

Illustrative images	Codes	Commercial models / Sizes	
		Length (L) x Width (W) x Height (H)	Volume
	018-25-B40-00-01	BLOCK - 10 x 10 x 03 mm	0,30 cm ³
	018-25-B40-00-02	BLOCK - 10 x 15 x 03 mm	0,45 cm ³
	018-25-B40-00-03	BLOCK - 10 x 10 x 05 mm	0,50 cm ³
	018-25-B40-00-04	BLOCK - 20 x 20 x 08 mm	3,20 cm ³
	018-25-B40-00-05	BLOCK - 20 x 20 x 10 mm	4,00 cm ³
	018-25-B40-00-06	BLOCK - 20 x 30 x 08 mm	4,80 cm ³
	018-25-B40-00-07	BLOCK - 20 x 30 x 10 mm	6,00 cm ³
	018-25-B40-00-08	WEDGE - 20 x 15 x 6 mm - 9°	1,35 cm ³
	018-25-B40-00-09	WEDGE - 20 x 15 x 8 mm - 14°	1,65 cm ³
	018-25-B40-00-10	WEDGE - 20 x 15 x 10 mm - 19°	1,95 cm ³
	018-25-B40-00-11	WEDGE - 20 x 15 x 12 mm - 24°	2,25 cm ³
	018-25-B40-00-12	WEDGE - 20 x 15 x 14 mm - 29°	2,55 cm ³
		018-25-B40-00-13	CYLINDER - 08 x 10 mm
018-25-B40-00-14		CYLINDER - 10 x 10 mm	0,80 cm ³
018-25-B40-00-15		CYLINDER - 12 x 10 mm	1,10 cm ³
018-25-B40-00-16		CYLINDER - 14 x 10 mm	1,50 cm ³
018-25-B40-00-17		CYLINDER - 08 x 20 mm	1,00 cm ³
018-25-B40-00-18		CYLINDER - 10 x 20 mm	1,60 cm ³
018-25-B40-00-19		CYLINDER - 12 x 20 mm	2,30 cm ³
018-25-B40-00-20		CYLINDER - 14 x 20 mm	3,10 cm ³

ANCILLARY COMPONENTS

Plenum® Oss 3D_β medical devices can be fixed with fixation screws in order to achieve primary stability of the graft to the recipient bed (ancillary components are not part of this product family, they have a separate ANVISA

registration and must be purchased separately. Maximum recommended torque limit for screw fixing is 10 N.cm.

The commercial models of fixing screws in the Plenum® Graft family are described below:

Commercial name	Technical name	Composition	Model / Dimension
Plenum® Graft Anvisa Registration 81684340018	Non-bioabsorbable, non-sterile craniofacial bone screw	The screws in the Plenum® Graft family are made from ELI - Ti6Al4V titanium alloy, in accordance with ASTM F136-13 (2017) and ABNT NBR ISO 5832-3:2017 standards.	CORTICAL SCREW RT Ø1,4X4mm CORTICAL SCREW RT Ø1,4X5mm CORTICAL SCREW RT Ø1,4X6mm CORTICAL SCREW RT Ø1,4X7mm CORTICAL SCREW RT Ø1,4X8mm CORTICAL SCREW RT Ø1,4X9mm CORTICAL SCREW RT Ø1,4X10mm CORTICAL SCREW RT Ø1,4X11mm CORTICAL SCREW RT Ø1,4X12mm CORTICAL SCREW RT Ø1,4X13mm CORTICAL SCREW RT Ø1,4X14mm CORTICAL SCREW RT Ø1,4X15mm CORTICAL SCREW RT Ø1,4X16mm CORTICAL SCREW RT Ø1,4X17mm CORTICAL SCREW RT Ø1,4X18mm CORTICAL SCREW RT Ø1,4X19mm CORTICAL SCREW RT Ø1,8X4mm CORTICAL SCREW RT Ø1,8X5mm CORTICAL SCREW RT Ø1,8X6mm CORTICAL SCREW RT Ø1,8X7mm CORTICAL SCREW RT Ø1,8X8mm CORTICAL SCREW RT Ø1,8X9mm CORTICAL SCREW RT Ø1,8X10mm CORTICAL SCREW RT Ø1,8X11mm CORTICAL SCREW RT Ø1,8X12mm CORTICAL SCREW RT Ø1,8X13mm CORTICAL SCREW RT Ø1,8X14mm CORTICAL SCREW RT Ø1,8X15mm CORTICAL SCREW RT Ø1,8X16mm CORTICAL SCREW RT Ø1,8X17mm CORTICAL SCREW RT Ø1,8X18mm CORTICAL SCREW RT Ø1,8X19mm

2. DEVICE COMPOSITION

The commercial models that are part of the Plenum® Oss 3D_β family are composed of β-tricalcium phosphate (≥ 95% β-TCP).

3. DEVICE PRESENTATION FORM

Plenum® Oss 3D_β is supplied unit-packed in the models described in Table 1, in sterile condition, and presented as follows: O1 (one) Plenum® Oss 3D_β and O3 (three) identification and traceability labels in the following packaging:

- Primary packaging in the form of a case with a transparent lid in crystal polystyrene (GPPS), a white base in polypropylene (PP) and an orange silicone cradle, containing O1 unit of the product in one of the presentations described in Table 1;
- Secondary packaging in the form of a PET blister and heat-sealed with Tyvek surgical grade paper, sterilized by Gamma Irradiation;
- External packaging in the form of a rigid carton containing the primary and secondary packaging and 3 identification and traceability labels;

Below are illustrative images of the primary and secondary packaging (Figure 1) and the outer packaging (Figure 2) of a model of Plenum® Oss 3D_β, in the form in which it will be marketed.


FIGURE 1: Illustrative images of primary and secondary packaging.


FIGURE 2: Illustrative images of external packaging.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Plenum® Oss 3D_β is a synthetic bone graft whose intended use is to fill and/or reconstruct bone cavities or defects in the skeletal system that are not intrinsic to the stability of the bone structure. It is indicated for guided tissue regeneration, augmentation/reconstruction and filling of bone defects such as atrophic alveolar ridge, empty spaces (gaps), bone cavities created surgically or through traumatic injuries and bone deformities.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Plenum® Oss 3D_β is a resorbable synthetic bone graft for use in the reconstruction and/or guided regeneration of bone defects. The porous macro and microstructures and chemical composition confer effective osteoconductivity, high hydrophilic properties, controlled resorption/dissolution of the crystalline phase constituting this bone graft, which is gradually absorbed by the body and replaced by newly formed bone tissue during the process of bone tissue repair or regeneration.

The compressive strength of these bone grafts ensures that they can be used for drilling and screw fixation. The adjacent contact of the bone graft (Plenum® Oss 3D_β) with the recipient bed increases the area of contact between them,

facilitating the revascularization and the migration of osteogenic cells into the bone graft. In general, resorbable synthetic materials for bone regeneration are the materials of choice, as these materials are free of possible contaminants or organic molecules from homogenous and heterogenous grafts, which could induce immune responses. In this way, synthetic materials minimize the risk of infection, inflammation and other post-operative complications.

6. INSTRUCTION FOR USE

Carefully read the instructions for use, recommendations, and guidelines before using the Plenum® Oss 3D_β synthetic bone substitute:

- Verify the device's expiration date before use. Do not use devices if the packaging has been damaged or if the expiration date has expired;
- Plenum® Oss 3D_β is supplied ready for use;
- Plenum® Oss 3D_β does not require prior hydration in saline solution or the patient's own blood to improve the efficiency of the product's handling;
- The surgeon (physician or dentist) is responsible for choosing the device model with the most appropriate characteristics and dimensions after an individual clinical analysis of each case. Always have more product available than you intend to use;
- The use of this bone graft with inappropriate techniques, procedures and surgical conditions could harm the patient, resulting in unsatisfactory reconstruction or tissue regeneration results;
- The techniques employed for the use and application of Plenum® Oss 3D_β depend on the preference of the surgeon (physician or dentist), and they are responsible for choosing the therapeutic approach that is used, model and dimensions of the device, application technique, as well as the criteria for monitoring and evaluating the results of surgery;
- The surgeon must take a careful surgical approach to obtain enough soft tissue for full, tension-free coverage of the Plenum® Oss 3D_β graft by first intention;
- Proper preparation of the site (host bone) at which the Plenum® Oss 3D_β will be applied is of the utmost importance. The device must fill the entire area to be treated, and any bleeding must be controlled so that the bone substitute adapts perfectly to the host bone and provides the expected results;
- The surgeon shall use only appropriate and sterile surgical instruments during procedures for manipulation and implantation of Plenum® Oss 3D_β;
- The surgeon must consider their clinical and surgical experience to choose the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients;
- Association with other regenerative products, i.e. complementary biomaterials (ancillary components) is the surgeon's criterion of choice and responsibility (these products are not part of this product family and must be purchased separately).

Handling and manipulation of Plenum® Oss 3D_β:

1. Remove the blister containing the product from the outer carton.

2. Open the blister on a sterile surgical drape and place the primary packaging (case) plus product on the operating table.

3. Open the sterile glove case by turning the lid anti-clockwise.


4. Once the case has been opened, remove the block using tweezers or the thumb and forefinger.

CUSTOMIZATION OF PLENUM® OSS 3D_β AND ADAPTATION TO THE RECEIVING BED:

- Plenum® Oss 3D_β are medical devices that can be adapted and adjusted, but the limits of mechanical conformation must be respected as shown in the table below:

Minimum thickness for screw fixation	3 mm	
Minimum graft edge thickness for screw fixation	3 mm	
Maximum number of fixing screws per commercial model	<p>Minimum: 1 screw (for fixings in the center of the bone graft);</p> <p>Maximum: 2 screws (for fixation at the ends of the bone graft); Note: The fixation screws must be positioned in such a way that they respect the minimum distance of 3 mm between them and between the ends of the graft (edge of the graft).</p>	

- Customize the Plenum® Oss 3D_β according to the anatomy of the recipient bed using sterile straight-piece drills and surgical/prosthetic discs (examples: 702 drills, Maxicut line models and diamond discs);
- Finish properly, leaving no sharp edges or corners to avoid perforating the flap;
- Adapt the graft to the recipient bed without leaving steps, leveling all sides of the graft to the bone defect.
- Maximum recommended torque limit for screw fastening: 10 N.cm

ATTENTION: The pictures are for illustrative purposes only. They do not represent the actual dimensions and characteristics of the product.

7. PRECAUTIONS

The surgeon should not use Plenum® Oss 3D_β before carefully reading all the instructions for use and the information on the packaging. All the information on use, warnings, precautions and recommendations mentioned in the instructions for use must be taken into account. Failure to observe these points could lead to possible adverse events for the patient;

- The handling of Plenum® Oss 3D_β should only be carried out by surgeons (doctors or dentists) who have been trained and qualified for this procedure, with experience and knowledge of facial bone reconstruction techniques, guided tissue regeneration and the use of bone substitutes;
- The product label shows the expiry date, product code and batch number, ANVISA registration number, technical manager and company details, as well as the respective symbology, allowing the product to be fully traceable;
- The Plenum® Oss 3D_β is available in various dimensions (models available in Table 01), and it is up to the surgeon (doctor or dentist) to choose the product with the most appropriate characteristics and dimensions, after a complete individual clinical analysis of each case;
- The wrong choice of bone graft to be used, as well as errors in the indication, manipulation, preparation of bone tissue or soft tissue, and installation of the product, can cause damage to the physical structure and contamination of the product and contribute to its failure;
- The decision on the graft model chosen for the reconstruction of the bone defect depends directly on the surgeon's preference in relation to the size of the bone defect, handling, usability and surgical technique to be employed, where the aim is to maximize contact between the bone tissue and the medical device to enhance graft integration and bone healing;
- Before starting the surgery, the surgeon must subject the patient to rigorous oral antisepsis to prevent the product from coming into contact with contaminating substances;
- Any product that is contaminated or damaged during handling and preparation must be disposed of in accordance with the guidelines described in the Product Disposal section;
- Any leftovers of Plenum® Oss 3D_β must not be used. Any product left over after the end of surgery must be disposed of in accordance with health service waste management regulations (hospitals, practices and clinics) and applicable environmental legislation;
- It is recommended to install dental implants and/or load the reconstruction site with Plenum® Oss 3D_β 8 months after surgery;
- Patients must be warned that a failure to comply with the post-operative instructions may lead to failure of the surgical procedure, requiring further surgery for revision or removal of the bone graft;
- Patients with inadequate oral hygiene, clinical conditions with hematological disorders and other uncontrolled systemic alterations can reduce the effectiveness of the regenerative process and lead to failure of the bone graft;
- Before the surgical procedure, the surgeon (doctor or dentist) must carry out rigorous preoperative planning, by means of clinical and imaging examinations (X-rays or CT scans). They should also consider the need for a detailed anamnesis and complementary tests on the patient's general health (blood count, coagulogram, calcium dosage, blood glucose, etc.), since this information directly influences the body's biological response to the use of the bone graft, Plenum® Oss 3D_β, and then define a surgical plan in advance, selecting the most suitable model for the tissue reconstruction/regeneration procedure;
- It is recommended that preoperative examinations be requested according

to the patient's anamnesis;

- The main and most likely causes of the lack of incorporation and resorption of Plenum® Oss 3D_β with the consequent loss of the graft are: infections, localized or systemic diseases, use of membranes, meshes or screens inappropriately, insufficient blood supply or nutrients, inadequate hygiene, trauma or overload at the surgical site, contamination of the product in the pre- and/or trans-operative period or in the post-operative period when there is exposure of the bone graft in the cranio-maxillofacial region and oral cavity, mechanical instability due to lack of adequate graft fixation, bone necrosis, lack of professional training, technical error in the installation of the product, and patients who do not follow the post-operative instructions. The surgeon must consider all these aspects in his preoperative planning;
- Plenum® Oss 3D_β is just an auxiliary material for bone grafting. The success of the therapy depends on the intervention technique and the patient's systemic response;
- The success of a surgical procedure for guided tissue regeneration or bone reconstruction using Plenum® Oss 3D_β associated with other regenerative products, bone grafts (autogenous, homogenous or heterogenous) or membranes is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/recipient bed incorporation.
- The patient must be informed of all the indications for use, contraindications, precautions, warnings, limitations, risks, complications and possible adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques using synthetic bone substitutes;
- The patient must be informed of the risks involved in applying excessive mechanical load or trauma affecting the operated region, which may result in failure, loss of mechanical stability, exposure and displacement of the bone graft;
- Post-operative care and correct hygiene/cleaning of the surgical site are extremely important. The patient should be instructed about the risks until the bone and/or gum tissue at the surgical site has completely repaired;
- The patient must be properly instructed on post-operative care. The patient's ability and willingness to follow the surgeon's instructions is one of the most important aspects of a surgical procedure;
- For elderly patients, those with motor problems and limitations, those with mental problems, drug addicts, or those with any other difficulty in following post-operative recommendations, this product may represent a greater risk of failure, as these patients may ignore the instructions and restrictions arising from this procedure;
- Plenum® Oss 3D_β has not been tested on breastfeeding women, so it is not known whether the metabolites of Plenum® Oss 3D_β pass into breast milk;
- The surgeon must clearly guide the patient, who must understand the importance of post-operative follow-up, in which the stability and integrity of the product must be verified, as well as incorporation and/or resorption at the appropriate time;
- It is important and fundamental that the surgeon explains all of the above to the patient, verbally and in writing, and that the patient has a complete and clear understanding of all these aspects;
- In order to protect the patient and their surgeon, it is important to ensure the correct and complete traceability of the Plenum® Oss 3D_β used, by noting the product code and batch in the patient's protocol, or by attaching the identification and traceability labels that come with the product to the patient's protocol.

8. WARNINGS

- Plenum® Oss 3D_β is supplied in sterile condition, sterilized by Gamma Irradiation. Do not use the product if its packaging has been tampered with or if its expiry date has passed;
- This is a SINGLE USE product, so reprocessing, resterilization and reuse are prohibited and may pose risks of failure and infection to the patient. If the product is exposed to one of these situations, it will result in contamination and permanent damage to the product and/or cross-infection of the patient;
- The product should be kept in its original packaging until use and should only be handled in a sterile surgical field;
- A patient who is an alcoholic or drug user may not follow the post-operative instructions correctly, resulting in failure of the treatment with the medical device;
- Using the medical device on a patient with a mental illness or who is unable to follow the recommendations and care during the post-operative phase can lead to treatment failure;
- Using the medical device on a patient who is unable to maintain good hygiene at the surgical site (motor or physical limitations) can lead to treatment failure;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques.

9. CONTRAINDICATIONS

The surgeon must prudently assess the patient's condition before surgery and the use of Plenum® Oss 3D_β, and check for any contraindications. Contraindications may include, but are not limited to:

- Patients with insufficient quantity and quality of gingival or soft tissue to provide adequate flap coverage without tension to avoid accidental exposure of the synthetic bone graft and render the procedure unfeasible;
- Patients with poor bone quality or patients with type IV bone due to osteoporosis or osteopenia, which can lead to early failure of the bone graft due to inefficient tissue repair or regeneration;
- Use in sites that suffer stress, direct mechanical strain, or parts of the graft that are exposed to constant deformation and/or tension forces, such as the use of mobile prostheses on the bone graft before the recommended period of integration of the bone graft into the recipient bed of 08 months post-operatively;
- Inflammation and/or infection near or at the surgical site;
- Recent history of systemic or localized infection;
- Fever;
- Leukocytosis;
- Patient with allergy or sensitivity to tricalcium β-phosphate;
- Patient with developmental disorders of the skeletal system;
- Pregnancy and lactation;
- Use in patients with immunological alterations and systemic diseases that result in poor healing or a lack of tissue integrity at the surgical site;
- Use in patients with vascular disorders or coagulation disorders;

- Calcium metabolism disorder;
- Treatment with steroids and other drugs that interfere with calcium metabolism;
- Immunosuppressive pathologies;
- Patients undergoing radiotherapy, chemotherapy or the use of drugs that compromise the vascularization of bone and soft tissue.

10. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure presents risks and the possibility of complications, such as infections, bleeding, allergic drug reactions and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Oss 3D_β may include, but are not limited to:

- Trans- and post-operative complications, such as: bleeding, pain, edema, abscess, fistula, dehiscence, nerve damage or compression (sensory or motor paresthesia and hyperesthesia, partial or total, transient or permanent), soft tissue damage, local or systemic infection, and other complications typical of a bone graft placement procedure;
- Sensitivity of the tissue at the surgical site;
- Pain or discomfort due to the presence of the bone graft or as a result of the surgical procedure;
- Lack of incorporation of the bone graft, or delayed healing, which can result in product failure;
- Loss of stability, displacement of the bone graft resulting from improper installation of the product or failure to incorporate it;
- Displacement or exposure of the bone graft in the oral cavity due to use in sites that may suffer stress, direct mechanical strain, exposure to constant deformation and/or tension forces, and inadequate fixation of the product;
- Allergic reaction to tricalcium β-phosphate;
- Intense fibrotic tissue reaction around the surgical area;
- Insufficient vascularization and bone formation can lead to loss of the grafted area;
- Premature suture dehiscence, edge coaptation tension can alter the healing processes and result in the failure of the surgical procedure;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques.

11. SYMBOLS

Symbols that represent standards and conformities associated with the product and its use.

	Catalog number
	Batch code.
	Unique device identifier
	Raw material
	Expiry date
	Country and date of manufacture
	Manufacturer
	Sterilization by irradiation
	Do not reuse
	Do not re-sterilize.
	Consult the instructions for use (plenum.bio/ifu)
	Do not use if packaging is damaged
	Keep dry
	Keep away from sunlight
	Single sterile barrier system with protective packaging inside

NOTE: The graphic symbols related to the labeling of the Plenum® Oss 3D_β, referenced above, meet the requirements established in ISO 15223-1:2021

12. STERILIZATION

- Plenum® Oss 3D_β is supplied in sterile condition, sterilized by Gamma Irradiation;
- Single-use medical product and cannot be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the packaging and the expiration date. Do not use the product if the packaging has been tampered with or if the expiry date has passed;
- Date of manufacture, expiry date, product code, batch number and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the product must be opened at the time of use and handled in appropriate and aseptic environments.

13. EXPIRATION DATE

- The product has a shelf life of 02 (two) years.
- The date of manufacture, expiry date, code and batch number of the product should be consulted on the product label;
- Before use, check the expiration date. Do not use the product if it has expired.

14. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of Plenum® Oss 3D_β and comply with health surveillance requirements, the surgeon and other staff must keep information regarding the device used in the patient's medical records. Additionally, this information must also be passed on to the device distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the hospital/healthcare facility to ensure the complete traceability of Plenum® Oss 3D_β by noting the code and batch of the device used in the patient's medical records, or by identifying it using one of the 3 traceability labels that come inside packaging.

The labels contain the following information necessary for the traceability of the device used:

- Manufacturer identification;
- Device name or commercial model;
- Device code;
- Batch code;
- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or by the patient to the regulatory authority and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The packaging must be intact upon receipt; do not use the product if the packaging is damaged or broken, or if the expiration date has passed. If the packaging is damaged, the product should be considered non-sterile and discarded;
- The device must be stored and transported in a clean environment, protected from sources of heat or humidity, away from the sun, free from the effects of the weather, in such a way as to prevent any damage or alteration to the packaging or its physical form that would make it impossible to use;
- The product should be stored on shelves or in drawers/cabinets in such a way that it is easy for the operator to see and handle. Do not store the product near light bulbs, so as not to affect the information on the packaging;
- The product should only be opened at the time of use in a sterile surgical field;
- The product should only be handled by qualified and trained personnel.

16. DISPOSAL OF DEVICE

- Discard the product if any of the successive packages that ensure sterility are damaged or opened;
- After use or any removal of the product from the patient, synthetic bone graft (Plenum® Oss 3D_β) that is not the subject of further analysis or study must be disposed of by healthcare establishments (hospitals, practices and clinics);
- It is recommended that Plenum® Oss 3D_β be completely disposed of so that it cannot be reused. Disposal of the product is the sole responsibility of the hospital/healthcare facility;
- After de-characterization, Plenum® Oss 3D_β must be identified as being unfit for use and disposed of in accordance with the rules for managing healthcare waste (hospitals, practices and clinics) in accordance with applicable environmental legislation;
- Leftover Plenum® Oss 3D_β that is not used during surgery should also be disposed of following the guidelines described above.

17. CUSTOMER SERVICE

If you would like to obtain the Instructions for Use in printed form and at no additional cost (including postage), please call (11) 3109- 9000, e-mail (relacionamento@plenum.bio) or visit the website (www. plenum.bio/ifu);

If the product presents any adverse event, with risk potential that could harm the patient, generates or has the potential to harm or threaten public health, is outside its specifications, or is generating any dissatisfaction, the healthcare professional must generate a technical report containing the details of the non-conformity, with all the data of the complainant, product code and batch, and patient data; notifying Plenum's Customer Relationship Service directly;

If necessary, send the cleaned, decontaminated and packaged product, duly identified and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Telephone/Fax: (11) 3109-9000; E-mail: relacionamento@plenum.bio

ANVISA also provides a channel for receiving notifications of incidents, adverse events and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA.

18. LIMITED WARRANTY

Products manufactured by Plenum have a factory guarantee against manufacturing defects. Defects in the product must be reported immediately to the manufacturer, observing the legal deadline, attaching a technical report to the product detailing the non-conformity found. This report must be generated and signed by a qualified healthcare professional and, in the event of a risk to the patient's health, must also comply with the relevant ANVISA regulations (see above under Customer Service).

The guarantee for products manufactured by Plenum is subject to strict observation and compliance with the indications for use, contraindications, precautions, warnings and other information contained in the instructions for

use. Failure to do so will exempt the manufacturer and distributor from any liability whatsoever, especially in the following cases:

- Reprocessing, resterilization or reuse of the product by the user;
- Improper use or use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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