PLENUM® OSS 3D<sub>HA</sub>  
FIT  
**Português**

PLENUM® OSS 3D<sub>HA</sub>  
FIT  
**English**

PLENUM® OSS 3D<sub>HA</sub>  
FIT  
**Español**

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TECHNICAL NAME: Bone grafts

TRADE NAME: Plenum<sup>®</sup> Oss 3D<sub>HA</sub> fit

MODEL: For a detailed description of the device, see Table 1

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E-mail: [customerservice@plenum.bio](mailto:customerservice@plenum.bio)

Website: [www.plenum.bio](http://www.plenum.bio)

ATTENTION: Read all the Instructions for Use carefully before use.

Comply with all indications, warnings and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to the occurrence of complications.

Product for Medical and Dental Use. Sterile Product. Sterilized by Gamma Irradiation. Expiration date: See device labeling. Single Use only.

Prohibited reprocessing, resterilize, and reuse.

ANVISA REGISTRATION No.: 81684340025

Device Code, Batch Number and Date of Manufacture and Sterilization: See device labeling.

REV.00

## SUMMARY






1. DETAILED DEVICE DESCRIPTION	3
2. PRODUCT COMPOSITION	4
3. FORM OF PRODUCT PRESENTATION	4
4. INDICATION, PURPOSE OR INTENDED USE OF THE PRODUCT	4
5. FUNDAMENTALS OF OPERATION AND INTENDED PERFORMANCE OF THE PRODUCT	4
6. INSTRUCTIONS FOR USE	4
7. PRECAUTIONS	6
8. WARNINGS	7
9. CONTRAINDICATIONS	7
10. POSSIBLE ADVERSE EFFECTS	7
11. SYMBOLS	8
12. STERILIZATION	8
13. EXPIRATION DATE	8
14. IDENTIFICATION AND TRACEABILITY	8
15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	9
16. DISPOSAL OF DEVICE	9
17. CUSTOMER SERVICE	9
18. LIMITED WARRANTY	9

# 1. DETAILED DEVICE DESCRIPTION

The Plenum® Oss 3D<sub>HA</sub> fit is a patient-specific (customized) bone graft based on synthetic bioceramics (hydroxyapatite) produced from the additive manufacturing process (3D printing). The virtual/digital planning and conception of the product is obtained from a planning through image processing software for design and 3D modeling, which treats three-dimensional (3D) medical images in DICOM (Digital Imaging and Communications in Medicine) format

from imaging exams (example: computed tomography), in which the bone graft is designed with complex architecture and faithful to the anatomy of the bone tissue to be reconstructed with interlocking porous macrostructure. After the medical device planning process, the virtual design is converted into a file in STL (Standard Triangle Language) extension, which associated with additive manufacturing technology, obtains the production of parts with complex geometries that corroborate personalized medical devices.

**TABLE 1:** Illustrative images, codes, commercial models, dimensional and anatomical limits and limit quantities of fixing screws per model of the Plenum Oss 3D<sub>HA</sub> fit.

Indication regions and anatomical boundaries	Illustrative figures (Patient-specific grafts)	Code	Model (Volume Dimensions)	Dimensional Limits on the x, y, and z Axes / Structural Cell	Maximum number of set screws per commercial model
<p><b>Person face drawing (AI-generated content may be incorrect)</b></p>  <p>ANATOMICAL LIMITS:  <b>Maxilla:</b> nasal cavity, maxillary sinus floor and frontal process.  <b>Mandible:</b> related to the body of the mandible - alveolar process (region of the dental alveoli) and base of the mandible (angle of the mandible, point at which the body joins)</p>		018-25-B70-00-01	0,35 cm <sup>3</sup> ≤ volume ≤ 1 cm <sup>3</sup>	<p><b>“Z”-Axis</b>                      Altura: 3 a 16 (mm)</p> <p><b>“Y”-Axis</b>                      Largura: 3 a 30 (mm)</p>	<p><b>Minimum:</b> 1 screw (for fixations in the center of the bone graft);  <b>Maximum:</b> 2 screws (for fixations at the ends of the bone graft).</p> <p>Note: The fixation screws must be positioned so that they maintain a minimum distance of 3 mm from each other and from the edges of the graft (graft border).</p>
		018-25-B70-00-02	1 cm <sup>3</sup> < volume ≤ 5 cm <sup>3</sup>	<p><b>“X”-Axis</b>                      Comprimento: 3 a 50 (mm)</p>	
		018-25-B70-00-03	5 cm <sup>3</sup> < volume ≤ 24 cm <sup>3</sup>	 <p><b>Edge:</b> 2,4mm  <b>Ø Pore:</b> 400µm</p>	

NOTE: The three-dimensional axes (x, y, and z) can be rotated depending on the normal direction of the recipient bed (depending on the geometry of the bone defect).

## Ancillary components

Plenum® Oss 3D<sub>HA</sub> fit medical devices, in order to have the intended clinical performance, must be fixed with fixation screws, in order to obtain the primary stability of the graft next to the recipient bed (the ancillary components are

not part of this family of products, they have separate ANVISA registration, and must be purchased separately).

The commercial models of attachment screws in the Plenum® Graft family are described below:

Commercial name	Technical name	Composition	Model/Dimensional	Screw Housing Hole		
				Ø for head housing	Ø for body housing 2.0 (mm)	Solid Contour Thickness
Plenum® GraftAn- visa Registration 81684340018	Non- bioabsorbable, non-sterile craniofacial bone screw	The screws of the Plenum® Graft family are manufactured in titanium alloy ELI-Ti6Al4V, according to ASTM F136-13(2017) and ABNT NBR ISO 5832-3:2017 standards	CORTICAL SCREW RT Ø1,4X4mm CORTICAL SCREW RT Ø1,4X5mm CORTICAL SCREW RT Ø1,4X6mm CORTICAL SCREW RT Ø1,4X7mm CORTICAL SCREW RT Ø1,4X8mm CORTICAL SCREW RT Ø1,4X9mm CORTICAL SCREW RT Ø1,4X10mm CORTICAL SCREW RT Ø1,4X11mm CORTICAL SCREW RT Ø1,4X12mm CORTICAL SCREW RT Ø1,4X13mm CORTICAL SCREW RT Ø1,4X14mm CORTICAL SCREW RT Ø1,4X15mm CORTICAL SCREW RT Ø1,4X16mm CORTICAL SCREW RT Ø1,4X17mm CORTICAL SCREW RT Ø1,4X18mm CORTICAL SCREW RT Ø1,4X19mm CORTICAL SCREW RT Ø1,8X4mm CORTICAL SCREW RT Ø1,8X5mm CORTICAL SCREW RT Ø1,8X6mm CORTICAL SCREW RT Ø1,8X7mm CORTICAL SCREW RT Ø1,8X8mm CORTICAL SCREW RT Ø1,8X9mm CORTICAL SCREW RT Ø1,8X10mm CORTICAL SCREW RT Ø1,8X11mm CORTICAL SCREW RT Ø1,8X12mm CORTICAL SCREW RT Ø1,8X13mm CORTICAL SCREW RT Ø1,8X14mm CORTICAL SCREW RT Ø1,8X15mm CORTICAL SCREW RT Ø1,8X16mm CORTICAL SCREW RT Ø1,8X17mm CORTICAL SCREW RT Ø1,8X18mm CORTICAL SCREW RT Ø1,8X19mm	4.5 mm	2.0 mm	0.2 mm

## 2. PRODUCT COMPOSITION


Plenum® Oss 3D<sub>HA</sub> fit is an absorbable synthetic bone graft based on hydroxyapatite (HA), manufactured according to ABNT NBR ISO 13175-3 - Implants for surgery - Calcium phosphates - Part 3: Bone substitutes of tricalcium β-phosphate and hydroxyapatite.

## 3. FORM OF PRODUCT PRESENTATION

Plenum® Oss 3D<sub>HA</sub> fit is supplied packaged in the models described in Table 01, in sterile condition, and presented as follows: 01 (one) Plenum® Oss 3D<sub>HA</sub> fit and 03 (three) identification and traceability labels in the following packages:

- Primary packaging in the form of a case composed of acrylonitrile-butadiene-styrene (ABS) frame and acetate film film, containing 01 unit of the product;
- Secondary packaging: the case containing the product is packaged in a PET blister and thermo-sealed with surgical grade Tyvek paper;
- Tertiary packaging in the form of a rigid carton paper box containing the primary and secondary packaging and 3 identification and traceability labels;

See below the illustrative images of the primary and secondary (Figure 1) and tertiary (Figure 2) packaging of a Plenum® Oss 3D<sub>HA</sub> fit model in the form in which it will be marketed.

**FIGURE 1:** Illustrative images of the primary and secondary packaging.

**FIGURE 2:** Illustrative images of tertiary packaging.  
NOTE: The Plenum® Oss 3D<sub>HA</sub> fit has no accessories.

## 4. INDICATION, PURPOSE OR INTENDED USE OF THE PRODUCT

The Plenum® Oss 3D<sub>HA</sub> fit is a synthetic bone graft, patient-specific model, whose intended use is to fill, augment or reconstruct bone defects in the bucco-maxillofacial region.

It is indicated for augmentation, reconstruction, filling and/or contouring of bones of the mid-third, alveolar ridges (maxilla and mandible). The Plenum® Oss 3D<sub>HA</sub> fit medical device is suitable for defects resulting from traumatic or vascular pathologies, bone deformities and the need for augmentation and/or remodeling for facial aesthetics and atrophic regions, which are not intrinsic to the stability of the bone structure.

## 5. FUNDAMENTALS OF OPERATION AND INTENDED PERFORMANCE OF THE PRODUCT

The Plenum® Oss 3D<sub>HA</sub> fit is a synthetic bone graft to be used in the reconstruction and/or guided regeneration of bone defects. The interconnected porous macrostructure and chemical composition of the product confer effective osteoconductivity, high hydrophilic property, controlled dissolution of the crystalline phase constituting this bone graft, being slowly replaced by newly formed bone tissue during the process of bone tissue repair or regeneration. Plenum® Oss 3D<sub>HA</sub> fit has the differential of being a bone graft that adapts perfectly to the recipient bed, especially with complex contours and geometries, increasing the contact area between them, facilitating vascularization and providing the migration of osteogenic cells into the bone graft. Plenum® Oss 3D<sub>HA</sub> fit works as a bone substitute favoring the three-dimensional reconstruction of bone defects, providing osteoconduction through its chemical composition, macro porosity (gyroid structure), perfect adaptation and integration to adjacent bone tissue, maintenance of the defect space to be reconstructed, characteristics necessary for a bone neoformation.

## 6. INSTRUCTIONS FOR USE


Carefully read all instructions for use, recommendations and guidelines before using the Plenum® Oss 3D<sub>HA</sub> fit synthetic bone graft;

- Before use, check the expiration date of the product. Do not use with the packaging violated or with the expiration date;
- The Plenum® Oss 3D<sub>HA</sub> fit is a custom bone graft and is supplied ready for use;
- The techniques for use and application of the Plenum® Oss 3D<sub>HA</sub> fit vary according to the surgeon's preference, and it is up to him to choose the therapeutic approach, application technique, planning and virtual modeling of the patient-specific product must be approved by the surgeon in charge;
- It is up to the surgeon to perform all preoperative exams and care to ensure that locally and systemically the patient is able to undergo a reconstructive surgical procedure;
- The surgeon should take a careful surgical approach to obtain enough soft tissue for total coverage and no tension to cover the Plenum® Oss 3D<sub>HA</sub> fit graft by first intention;
- Proper preparation of the site (receiving bed), where the Plenum® Oss 3D<sub>HA</sub> fit will be applied, is very important. The surgical bed area (recipient) should be scarified to allow blood perfusion and nutrition of the Plenum® Oss 3D<sub>HA</sub> fit graft,

controlled bleeding is essential to nourish the graft and the results are achieved;

- In the procedure of manipulation and fixation of the Plenum® Oss 3D<sub>HA</sub> fit, the surgeon should use only appropriate surgical instruments for bone graft fixation and sterile;
- The surgeon must consider his clinical and surgical experience to choose the most appropriate technique for each patient (individualized procedure). Therefore, a single surgical technique is not recommended to perform the procedure;
- Association with other regenerative products, i.e., complementary biomaterials (ancillary components) is the surgeon's criterion of choice and responsibility (these products are not part of the Plenum® Oss 3D<sub>HA</sub> fit family and must be purchased separately).

### Handling of Plenum® Oss 3D<sub>HA</sub> fit

1. Remove the blister containing the product from the outer packaging (carton).


2. Open the blister in sterile drape and place the primary packaging (case) plus product set on the operating table.

3. Open the case by pressing the indication where it says "PUSH". After the case is opened, carefully remove the graft with the help of tweezers or with your thumb and index fingers.

### Plenum® Oss 3D<sub>HA</sub> fit Installation Schematic Sequence

1. Initial aspect of the clinical defect.

2. Tissue detachment for exposure of the surgical bed.

3. Previous use of the surgical guide\* (bone graft analogue) to verify the adaptation of the Plenum® Oss 3D<sub>HA</sub> fit. (\*This guide is provided separately and non-sterile and must be sterilized by moist heat - autoclave, before the surgical procedure).

4. Perforation of the surgical bed to allow blood perfusion and nutrition of the Plenum® Oss 3D<sub>HA</sub> fit.

5. Fit Plenum® Oss 3D<sub>HA</sub> fit over the surgical bed and stabilization with screw.

6. Final appearance of the Plenum® Oss 3D<sub>HA</sub> fit adapted and fixed onto the surgical bed.

**ATTENTION:** The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

## 7. PRECAUTIONS

The surgeon should not use the Plenum® Oss 3D<sub>HA</sub> fit until he has carefully read all instructions for use and the information contained in the package. All usage information, warnings, precautions, and recommendations mentioned in the instructions for use should be considered. Failure to observe these points may lead to the occurrence of possible adverse events for the patient.

- The handling of the Plenum® Oss 3D<sub>HA</sub> fit should be performed exclusively by surgeons (physicians or dentists) previously trained and qualified for this procedure, with experience and knowledge about bone tissue reconstruction techniques, guided tissue regeneration and the use of bone substitutes;
- The product label shows the expiration date, the code and batch number of the product, the technical responsible and the company's data, as well as the respective symbology, allowing full traceability of the product;
- The Plenum® Oss 3D<sub>HA</sub> fit, as it is a personalized bone graft, requires virtual planning for subsequent manufacture of the patient-specific medical device, which is the responsibility and competence of the surgeon (doctor or dentist) to choose the appropriate product model with the most appropriate characteristics and dimensions, after the complete individual clinical analysis and approval of the virtual planning and design of the product model of each clinical case;
- Before starting the surgery, the surgeon must submit the patient to a rigorous oral antisepsis, to prevent the product from coming into contact with contaminating substances;
- Any product that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item Product Disposal;
- Plenum® Oss 3D<sub>HA</sub> fit is a customized product, patient-specific model, therefore, it is forbidden to use the product on another patient. The product that may be left over after the end of a surgery must be disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;
- It is recommended to install dental implants and/or apply loading at the reconstruction site with Plenum® Oss 3D<sub>HA</sub> fit after 08 months postoperatively;
- The patient should be warned that failure to comply with the instructions of the postoperative period may lead to the failure of the surgical procedure and the bone graft, requiring new surgery for revision or removal of the bone graft;
- Patients with inadequate oral hygiene, clinical conditions with hematological disorders and other uncontrolled systemic changes, can reduce the effectiveness of the regenerative process and lead to bone graft failure;
- Before the surgical procedure, the surgeon (doctor or dentist) must carry out a rigorous preoperative planning, through clinical examinations and images (radiographs or CT scans). They should also consider the need for a detailed anamnesis and complementary tests on the patient's general health (blood count, coagulogram, calcium dosage, blood glucose, etc.), since this information directly influences the body's biological response to the use of the bone graft, Plenum® Oss 3D<sub>HA</sub>, and then previously define a surgical plan and the precise and appropriate virtual planning for the tissue reconstruction/regeneration procedure;
- It is recommended to request preoperative tests according to the patient's anamnesis;
- The use of this personalized bone graft with inappropriate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results of tissue reconstruction or regeneration;
- The main and most probable causes of the lack of incorporation and

accelerated degradation of Plenum® Oss 3D<sub>HA</sub> with the consequent loss of the graft are: infections, localized or systemic diseases, insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgical site, contamination of the product in the pre and/or intraoperative or postoperative period when there is exposure of the bone graft in the bucco-maxillofacial region, mechanical instability due to lack of adequate graft fixation, bone necrosis, lack of professional training, technical error in the installation of the product, and patients who do not follow the postoperative instructions. The surgeon should consider all these aspects in his preoperative planning;

- Plenum® Oss 3D<sub>HA</sub> is only an auxiliary material for bone grafting. The success of therapy depends on the intervention technique and the patient's systemic response;
- The success of a guided tissue regeneration or bone reconstruction surgical procedure using Plenum® Oss 3D<sub>HA</sub> associated with other regenerative products, bone grafts (autogenous, homogeneous or heterogenous) or membranes is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/recipient bed incorporation;
- The patient should be informed of all indications for use, contraindications, precautions, warnings, limitations, risks, complications, and possible adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques with the use of synthetic bone substitutes;
- The patient should be informed of the risks involved in the application of excessive mechanical load or trauma affecting the operated region, which may imply failure, loss of mechanical stability, fracture, exposure and displacement of the bone graft;
- Postoperative care and correct hygiene/cleanliness of the operated site are extremely important. The patient should be instructed about the existing risks until there is complete repair of the bone and/or gingival tissue at the surgical site;
- The patient should be adequately oriented regarding postoperative care. The patient's ability and willingness to follow the surgeon's instructions is one of the most important aspects in a surgical procedure;
- For elderly patients, with motor problems and limitations, with mental problems, drug addicts, or with any other difficulty in following postoperative recommendations, this product may represent a greater risk of failure, as these patients may ignore the instructions and restrictions arising from this procedure;
- Plenum® Oss 3D<sub>HA</sub> fit has not been tested in lactating women, therefore it is not known whether the metabolites of Plenum® Oss 3D<sub>HA</sub> fit pass into breast milk;
- The surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period;
- It is important and fundamental that the surgeon explains, verbally and in writing, all the above topics to his patient, and that he has a full understanding of all these aspects;
- In order to protect the patient and their surgeon, it is important to ensure the correct and total traceability of the Plenum® Oss 3D<sub>HA</sub> fit used, by annotating the product code and batch in the patient protocol, or by fixing the identification and traceability labels that accompany the product in the patient protocol.

## 8. WARNINGS

- Plenum® Oss 3D<sub>HA</sub> fit is supplied in sterile condition, sterilized by Gamma Irradiation. Do not use the product if the packaging is violated or expired;
- SINGLE-USE product, therefore, reprocessing, resterilization and reuse are prohibited and may bring risks of failure and infection to the patient. If the product is exposed to one of these situations, it will lead to contamination and permanent damage to the product and/or cross-infection in the patient;
- The product should be kept in the original packaging until the moment of use, and should only be handled in a sterile surgical drape;
- Attention: As it is a patient-specific medical device, that is, a product dependent on the planning and analysis of imaging exams, full attention should be paid if the product is not used after 90 days from the date of manufacture. In these cases, the surgeon must request a review of the planning carried out, with a new imaging exam for the current date; if this revision shows the lack of adaptation of the product to the receiving bed and/or change in geometry, a new product must be made; since the accuracy of the product's adaptation to the receiving bed and the predictability of clinical performance may be compromised;
- An alcoholic patient or drug user may not correctly follow the postoperative instructions, and may cause treatment failure with the medical device;
- Use of the medical device in patients with mental illness or unable to follow the recommendations and care of the postoperative phase may lead to treatment failure;
- Use of the medical device in a patient unable to maintain good hygiene at the surgical site (motor or physical limitations) may result in treatment failure;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques.

## 9. CONTRAINDICATIONS

The surgeon should assess the patient's condition prudently before surgery and the use of Plenum® Oss 3D<sub>HA</sub> fit, and check for any contraindications. Contraindications may include, but are not limited to:

- Patient who has insufficient quantity and quality of gingival or soft tissue in general, who is unable to provide adequate flap coverage without tension, to avoid accidental exposure of the synthetic bone graft and compromise the procedure;
- Patients with poor bone quality or patients with type IV bone due to osteoporosis or osteopenia, which can cause early bone graft failure due to inefficient tissue repair or regeneration;
- Use in places that suffer stress, direct mechanical stress, or parts of the graft that are exposed to constant deformation forces and/or tension, such as the use of mobile prostheses on the bone graft before the recommended period of integration of the bone graft into the recipient bed of 08 months postoperatively;
- Inflammation and/or infection near or at the surgical site;
- Recent history of systemic or localized infection;
- Fever;
- Leukocytosis;

- Patient with allergy or sensitivity to hydroxyapatite;
- Patient with developmental disorders of the skeletal system;
- Pregnancy and lactation period;
- Use in patients with immunological alterations and systemic diseases that result in poor healing or in the absence of integrity of the existing tissue at the surgical site;
- Use in patients with vascular disorders or coagulation disorders;
- Calcium metabolism disorder;
- Treatment with steroids and other drugs that interfere with calcium metabolism;
- Immunosuppressive pathologies;
- Patient submitted to radiotherapy, chemotherapy or use of drugs that compromise the vascularization of tissues, bone and soft.

## 10. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure has risks and possibilities of complications. Possible adverse effects associated with the surgical procedure or the use of Plenum® Oss 3D<sub>HA</sub> fit may include, but are not limited to:

- Transoperative and postoperative complications, such as: hemorrhage, hematoma, pain, edema, abscess, fistula, dehiscence, nerve injury or compression (sensory or motor paresthesia and hyperesthesia, partial or total, transient or permanent), soft tissue injury, injury to adjacent teeth, local or systemic infection, and other complications typical of a bone graft placement procedure;
- Sensitivity of the tissue at the surgical site;
- Pain or discomfort due to the presence of the bone graft or due to the surgical procedure;
- Lack of incorporation of the bone graft, or even delay in healing, which can result in product failure;
- Loss of stability, displacement of the bone graft resulting from inadequate installation of the product or non-incorporation;
- Displacement or exposure of the bone graft at the grafted site due to use in places that may suffer stress, load, direct mechanical stress, exposure to constant deformation forces and/or tension, and inadequate fixation of the product;
- Allergic reaction to hydroxyapatite;
- Intense fibrotic tissue reaction around the surgical area;
- Insufficient vascularity and insufficient bone formation can cause loss of the grafted area;
- Premature dehiscence of sutures, tension in the coaptation of the surgical edges, or premature exposure of the Plenum® Oss 3D<sub>HA</sub> fit and the grafted area can alter the scarring processes and result in failure of the surgical procedure;

The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques with the use of synthetic bone substitutes.

## 11. SYMBOLS

Symbols that represent standards and conformities associated with the product and its use.

	Catalog number
	Use-by date
	Batch code
	Unique device identifier
	Device manufacturing material
	Do not re-use
	Country and date of manufacture
	Manufacturer
	Sterilized using Irradiation
	Do not re-sterilize
	Consult electronic instructions for use
Rev. XX	Revision No. of Instruction for Use
	Consult the instructions for use in electronic format (QR code)
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight
	Single sterile barrier system with protective packaging outside
	Patient ID
Reg. ANVISA	ANVISA device registration number
Proibido reprocessar	Reprocessing prohibited

NOTE: The graphic symbols related to the labeling of the Plenum® Oss 3D<sub>HA</sub> fit, referenced above, meet the requirements set forth in the current ISO 15223-1 standard.

## 12. STERILIZATION

- Plenum® Oss 3D<sub>HA</sub> fit is supplied in sterile condition, sterilized by Gamma Irradiation;
- Single-use medical product and cannot be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the package and the expiration date. Do not use the product if the packaging is violated or the product has expired;
- Date of manufacture, expiration date, product code, batch number and other information regarding identification and traceability are on the product labeling;
- To maintain sterility, the product should be opened only at the time of use and handled in appropriate and aseptic environments.

## 13. EXPIRATION DATE

- The validity period of the product is 01 (one) year.
- The date of manufacture, expiration date, code, and batch number of the product must be consulted on the product labeling;
- Before use check the validity, do not use the product if it has expired.

## 14. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Oss 3D<sub>HA</sub> fit and comply with the requirements of current legislation, the surgeon in charge and his team must keep information regarding the product used in the patient's medical record. In addition, this information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the hospital institution/health service establishment to ensure the complete traceability of the Plenum® Oss 3D<sub>HA</sub> fit by annotating the code and batch of the product used in the patient's medical record, or by means of one of the 3 traceability labels that accompany the package.

In addition, this information must be provided to the patient so that he has information about the product used in his procedure.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or business model of the product;
- Product code;
- Product batch number;
- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or by the patient himself to the National Health Surveillance Agency - ANVISA and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

## 15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the product if the packaging is damaged or violated, or with the expiration date expired. If the packaging is damaged, the product should be considered as non-sterile and discarded;
- The product must be stored and transported in a clean environment, protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any damage or alteration in the packaging or its physical form that makes it impossible to use;
- The storage of the product on shelves or drawers/cabinet should be in such a way that it facilitates the operator's viewing and handling. Do not store the product near lamps, so that the information contained in the package is not affected;
- The product should be opened only at the time of use in a sterile surgical field;
- The handling of the product should only be done by qualified and trained personnel to perform this procedure.

## 16. DISPOSAL OF DEVICE

- Discard the product if any of the successive packages that ensure sterility are damaged or opened;
- After use or eventual removal of the product from the patient, the synthetic bone graft (Plenum® Oss 3D<sub>HA</sub> fit) that is not the subject of further analysis or studies should be discarded by health service establishments (hospitals, offices and clinics);
- It is recommended that Plenum® Oss 3D<sub>HA</sub> fit be completely de-characterized, so that it cannot be reused. The decharacterization of the product is the sole responsibility of the hospital institution/health service establishment;
- After decharacterization, the Plenum® Oss 3D<sub>HA</sub> fit must be identified as being unfit for use and disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;
- Plenum® Oss 3D<sub>HA</sub> fit that is left over and not used during surgery should also be discarded following the guidelines outlined above.

## 17. CUSTOMER SERVICE

If you want to obtain the Instructions for Use in printed format and at no additional cost (including shipping), request it by calling (11) 3109-9000, or through the electronic address ([relacionamento@plenum.bio](mailto:relacionamento@plenum.bio)) or through the website ([www.plenum.bio/ifu](http://www.plenum.bio/ifu));

If the product presents any adverse event, with potential risk that may injure the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-conformity, with all the claimant's data, product code and lot, and patient data; by directly notifying Plenum's Customer Relationship Service;

If necessary, send the clean, decontaminated and packaged product, duly identified,

and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Phone: (11) 3109-9000; Email: [relacionamento@plenum.bio](mailto:relacionamento@plenum.bio)

ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints, related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at the electronic address [www.anvisa.gov.br](http://www.anvisa.gov.br), link NOTIVISA.

## 18. LIMITED WARRANTY

The products manufactured by Plenum have a factory warranty against manufacturing defects. The occurrence of defects in the product must be immediately reported to the manufacturer, observing the legal deadline, attaching a technical report with the details of the non-conformity found. This report must be generated and signed by a qualified and trained health professional, also observing the relevant ANVISA rules in case of risk to the patient's health (see the Customer Service topic above).

The Warranty of products manufactured by Plenum is subject to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. Failure to comply with this and follow exempts the manufacturer and distributor from any type of liability, especially in the following cases:

- Reprocessing, re-sterilization or reuse of the product by the user;
- Inappropriate use or other use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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