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PLENUM® OSS 3D_{HA}
English

PLENUM® OSS 3D_{HA}
Español

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TECHNICAL NAME: Bone grafts

TRADE NAME: Plenum[®] Oss 3D_{HA}

MODEL: Detailed description of the device see Table 1

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Read all the Instructions for Use carefully before use. Comply with all indications, warnings and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to the occurrence of complications.

Product for Medical and Dental Use. Sterile Product. Sterilized by Gamma Irradiation. Expiration date: See device labeling. Single Use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA REGISTRATION No.: 81684340024

Device Code, Batch Number and Date of Manufacture and Sterilization: See device labeling.

REV.00

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. PRODUCT COMPOSITION	4
3. DEVICE PRESENTATION FORM	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLE OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	5
7. PRECAUTIONS	6
8. WARNINGS	7
9. CONTRAINDICATIONS	7
10. POSSIBLE ADVERSE EFFECTS	8
11. SYMBOLS DESCRIPTION	8
12. STERILIZATION	9
13. EXPIRATION DATE	9
14. IDENTIFICATION AND TRACEABILITY	9
15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	9
16. DISPOSAL OF DEVICE	9
17. CUSTOMER SERVICES	9
18. LIMITED WARRANTY	10




1. DETAILED DEVICE DESCRIPTION

Plenum® Oss 3D_{HA} is a resorbable bone graft based on hydroxyapatite produced from the additive manufacturing process (3D printing), which allows a standardized interconnected porous macrostructure with pores of 400 µm,

favoring a replacement of the graft with predictable and concomitant neoformed bone tissue a slow resorption during the tissue repair or regeneration process.

Plenum® Oss 3D_{HA} is a single-use product, supplied in sterile condition and is available in the codes, dimensions and quantities described in Table 1.

TABLE 1: Illustrative image, codes and dimensions of the commercial models of Plenum® Oss 3D_{HA}

Illustrative Images	Codes	Commercial Models / Sizes	
		Length (L) x Width (W) x Height (H)	Volume
	018-25-B30-00-01	BLOCK - 10 x 10 x 03 mm	0,30 cm ³
	018-25-B30-00-02	BLOCK - 10 x 15 x 03 mm	0,45 cm ³
	018-25-B30-00-03	BLOCK - 10 x 10 x 05 mm	0,50 cm ³
	018-25-B30-00-04	BLOCK - 20 x 20 x 08 mm	3,20 cm ³
	018-25-B30-00-05	BLOCK - 20 x 20 x 10 mm	4,00 cm ³
	018-25-B30-00-06	BLOCK - 20 x 30 x 08 mm	4,80 cm ³
	018-25-B30-00-07	BLOCK - 20 x 30 x 10 mm	6,00 cm ³
	018-25-B30-00-08	WEDGE - 20 x 15 x 6 mm - 9°	1,35 cm ³
	018-25-B30-00-09	WEDGE - 20 x 15 x 8 mm - 14°	1,65 cm ³
	018-25-B30-00-10	WEDGE - 20 x 15 x 10 mm - 19°	1,95 cm ³
	018-25-B30-00-11	WEDGE - 20 x 15 x 12 mm - 24°	2,25 cm ³
	018-25-B30-00-12	WEDGE - 20 x 15 x 14 mm - 29°	2,55 cm ³
		018-25-B30-00-13	CYLINDER - 08 x 10 mm
018-25-B30-00-14		CYLINDER - 10 x 10 mm	0,80 cm ³
018-25-B30-00-15		CYLINDER - 12 x 10 mm	1,10 cm ³
018-25-B30-00-16		CYLINDER - 14 x 10 mm	1,50 cm ³
018-25-B30-00-17		CYLINDER - 08 x 20 mm	1,00 cm ³
018-25-B30-00-18		CYLINDER - 10 x 20 mm	1,60 cm ³
018-25-B30-00-19		CYLINDER - 12 x 20 mm	2,30 cm ³
018-25-B30-00-20		CYLINDER - 14 x 20 mm	3,10 cm ³

ANCILLARY COMPONENTS

Plenum® Oss 3D_{HA} medical devices can be fixed with fixation screws, in order to obtain the primary stability of the graft next to the recipient bed (the ancillary components are not part of this family of products, they have

separate ANVISA registration, and must be purchased separately). Maximum recommended torque limit for screw fixing is 10 N.cm.

The commercial models of attachment screws in the Plenum® Graft family are described below:

Commercial Name	Technical Name	Composition	Model / Dimension
Plenum® Graft ANVISA Registration no.: 81684340018	Non-absorbable, non-sterile craniofacial bone screw	The screws of the Plenum® Graft family are manufactured in titanium alloy ELI-Ti6Al4V, according to ASTM F136-13(2017) and ABNT NBR ISO 5832-3:2017 standards	CORTICAL SCREW RT Ø1,4X4mm CORTICAL SCREW RT Ø1,4X5mm CORTICAL SCREW RT Ø1,4X6mm CORTICAL SCREW RT Ø1,4X7mm CORTICAL SCREW RT Ø1,4X8mm CORTICAL SCREW RT Ø1,4X9mm CORTICAL SCREW RT Ø1,4X10mm CORTICAL SCREW RT Ø1,4X11mm CORTICAL SCREW RT Ø1,4X12mm CORTICAL SCREW RT Ø1,4X13mm CORTICAL SCREW RT Ø1,4X14mm CORTICAL SCREW RT Ø1,4X15mm CORTICAL SCREW RT Ø1,4X16mm CORTICAL SCREW RT Ø1,4X17mm CORTICAL SCREW RT Ø1,4X18mm CORTICAL SCREW RT Ø1,4X19mm CORTICAL SCREW RT Ø1,8X4mm CORTICAL SCREW RT Ø1,8X5mm CORTICAL SCREW RT Ø1,8X6mm CORTICAL SCREW RT Ø1,8X7mm CORTICAL SCREW RT Ø1,8X8mm CORTICAL SCREW RT Ø1,8X9mm CORTICAL SCREW RT Ø1,8X10mm CORTICAL SCREW RT Ø1,8X11mm CORTICAL SCREW RT Ø1,8X12mm CORTICAL SCREW RT Ø1,8X13mm CORTICAL SCREW RT Ø1,8X14mm CORTICAL SCREW RT Ø1,8X15mm CORTICAL SCREW RT Ø1,8X16mm CORTICAL SCREW RT Ø1,8X17mm CORTICAL SCREW RT Ø1,8X18mm CORTICAL SCREW RT Ø1,8X19mm

2. PRODUCT COMPOSITION

The commercial models that are part of the Plenum® Oss 3D_{HA} family are hydroxyapatite-based ($\geq 95\%$ HA).

3. DEVICE PRESENTATION FORM

Plenum® Oss 3D_{HA} is supplied individually packaged in the models described in Table 1 (see item 1.1), in sterile condition, and presented as follows:

- Primary packaging in the form of a case with a transparent lid in crystal polystyrene (GPPS), white base in polypropylene (PP) and cradle in orange silicone, containing 01 unit of the product in one of the presentations described in Table 1;
- Secondary packaging in the form of a PET blister and thermo-sealed with Tyvek-type surgical grade paper, sterilized by Gamma Irradiation;
- Outer packaging in the form of a rigid carton paper box containing the primary and secondary packaging, and 3 identification and traceability labels;

See below the illustrative images of the primary and secondary packaging (Figure 1) and external packaging (Figure 2) of a Plenum® Oss 3D_{HA} model, in the form in which it will be marketed.


FIGURE 1: Illustrative images of the primary and secondary packaging.

FIGURE 2: Illustrative images of the outer packaging.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Plenum® Oss 3D_{HA} is a synthetic bone graft, whose intended use is the filling and/or reconstruction of bone cavities or skeletal system defects that are not intrinsic to the stability of the bone structure. It is indicated for guided tissue regeneration, augmentation/reconstruction and filling of bone defects such as atrophic alveolar ridges, empty spaces (gaps), bone cavities created surgically or through traumatic injuries and bone deformities, and for bone augmentation and/or remodeling needs for facial aesthetics.

5. PRINCIPLE OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Plenum® Oss 3D_{HA} is a resorbable synthetic bone graft to be used in the reconstruction and/or guided regeneration of bone defects. The interconnected porous macrostructure and chemical composition confer effective osteoconductivity, high hydrophilic property, controlled resorption/dissolution of the crystalline phase constituting this bone graft, which is absorbed by the body gradually and replaced by newly formed bone tissue during the process of bone tissue repair or regeneration.

The compressive strength property of these bone grafts ensures a mode of use that allows drilling by drills and fixation by means of screws. The adjacent contact of the bone graft (Plenum® Oss 3D_{HA}) with the recipient bed increases the contact area between them, facilitating revascularization and providing the migration of osteogenic cells into the bone graft. In general, resorbable synthetic materials for bone regeneration are the materials of choice; since


these materials are free of possible contaminants or organic molecules from homogeneous and heterogenous grafts, which can induce immune responses. In this way, synthetic materials minimize the risk of infection, inflammation, and other postoperative complications.

6. INSTRUCTIONS FOR USE

Carefully read all instructions for use, recommendations and guidelines before using the Plenum® Oss 3D_{HA} synthetic bone graft;

- Before use, check the expiration date of the product. Do not use with the packaging violated or with the expiration date;
- Plenum® Oss 3D_{HA} is a bone graft and supplied ready for use;
- It does not require prior hydration in saline solution or the patient's own blood for a gain in effectiveness in handling the product;
- The surgeon must understand and plan all the steps of the procedure and verify the amount he intends to use. Always have the product available in quantities higher than the one you intend to use;
- The use of this bone graft with inappropriate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results of tissue reconstruction or regeneration;
- The techniques for using and applying Plenum® Oss 3D_{HA} vary according to the surgeon's preference, and it is up to him to choose the therapeutic approach, application technique, prior planning of the surgical procedure;
- The surgeon should take a careful surgical approach to obtain enough soft tissue for total coverage and no tension to cover the Plenum® Oss 3D_{HA} graft by first intention;
- Proper preparation of the site (recipient bed) where Plenum® Oss 3D_{HA} will be applied is very important. The entire area of the surgical bed (recipient) should be prepared to allow blood perfusion and nutrition of the Plenum® Oss 3D_{HA} graft, controlled bleeding is essential to nourish the graft and the results are achieved;
- In the procedure of manipulation and fixation of the Plenum® Oss 3D_{HA}, the surgeon should use only appropriate surgical instruments for bone graft fixation and sterile;
- The surgeon must consider his clinical and surgical experience to choose the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients;
- Association with other regenerative products, i.e., complementary biomaterials (ancillary components) is the surgeon's criterion of choice and responsibility (these products are not part of this product family and must be purchased separately).


Handling and Manipulation of Plenum® Oss 3D_{HA}:

1. Remove the blister containing the product from the outer packaging (carton).

2. Open the blister in sterile drape and place the primary packaging (case) plus product set on the operating table.

3. Open the sterile sleeve case by turning the lid counterclockwise.


4. Once the case is open, remove the pad with tweezers or your thumb and index fingers.

CUSTOMIZATION OF THE BLOCK AND ADAPTATION TO THE RECEIVING BED:

The Plenum® Oss 3D_{HA} are medical devices that can be adapted and adjusted, but the limits of mechanical conformation must be respected according to the table below.

Minimum thickness for screw fastening	3 mm	
Minimum graft edge thickness, for screw fixation	3 mm	
Maximum number of set screws per commercial model	<p>Minimum: 1 screw (for fixations in the center of the bone graft);</p> <p>Maximum: 2 screws (for fixations at the ends of the bone graft); Note: The fixation screws must be positioned so that they respect the minimum distance of 3 mm, between them and between the ends of the graft (graft edge).</p>	

- Customize the Plenum® Oss 3D_{HA} according to the anatomy of the recipient bed with the help of drills and sterile surgical/prosthetic discs for straight pieces (examples: 702 drills, Maxicut line models and diamond discs);
- Give adequate finishing, leaving no sharp edges and corners to avoid perforation of the flap;
- Adapt the graft to the recipient bed without leaving steps, leveling all sides of the graft to the bone defect.
- Maximum recommended torque limit for bolt fastening: 10 N.cm

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

7. PRECAUTIONS

- Surgeons must not use Plenum® Oss 3D_{HA} until he has carefully read all instructions for use and the information contained in the package. All usage information, warnings, precautions, and recommendations mentioned in the instructions for use should be considered. Failure to observe these points may lead to the occurrence of possible adverse events for the patient;
- The handling of the Plenum® Oss 3D_{HA} should be performed exclusively by surgeons (physicians or dentists) previously trained and qualified for this procedure, with experience and knowledge about bone tissue reconstruction techniques, guided tissue regeneration and the use of bone substitutes;
- The product label presents the expiration date, the code and batch number of the product, ANVISA registration number, technical responsible and the company's data, as well as the respective symbology, allowing full traceability of the product;
- The Plenum® Oss 3D_{HA} is available in several dimensions (models available in Table 01), and it is the surgeon's responsibility (physician or dentist) to choose the product with the most appropriate characteristics and dimensions, after the complete individual clinical analysis of each case;
- The wrong choice of the bone graft that will be used, as well as errors in the indication, manipulation, preparation of bone tissue or soft tissues, and installation of the product, can cause damage to the physical structure and contamination of the product and contribute to its failure;
- The decision on the graft model chosen for bone defect reconstruction depends directly on the surgeon's preference regarding the size of the bone defect, handling, usability, and surgical technique to be employed, where the objective is to maximize the contact between bone tissue and the medical device to enhance graft integration and bone healing;
- Before starting the surgery, the surgeon must submit the patient to a rigorous antisepsis of the region, to prevent the product from coming into contact with contaminating substances;
- Any product that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item Product Disposal;
- Any leftovers of Plenum® Oss 3D_{HA} are prohibited. The product that may be left over after the end of a surgery must be disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;
- It is recommended the installation of dental implants and/or application of loading at the reconstruction site with the Plenum® Oss 3D_{HA} after 08 months postoperatively;
- The patient should be warned that failure to comply with the instructions of the postoperative period may lead to the failure of the surgical procedure, requiring new surgery for revision or removal of the bone graft;
- Patients with inadequate oral hygiene, clinical conditions with hematological disorders and other uncontrolled systemic changes, can reduce the effectiveness of the regenerative process and lead to bone graft failure;
- Before the surgical procedure, the surgeon (doctor or dentist) must carry out a rigorous preoperative planning, through clinical examinations and images (radiographs or CT scans). It should also consider the need for a detailed anamnesis and complementary tests on the patient's general health (blood count, coagulogram, calcium dosage, blood glucose, etc.), since this information directly influences the body's biological response to the use of the bone graft, Plenum® Oss 3D_{HA}, and then previously define a surgical plan, with the selection of the most appropriate model for the tissue reconstruction/regeneration procedure;

- It is recommended to request preoperative tests according to the patient's anamnesis;

- The main and most likely causes of the lack of incorporation and resorption of Plenum® Oss 3D_{HA} with the consequent loss of the graft are: infections, localized or systemic diseases, inappropriate use of membranes, meshes or meshes, insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgical site, contamination of the product in the pre and/or intraoperative period or in the postoperative period when there is exposure of the bone graft, mechanical instability due to lack of adequate graft fixation, bone necrosis, lack of professional training, technical error in the installation of the product, and patients who do not follow the postoperative instructions. The surgeon should consider all these aspects in his preoperative planning;

- Plenum® Oss 3D_{HA} is only an auxiliary material for bone grafting. The success of therapy depends on the intervention technique and the patient's systemic response;

- The success of a guided tissue regeneration or bone reconstruction surgical procedure using Plenum® Oss 3D_{HA} associated with other regenerative products, bone grafts (autogenous, homogeneous or heterogenous) or membranes is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/recipient bed incorporation.

- The patient should be informed of all indications for use, contraindications, precautions, warnings, limitations, risks, complications, and possible adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques with the use of synthetic bone substitutes;

- The patient should be informed of the risks involved in the application of excessive mechanical load or trauma affecting the operated region, which may imply failure, loss of mechanical stability, exposure and displacement of the bone graft;

- Postoperative care and correct hygiene/cleanliness of the operated site are extremely important. The patient should be instructed about the existing risks until there is complete repair of the bone and/or gingival tissue at the surgical site;

- The patient should be adequately oriented regarding postoperative care. The patient's ability and willingness to follow the surgeon's instructions is one of the most important aspects in a surgical procedure;

- For elderly patients, with motor problems and limitations, with mental problems, drug addicts, or with any other difficulty in following postoperative recommendations, this product may represent a greater risk of failure, as these patients may ignore the instructions and restrictions arising from this procedure;

- Plenum® Oss 3D_{HA} has not been tested in lactating women, therefore it is not known whether the metabolites of Plenum® Oss 3D_{HA} pass into breast milk;

- The surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period, in which the stability and integrity of the product must be verified, as well as the incorporation and/or resorption at the appropriate time;

- It is important and fundamental that the surgeon explains, verbally and in writing, all the above topics to his patient, and that he has a complete and clear understanding of all these aspects;

- In order to protect the patient and their surgeon, it is important to ensure the correct and total traceability of the Plenum® Oss 3D_{HA} used, by annotating the product code and lot in the patient protocol, or by fixing the identification and traceability labels that accompany the product in the patient protocol.

8. WARNINGS

- Plenum® Oss 3D_{HA} is supplied in sterile condition, sterilized by Gamma Irradiation. Do not use the product if the packaging is violated or expired;

- SINGLE-USE product, therefore, reprocessing, resterilization and reuse are prohibited and may bring risks of failure and infection to the patient. If the product is exposed to one of these situations, it will lead to contamination and permanent damage to the product and/or cross-infection in the patient;

- The product should be kept in the original packaging until the moment of use, and should only be handled in a sterile surgical drape;

- An alcoholic patient or drug user may not correctly follow the postoperative instructions, and may cause treatment failure with the medical device;

- Use of the medical device in patients with mental illness or unable to follow the recommendations and care of the postoperative phase may lead to treatment failure;

- Use of the medical device in a patient unable to maintain good hygiene at the surgical site (motor or physical limitations) may result in treatment failure;

- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques.

9. CONTRAINDICATIONS

- The surgeon should prudently assess the patient's condition prior to surgery and the use of Plenum® Oss 3D_{HA}, and check for any contraindications. Contraindications may include, but are not limited to:

- Patient who has insufficient quantity and quality of gingival or soft tissue in general to not provide adequate flap coverage and without tension to avoid accidental exposure of the synthetic bone graft and make the procedure unfeasible;

- Patients with poor bone quality or patients with type IV bone due to osteoporosis or osteopenia, which can cause early bone graft failure due to inefficient tissue repair or regeneration;

- Use in places that suffer stress, direct mechanical stress, or parts of the graft that are exposed to constant deformation forces and/or tension, such as the use of mobile prostheses on the bone graft before the recommended period of integration of the bone graft into the recipient bed of 08 months postoperatively;

- Inflammation and/or infection near or at the surgical site;

- Recent history of systemic or localized infection;

- Fever;

- Leukocytosis;

- Patient with allergy or sensitivity to hydroxyapatite;

- Patient with developmental disorders of the skeletal system;

- Pregnancy and lactation period;

- Use in patients with immunological alterations and systemic diseases that result in poor healing or the absence of integrity of the existing tissue at the surgical site;

- Use in patients with vascular disorders or coagulation disorders;

- Calcium metabolism disorder;
- Treatment with steroids and other drugs that interfere with calcium metabolism;
- Immunosuppressive pathologies;
- Patient submitted to radiotherapy, chemotherapy or use of drugs that compromise the vascularization of tissues, bone and soft.

10. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure has risks and possibilities of complications, such as infections, bruises, hemorrhages, allergic drug reactions, and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Oss 3D_{HA} may include, but are not limited to:

- Transoperative and postoperative complications, such as: hemorrhage, hematoma, pain, edema, abscess, fistula, dehiscence, nerve injury or compression (sensory or motor paresthesia and hyperesthesia, partial or total, transient or permanent), soft tissue injury, local or systemic infection, and other complications typical of a bone graft placement procedure;
- Sensitivity of the tissue at the surgical site;
- Pain or discomfort due to the presence of the bone graft or due to the surgical procedure;
- Lack of incorporation of the bone graft, or even delay in healing, which can result in product failure;
- Loss of stability, displacement of the bone graft resulting from inadequate installation of the product or non-incorporation;
- Displacement or exposure of the bone graft due to use in places that may suffer stress, direct mechanical stress, exposure to constant deformation forces and/or tension, and inadequate fixation of the product;
- Allergic reaction to hydroxyapatite;
- Intense fibrotic tissue reaction around the surgical area;
- Insufficient vascularity and insufficient bone formation can cause loss of the grafted area;
- Premature dehiscence of sutures, tension in the coaptation of the surgical edges, or premature exposure of the Plenum® Oss 3D_{HA} and the grafted area can alter the scarring processes and result in failure of the surgical procedure;

The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the surgeon of all these complications, and of the risks arising from guided tissue regeneration surgeries and bone grafting techniques.

11. SYMBOLS DESCRIPTION

The symbols described in Table 4 appear on the instruction for use, on the packaging or on the label of the product. The symbols represent standards and conformities associated with the product and its use.

TABLE 4: Symbols that represent standards and conformities associated with the product and its use.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)

	Consult the instructions for use in electronic format (QR code)
Reg. ANVISA	ANVISA device registration number

NOTE: The graphic symbols related to the labeling of the Plenum® Oss 3D_{HA}, referenced above, meet the requirements set forth in the current ISO 15223-1 standard.

12. STERILIZATION

- Plenum® Oss 3D_{HA} is supplied in sterile condition, sterilized by Gamma Irradiation;
- Single-use medical product and cannot be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the package and the expiration date. Do not use the product if the packaging is violated or the product has expired;
- Date of manufacture, expiration date, product code, batch number and other information regarding identification and traceability are on the product labeling;
- To maintain sterility, the product should be opened only at the time of use and handled in appropriate and aseptic environments.

13. EXPIRATION DATE

- Check the expiration date of the product on the label.
- The date of manufacture, expiration date, code, and batch number of the product must be consulted on the product labeling;
- Before use check the validity, do not use the product if it has expired.

14. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Oss 3D_{HA} and comply with the requirements of current legislation, the surgeon and his team must keep information regarding the product used in the patient's medical record. In addition, this information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the hospital/health service establishment to ensure the complete traceability of the Plenum® Oss 3D_{HA}, by noting in the patient's medical record the code and batch of the product used, or by means of one of the 3 traceability labels that accompany the package.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or business model of the product;
- Product code;
- Product batch number;

- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or by the patient himself to the National Health Surveillance Agency - ANVISA and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the product if the packaging is damaged or violated, or with the expiration date expired. If the packaging is damaged, the product should be considered as non-sterile and discarded;
- The product must be stored and transported in a clean environment, protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any damage or alteration in the packaging or its physical form that makes it impossible to use;
- The storage of the product on shelves or drawers/cabinet should be in such a way that it facilitates the operator's viewing and handling. Do not store the product near lamps, so that the information contained in the package is not affected;
- The product should be opened only at the time of use in a sterile surgical field;
- The handling of the product should only be done by qualified and trained personnel to perform this procedure.

16. DISPOSAL OF DEVICE

- Discard the product if any of the successive packages that ensure sterility are damaged or opened;
- After the use or eventual removal of the product from the patient, the synthetic bone graft (Plenum® Oss 3D_{HA}) that is not the subject of further analysis or studies should be discarded by health service establishments (hospitals, offices and clinics);
- It is recommended that Plenum® Oss 3D_{HA} be completely de-characterized, so that it cannot be reused. The decharacterization of the product is the sole responsibility of the hospital institution/health service establishment;
- After decharacterization, the Plenum® Oss 3D_{HA} must be identified as being unfit for use and disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation;
- Plenum® Oss 3D_{HA} that is left over and not used during surgery should also be discarded following the guidelines outlined above.

17. CUSTOMER SERVICES

If you want to obtain the Instructions for Use in printed format and at no additional cost (including shipping), request it by phone (11) 3109-9000, electronic address (relacionamento@plenum.bio) or through the website (www.plenum.bio/ifu);

If the product presents any adverse event, with potential risk that may injure the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-conformity, with all the claimant's data, product code and lot, and patient data; by directly notifying Plenum's Customer Relationship Service;

If necessary, send the clean, decontaminated and packaged product, duly identified, and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Phone/Fax: (11) 3109-9000; Email: relacionamento@plenum.bio

ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints, related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at the electronic address www.anvisa.gov.br, link NOTIVISA.

18. LIMITED WARRANTY

The products manufactured by Plenum have a factory warranty against manufacturing defects. The occurrence of defects in the product must be immediately reported to the manufacturer, observing the legal deadline, attaching a technical report with the details of the non-conformity found. This report must be generated and signed by a qualified and trained health professional, also observing the relevant ANVISA rules in case of risk to the patient's health (see the Customer Service topic above).

The Warranty of products manufactured by Plenum is subject to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. Failure to comply with this and follow exempts the manufacturer and distributor from any type of liability, especially in the following cases:

- Reprocessing, re-sterilization or reuse of the product by the user;
- Inappropriate use or other use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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