PLENUM[®]
TISSUE ORTHO[™]
Português

PLENUM[®]
TISSUE ORTHO[™]
English

PLENUM[®]
TISSUE ORTHO[™]
Español

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Plenum[®] Tissue ortho[™]

INSTRUCTIONS FOR USE

TECHNICAL NAME: Regenerative membrane

TRADE NAME: Plenum[®] Tissue ortho[™]

MODEL: For a detailed description of the product see Table 1.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Ethylene Oxide. Expiration date: See device labeling. Single-use device. Prohibited reprocessing, resterilize, and reuse.

ANVISA DEVICE REGISTRATION NUMBER: 81684340023

Device code, batch code, and date of manufacture and sterilization: See device labeling.

SUMMARY

1. DETAILED PRODUCT DESCRIPTION	4
2. PRODUCT COMPOSITION	4
3. FORM OF PRODUCT PRESENTATION	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	5
7. HANDLING AND MANIPULATION OF PLENUM® TISSUE ORTHO™	5
8. PRECAUTIONS	6
9. WARNINGS	6
10. CONTRAINDICATIONS	7
11. POSSIBLE ADVERSE EFFECTS	7
12. SYMBOLS	7
13. STERILIZATION	8
14. EXPIRATION DATE	8
15. IDENTIFICATION AND TRACEABILITY	8
16. STORAGE, TRANSPORT AND HANDLING CONDITIONS	8
17. DISPOSAL OF DEVICE	9
18. CUSTOMER SERVICE	9
19. LIMITED WARRANTY	9

1. DETAILED PRODUCT DESCRIPTION

Plenum® Tissue ortho™ is a white absorbable synthetic membrane composed of polydioxanone, to be used in tissue regeneration. The microstructure of Plenum® Tissue ortho™ ensures high swelling efficiency, acts as a mechanical barrier promoting space maintenance, and efficient guided tissue regeneration. Plenum® Tissue ortho™ is a membrane with sub- and micrometric fibers, morphology and surface topography similar to the extracellular matrix, whose physicochemical and morphological characteristics facilitate the diffusion of biological fluids and cell adhesion.

Plenum® Tissue ortho™ is a single-use product, supplied in sterile condition and is available in the models, codes, and dimensions listed below in Table 1.

TABLE 1. Illustrative image, codes, trade name and dimensions of Plenum® Tissue ortho™.

<p>Illustrative Image</p> 	
<p>Code</p>	<p>Model / Dimensions (thickness × width × length)</p>
<p>018-25-B22-80-45</p>	<p>0.25 × 40 × 50 mm</p>
<p>018-25-B22-80-55</p>	<p>0.25 × 50 × 50 mm</p>
<p>018-25-B22-80-66</p>	<p>0.25 × 60 × 60 mm</p>
<p>018-25-B22-81-45</p>	<p>0.5 × 40 × 50 mm</p>
<p>018-25-B22-81-55</p>	<p>0.5 × 50 × 50 mm</p>
<p>018-25-B22-81-66</p>	<p>0.5 × 60 × 60 mm</p>

2. PRODUCT COMPOSITION

Plenum® Tissue ortho™ is a 100% synthetic polydioxanone-based membrane

3. FORM OF PRODUCT PRESENTATION

Plenum® Tissue ortho™ is not provided with accessories.

Plenum® Tissue ortho™ in the models described in table 01 is supplied individually packaged, in the sterile condition, and presented as follows:

- Primary packaging: polyethylene terephthalate (PET) blister and heat-sealed with surgical grade paper;
- Secondary packaging: PET/aluminum/polypropylene (PP) envelope containing 01 unit of the medical device, sterilized by Ethylene Oxide - EtO;
- Tertiary packaging: rigid carton paper box containing the primary and secondary packaging and 3 identification and traceability labels;

See below the illustrative images of the primary and secondary (Figure 1) and tertiary (Figure 2) packaging of a Plenum® Tissue ortho™ model, in the form in which it will be marketed:


FIGURE 2. Illustrative images of tertiary packaging.

FIGURE 2. Illustrative images of tertiary packaging.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Plenum® Tissue ortho™ is indicated for tissue regeneration for the treatment of chondral or osteochondral lesions or defects of the different joints of the human body.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Plenum® Tissue ortho™ is a white absorbable synthetic membrane composed of polydioxanone. The microstructure of Plenum® Tissue ortho™ ensures high swelling efficiency, promoting space maintenance and efficient tissue regeneration. This membrane has sub- and micrometric fibers, morphology and surface topography similar to the extracellular matrix, whose physicochemical and morphological characteristics facilitate the diffusion of biological fluids and cell adhesion. The physicochemical, mechanical and biological properties allow this membrane to have a dual function, performing the function of a mechanical barrier due to its microstructure and resistance, and as a matrix, providing a stable environment for cell differentiation and tissue formation.

Membrane degradation occurs by hydrolysis (breakdown of the polymer chain in the presence of water), and the physical structure and composition allow a controlled resorption within an estimated period of 6 months, this period varies according to the dimensions of the membrane, and the by-products of degradation induce minimal or no adjacent foreign body reaction.

The intended performance of the Plenum® Tissue ortho™ medical device is that it functions as a resorbable synthetic regenerative membrane for the treatment of cartilage lesions and/or defects (chondral or osteochondral), providing space maintenance and tissue repair/regeneration.

6. INSTRUCTIONS FOR USE

Carefully read all instructions for use, recommendations and guidelines before using the Plenum® Tissue ortho™ membrane;

- Before use, check the expiration date of the product. Do not use if the packaging is violated or expired;
- Plenum® Tissue ortho™ is supplied ready for use, is available in different sizes and thicknesses, requiring no prior hydration. The appropriate membrane model, once selected, should be applied directly to the area to be treated. The microstructure of the membrane has excellent adaptability to the recipient bed, and the membrane is hydrated by the patient's own blood and biological fluids, also allowing the use of fixation systems (suture thread, surgical glues, tacks or screws). The mechanical barrier formed by the membrane promotes efficient maintenance of the space to be repaired or regenerated;
- Plenum® Tissue ortho™ is available in various models and sizes, and it is the surgeon's responsibility to choose the product model with the most appropriate characteristics and dimensions, after the individual clinical analysis of each case;
- The techniques for the use and application of Plenum® Tissue ortho™ vary according to the surgeon's preference, and it is up to him to choose the therapeutic approach, model and dimensions of the product, application technique, as well as the criteria for monitoring and evaluating the results of the surgery;
- Proper preparation of the site (host bone) where Plenum® Tissue ortho™ will be applied is very important. The membrane must cover the entire area to be treated. The cavity and/or tissue to be repaired must have the recipient bed already prepared, as well as with the bleeding controlled, so that the membrane has a perfect adaptation to the recipient bed and achieves the results;
- In the procedure of manipulation and implantation of Plenum® Tissue ortho™, the surgeon should use only appropriate and sterile surgical instruments;
- Before the surgical procedure, the surgeon must carry out rigorous preoperative planning, through clinical examinations and imaging (radiographs, CT scans, or magnetic resonance imaging). You should also consider the need for a detailed anamnesis and complementary tests on the patient's general health (complete hemogram, coagulogram, calcium dosage, etc.), since this information directly influences the body's biological response to the application of Plenum® Tissue ortho™, and then previously define a surgical plan, with the selection of the model and size of the membrane most appropriate to the patient's condition;
- The surgeon must indicate the surgical technique using their clinic and surgical experience and select the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients. Surgical protocols must be carried out according to references and previous experiences of the surgeon, always considering the most appropriate choice of model and membrane size, the technique, the installation sequence, and the use of membrane fixation methods (suture thread, surgical glues, tacks or screws). These decisions must be based on conventional and established therapeutic techniques for guided tissue regeneration and bone grafting;
- The patient should be adequately oriented regarding postoperative care. The patient's ability and willingness to follow the surgeon's instructions is one of the most important aspects of the success of a surgical procedure;
- For elderly patients, those with mental problems, drug addicts, or with any other difficulty in following postoperative recommendations, this product may represent a higher risk of failure, as these patients may ignore

the instructions and restrictions arising from this procedure;


- The surgeon should adopt a correct clinical/surgical planning described in the specific literature on tissue regeneration for the pre- and immediate and late postoperative periods;
- The surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period, in which the stability, integrity of the membrane, as well as integration and absorption at the appropriate time must be verified;
- In order to protect the patient and his surgeon, it is important to ensure the correct and total traceability of the Plenum® Tissue ortho™ used, by annotating the product code and batch in the patient protocol, or by fixing the identification and traceability labels that accompany the product in the patient protocol.

7. HANDLING AND MANIPULATION OF PLENUM® TISSUE ORTHO™

ON THE SIDE OF THE SECONDARY PACKAGING IT IS MARKED OPEN HERE IN THE PLACE WHERE THERE IS THE PERFORATION TO FACILITATE THE OPENING OF THE PACKAGING.

ON THE SIDE OF THE PRIMARY PACKAGE THERE IS A FLAP ON THE SURGICAL PAPER TO FACILITATE OPENING.

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

8. PRECAUTIONS

The surgeon should not use Plenum® Tissue ortho™ until he has carefully read all instructions for use and the information contained in the package. All usage information, precautions, and recommendations mentioned in these instructions for use should be considered. Failure to observe these points may lead to the occurrence of adverse events for the patient.

- In the event that the product is exposed to a temperature outside the specified maximum limit, there is no guarantee of product performance. This condition can cause the beginning of product degradation (the regenerating membrane becomes rigid and brittle), so the appearance of the membrane must be carefully observed before use;
- If physical aspects of membrane degradation are observed during handling, the product should not be used;
- Any product that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the Product Disposal;
- The surgeon must understand and plan all the stages of the surgical procedure, and verify that the membrane model he intends to use meets the requirements of the surgical protocol. Always have membranes available in quantities greater than what you intend to use;
- It is recommended to request preoperative tests according to the patient's anamnesis;
- The surgeon must be aware of the patient's clinical conditions. In the preoperative examination, the target site where the membrane will be installed must be rigorously verified through clinical and imaging examinations. The surgeon should consider the general clinical status of the patient, analyze local or systemic factors that may interfere with the process of bone healing, soft tissues, and membrane integration and absorption;
- Before starting surgery, the surgeon must submit the patient to a rigorous antisepsis at the surgical procedure site to prevent the membrane from coming into contact with contaminating substances;
- Patients with clinical conditions with hematological disorders and other uncontrolled systemic alterations have restrictions on the use of this product;
- Any leftovers of Plenum® Tissue ortho™ are prohibited. The membrane that may be left after the end of a surgery must be discarded, in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;
- The wrong choice of the membrane that will be used, as well as errors in the indication, manipulation, preparation of bone tissue or soft tissues, and application of the product, can cause damage to the physical structure and contamination of the membrane and contribute to the failure of its functioning;
- The use of the Plenum® Tissue ortho™ membrane with inappropriate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results in tissue repair or regeneration;
- The patient should be warned that failure to comply with the instructions of the postoperative period may lead to failure of the surgical procedure, requiring new surgery to review or remove the membrane;
- The main and most probable causes of the lack of integration and absorption of the membrane with the consequent loss of surgery are: infections, localized or systemic diseases (not controlled), low quality of the bone substitute used associated with the membrane, insufficient blood or nutrient supply, trauma or overload at the surgical site, undersized membrane,

contamination of the membrane in the intraoperative or postoperative period of the membrane in the recipient bed, mechanical instability of the membrane, bone necrosis, lack of professional training, technical error in the application of the membrane, and patients who do not follow postoperative instructions. The surgeon should consider all these aspects in his preoperative planning;

- Plenum® Tissue ortho™ has not been tested in lactating women, therefore it is not known whether the metabolites of Plenum® Tissue ortho™ pass into breast milk;
- Plenum® Tissue ortho™ is only an auxiliary material for tissue regeneration. Use alone does not guarantee the success of the treatment. The success of the therapy depends on the intervention technique and the patient's systemic capacity to regenerate the tissues;
- The success of a tissue regeneration surgical procedure using Plenum® Tissue ortho™ associated with bone substitutes (alloplastic materials or bone grafts) is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/recipient bed integration, as these extrinsic factors can be affected by numerous biological factors, clinical and biomechanical. These factors can cause failure and early loss of the membrane and compromise the surgical procedure. Therefore, faithful obedience to all indications for use, contraindications, precautions, warnings and other observations contained in these instructions for use are essential for the success of the surgical procedure;
- The patient should be informed of all the limitations, risks, complications, warnings, and possible adverse effects arising from the tissue regeneration procedure and techniques with the use of resorbable membranes;
- The patient should be informed of the risks involved in procedures with the use of regenerating membranes where the operated site will suffer stress and direct mechanical stress, for example, exposure to constant deformation forces and/or tension before the rest/rehabilitation period determined by the physician;
- The patient should be informed of the risks involved in the application of excessive mechanical load or trauma affecting the operated region, which may imply failure, loss of mechanical stability, exposure, and displacement of the membrane;
- It is important and fundamental that the surgeon explains all the above topics to his patient, and that he has a complete and clear understanding of all these aspects.

9. WARNINGS

- Plenum® Tissue ortho™ is supplied in a sterile condition, sterilized by Ethylene Oxide – EtO. Do not use the product if the packaging is violated or expired;
- The product label presents the expiration date, the code and batch number of the product, maximum and minimum limit of storage and transport temperature, ANVISA registration number, technical responsible and the company's data, as well as the respective symbology, allowing full traceability of the membrane;
- SINGLE-USE product, therefore, reprocessing, resterilization and reuse are prohibited and may bring risks of failure and infection to the patient. If the product is exposed to one of these situations, it will lead to contamination and permanent damage to the product and/or cross-infection in the patient;
- The product should be kept in the original packaging until the time of use, and should only be handled in sterile surgical drapes;

- This product is restricted to medical use. Plenum® Tissue ortho™ should only be used by trained, qualified surgeons with experience and knowledge of the therapeutic techniques of tissue reconstructions and regenerations using regenerative medical devices;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from tissue regeneration surgeries with the use of resorbable membranes.

10. CONTRAINDICATIONS

The surgeon should prudently assess the patient's condition prior to surgery and the use of Plenum® Tissue ortho™, and check for any contraindications listed below. Contraindications may include, but are not limited to:

- Patient with poor bone quality due to osteoporosis or osteopenia, which can cause early membrane failure due to inefficient tissue regeneration;
- Infection near or at the site of surgery;
- Recent history of systemic or localized infection;
- Fever;
- Leukocytosis;
- Patient with allergy or sensitivity to the polymer polydioxanone;
- Pregnancy and lactation period;
- Use in patients with immunological alterations and systemic diseases that result in poor healing or the absence of integrity of the existing tissue at the surgical site;
- Use in patients with vascular disorders or coagulation disorders;
- Use in an alcoholic patient or drug user who cannot correctly follow postoperative instructions;
- Use in patients with mental illness or unable to follow the recommendations and care of the postoperative phase;
- Use in patients unable to maintain good hygiene at the surgical site;
- Immunosuppressive pathologies;
- Patient undergoing radiotherapy or use of drugs that compromise the vascularization of tissues, bone and soft.

11. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure has risks and possibilities of complications, such as infections, hemorrhages, allergic drug reactions and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Tissue ortho™ may include, but are not limited to:

- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from tissue regeneration surgeries with the use of resorbable membranes;
- Transoperative and postoperative complications, such as: hemorrhage, pain, edema, abscess, fistula, dehiscence, hematoma, nerve injury or compression (sensory or motor paresthesia, partial or total, transient or permanent), soft tissue injury, local or systemic infection, and other

complications typical of a procedure for placing resorbable membranes;

- Sensitivity of the tissue at the surgical site;
- Pain or discomfort due to the surgical procedure;
- Lack of membrane integration and absorption, or even delayed healing, which can result in product failure;
- Loss of stability, displacement or failure of the membrane resulting from non-integration and early absorption, or even caused by inadequate application of the device;
- Displacement or exposure of the membrane in the host site due to use at sites that may be subject to stress, direct mechanical force, exposure to constant deformation and/or tension, and inadequate fastening of the device;
- Allergic reaction to polydioxanone;
- Intense fibrotic tissue reaction around the surgical area;
- Tissue necrosis, insufficient vascularization, and insufficient tissue formation, which can cause loss of the grafted area and membrane;
- Premature dehiscence of sutures, tension in the coaptation of the surgical edges can alter the healing processes and result in the failure of the surgical procedure.

12. SYMBOLS

Table 2 describes the symbols that may be on product labeling. Please check the physical labeling for the symbols applied to the product. The symbols express information related to the identification and use of the product.

TABLE 2. Symbols used for labeling to express information related to the identification and use of the device

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system

	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

13. STERILIZATION

- Plenum® Tissue ortho™ is supplied in a sterile condition, sterilized by Ethylene Oxide - EtO;
- Single-use medical device and cannot be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the package and the expiration date. Do not use the device if the packaging is violated or the device has expired;
- Date of manufacture, sterilization expiration date, code, batch number, and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments.

14. EXPIRATION DATE

- Indicated on device labeling;
- Before use, check the expiration date, do not use the device if it has expired.

15. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of Plenum® Tissue ortho™ and to comply with the requirements of current legislation, the surgeon or someone from his team must keep the information regarding the membrane used in the patient's record. In addition, this information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the surgeon in charge and/or the hospital/health service establishment to ensure the complete traceability of Plenum® Tissue ortho™, by annotating in the patient's medical record the code and batch of the product used, or by means of one of the 3 traceability labels that accompany the packaging, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that he has information about the product used in his procedure.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or business model of the product;
- Product code;
- Product batch number;
- ANVISA registration number.
- Traceability information is necessary for notification by the health service and/or by the patient himself to the National Health Surveillance Agency - ANVISA and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

16. STORAGE, TRANSPORT AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the product if the packaging is damaged or violated, or with the expiration date expired. If the packaging is damaged, the product should be considered as non-sterile and discarded;
- The product must be stored and transported in a clean environment (maximum temperature limit of 25°C and minimum of 4 °C), protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any damage or alteration in the packaging or its physical form that makes it impossible to use;
- It is the sole responsibility of the users of the product, health service establishments (hospitals, clinics and offices) to store, preserve and handle the product in a controlled environment at a maximum temperature of 25 °C and a minimum of 4 °C (24 hours/7 days);
- The device must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the device near lamps, so that the information contained in the packaging is not affected;
- The device should only be opened at the time of use in a sterile surgical field;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

17. DISPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After the use of the device or their removal from the patient, the membranes that are not subject to analysis or subsequent studies must be disposed of by the health service establishment (hospitals, offices, and clinics);
- It is recommended that the membrane be completely disassembled to ensure that it cannot be reused. The health service institution/establishment is solely responsible for disassembling the device;
- After disassembly, the Plenum® Tissue ortho must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- Remnant residuals of the Plenum® Tissue ortho not used during surgery must also be discarded following the guidelines above.

18. CUSTOMER SERVICE

- If the device presents any adverse event with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, and batch, and patient data; directly notifying the Plenum Customer Relations Service;
- If necessary, send the cleaned, decontaminated packaged device, properly identified and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Telephone: +55 11 3109- 9000; E-mail: customerservice@plenum.bio;
- ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints related to the use of devices and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA

19. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the Customer Services section above) must be observed in cases involving a risk to the patient's health.

The Warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing, resterilization or reuse of the device by the user;
- Improper use or use other than in a surgical setting;
- Device not used according to instructions for use.

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