PLENUM[®]
BONE ORTHO[™]
Português

PLENUM[®]
BONE ORTHO[™]
English

PLENUM[®]
BONE ORTHO[™]
Español

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TECHNICAL NAME: Bone grafts

TRADE NAME: Plenum[®] Bone ortho[™]

MODEL: For a detailed description of the device, see Table 1.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA DEVICE REGISTRATION N°.: 81684340026

Device Code, batch code, and date of manufacture and sterilization: See device labeling.

REV.00

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. PRODUCT COMPOSITION	3
3. DEVICE PRESENTATION FORM	3
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	3
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	4
7. PRECAUTIONS	5
8. WARNINGS	6
9. CONTRAINDICATIONS	6
10. POSSIBLE ADVERSE EFFECTS	6
11. SYMBOLS DESCRIPTION	7
12. STERILIZATION	7
13. EXPIRATION DATE	8
14. IDENTIFICATION AND TRACEABILITY	8
15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	8
16. DISPOSAL OF DEVICE	8
17. CUSTOMER SERVICES	8
18. LIMITED WARRANTY	9

1. DETAILED DEVICE DESCRIPTION

The bone graft Plenum® Bone ortho™ is a synthetic biphasic bioceramic composed of hydroxyapatite and tricalcium β -phosphate (HA: β -TCP) in a 70:30% ratio. The macro and microstructure, surface area and porosity of the Plenum® Bone ortho™ confer effective osteoconductivity, high hydrophilic property, controlled absorption of the constituent phases of the bioceramics, which are absorbed by the body, and replaced by newly formed tissue during the process of repair or regeneration of bone tissue. The composition of Plenum® Bone ortho™ reflects in the behavior of controlled bone replacement, i.e., a faster dissolution of the crystalline phase of β -TCP in relation to the crystalline phase of HA. Plenum® Bone ortho™ is a single-use product, supplied in sterile condition and is available in the codes, dimensions, and quantities described in Table 1.

Code	Model/Presentation (granule size/ quantity)	Model Content
018-60-B10-10-00	1,0 - 2,0 mm - 5 cc	1 unit of 5 cc
018-60-B10-10-01	1,0 - 2,0 mm - 10 cc	1 unit of 10 cc
018-60-B10-10-02	1,0 - 2,0 mm - 20 cc	1 unit of 20 cc
018-60-B10-10-03	1,0 - 2,0 mm - 30 cc	1 unit of 30 cc
018-60-B10-40-00	2,0 - 4,0 mm - 5 cc	1 unit of 5 cc
018-60-B10-40-01	2,0 - 4,0 mm - 10 cc	1 unit of 10 cc
018-60-B10-40-02	2,0 - 4,0 mm - 20 cc	1 unit of 20 cc
018-60-B10-40-03	2,0 - 4,0 mm - 30 cc	1 unit of 30 cc
018-60-B10-60-00	4,0 - 5,6 mm - 5 cc	1 unit of 5 cc
018-60-B10-60-01	4,0 - 5,6 mm - 10 cc	1 unit of 10 cc
018-60-B10-60-02	4,0 - 5,6 mm - 20 cc	1 unit of 20 cc
018-60-B10-60-03	4,0 - 5,6 mm - 30 cc	1 unit of 30 cc

2. PRODUCT COMPOSITION

Plenum® Bone ortho™ is an absorbable synthetic bone substitute based on biphasic bioceramic composed of hydroxyapatite and tricalcium β -phosphate (HA: β -TCP) in the 70:30% ratio, following the guidelines of the current ABNT NBR ISO 13175-3 standard.

3. DEVICE PRESENTATION FORM

The Plenum® Bone ortho™ has no accessories.

The Plenum® Bone ortho™ is supplied packaged in the models described in table 01, in sterile condition, and presented as follows:

- Primary packaging: polyethylene terephthalate (PET) blister and thermo-sealed with Tyvek-type surgical grade paper;
- Secondary packaging: PET blister and thermo-sealed with surgical grade Tyvek paper and sterilized by Gamma Irradiation;
- Tertiary packaging: rigid carton paper envelope containing the primary and secondary packaging and 3 identification and traceability labels;

See below the illustrative images of the primary and secondary (Figure 1) and tertiary (Figure 2) packaging of a Plenum® Bone ortho™ model in the form in which it will be marketed.


FIGURE 1. Illustrative images of the primary and secondary packaging.


FIGURE 2. Illustrative images of tertiary packaging.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Plenum® Bone ortho™ is a synthetic bone graft, which is intended for bone filling and/or reconstruction, being indicated for the regeneration of bone defects resulting from congenital, post-traumatic, post-surgical problems, in grafting procedures in the areas of traumatology, orthopedic and spine surgeries (cage filling), oncology, and other degenerative processes of the skeletal system, which are not intrinsic to the stability of the bone structure, being gradually replaced by newly formed bone.

Plenum® Bone ortho™ can also be used in procedures as an auxiliary product in hip, knee, shoulder and other joint arthroplasties, associated with membranes, plates, meshes, and articulated and non-articulated prostheses. The Plenum®

Bone ortho™ is partially reabsorbed by the body, and replaced by newly formed tissue during the process of bone tissue repair or regeneration.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The Plenum® Bone ortho™ has a macro and microstructure that mimics the trabecular bone tissue providing surface area, percentage of porosity and unique physicochemical characteristics that promote excellent osteoconductive property, high hydrophilicity and adaptability to the recipient bed. The porosity of Plenum® Bone ortho™ ensures efficient osteogenic cell growth by promoting bone regeneration at the bone defect site, since microporosity and macroporosity facilitate the diffusion of biological fluids and cell adhesion. In addition, the chemical composition of Plenum® Bone ortho™ (hydroxyapatite and tricalcium β -phosphate in a 70:30% ratio) confers high hydrophilic properties to the product, which makes it easily manipulable (with fluids such as saline solution or the patient's own blood) in various applications for bone reconstruction. The replacement of the Plenum® Bone ortho™ with the newly formed bone tissue occurs during the process of bone tissue repair or regeneration. In general, absorbable alloplastic materials for bone repair are the materials of choice; since these materials are free of possible contaminants or organic molecules from homogeneous and heterogenous grafts, which can induce immune responses. In this way, alloplastic materials minimize the risk of infection, inflammation, and other postoperative complications. The composition of Plenum® Bone ortho™ reflects in the controlled absorption behavior. As a faster absorption of β -TCP promotes a continuous increase in the porosity of the material, which allows an integration of the material/tissue due to cell adhesion and proliferation and growth of blood vessels (subsequent deposition of the newly formed bone), HA provides stability of the newly formed bone volume for a longer period.

6. INSTRUCTIONS FOR USE

- Carefully read all instructions for use, recommendations and guidelines before using the Plenum® Bone ortho™ synthetic bone graft;
- Before use, check the expiration date of the product. Do not use with the packaging violated or with the expiration date;
- The Plenum® Bone ortho™ is supplied ready for use. It is available in different quantities and in three particle size ranges (1 -2 mm, 2 - 4 mm and 4 - 5.6 mm), does not require prior hydration, but hydration is recommended (with fluids such as saline solution or the patient's own blood) for a gain in effectiveness in handling the product;
- Plenum® Bone ortho™ is available in various models and sizes of granules (see Table 1), and it is the surgeon's responsibility to choose the product model with the most appropriate characteristics and dimensions, after the individual clinical analysis of each case;
- The techniques for using and applying the Plenum® Bone ortho™ vary according to the surgeon's preference, and it is up to the surgeon to choose the therapeutic approach, model and dimensions of the product, application technique, use of complementary biomaterials, as well as the criteria for monitoring and evaluating the results of the surgery;
- Proper preparation of the site (receiving bed) where the Plenum® Bone ortho™ will be applied is very important. The entire area to be treated must be filled with the product, as well as bleeding must be controlled, so that the bone graft has a perfect adaptation to the recipient bed and achieves the

expected results;

- In the procedure of manipulation and implantation of the Plenum® Bone ortho™, the surgeon should use only appropriate and sterile surgical instruments;
- Before the surgical procedure, the surgeon must carry out rigorous preoperative planning, through clinical examinations and imaging (radiographs or CT scans). You should also consider the need for a detailed anamnesis and complementary tests on the patient's general health (such as complete blood count, coagulogram, calcium dosage), since this information directly influences the body's biological response to the use of the bone graft, Plenum® Bone ortho™, and then previously define a surgical plan, with the selection of the size of the granules and the most appropriate quantity for the patient's condition;
- The surgeon must consider his clinical and surgical experience to choose the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients;
- The surgical protocol should be carried out according to the surgeon's previous references and experiences, always considering the most appropriate choice of granule size and quantity of the product, technique, grafting sequence, and the use of membranes, plates, meshes, and other metallic and/or polymeric devices when necessary, always based on conventional and established therapeutic techniques of bone grafting;
- Postoperative care and correct hygiene/cleanliness of the operated site are extremely important. The patient should be instructed about the existing risks until there is complete repair of the bone and/or gingival tissue at the surgical site;
- The patient should be adequately oriented regarding postoperative care. The patient's ability and willingness to follow the surgeon's instructions is one of the most important aspects in a surgical procedure;
- For elderly patients, those with mental problems, drug addicts, or with any other difficulty in following postoperative recommendations, this product may represent a higher risk of failure, as these patients may ignore the instructions and restrictions arising from this procedure;
- The surgeon should adopt a correct clinical/surgical planning described in the specific literature on bone regeneration and bone grafting techniques for the immediate and late pre- and postoperative periods;
- The surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period, in which the stability and integrity of the product must be verified, as well as the incorporation and/or absorption at the appropriate time;
- In order to protect the patient and his surgeon, it is important to ensure the correct and total traceability of the Plenum® Bone ortho™ used, by annotating the product code and batch in the patient's medical record, or by fixing the identification and traceability labels that accompany the product in the patient's medical record.

FIGURE 3. Illustrative images of the handling and manipulation of the Plenum® Bone ortho™.

ATTENTION. The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

7. PRECAUTIONS

- Surgeons must not use the Plenum® Bone ortho™ until they have carefully read all instructions for use and the information contained in the package. All usage information, precautions, and recommendations mentioned in these instructions for use should be considered. Failure to observe these points may lead to the occurrence of adverse events for the patient;
- The handling of the Plenum® Bone ortho™ should be performed exclusively by surgeons previously trained and qualified for this procedure, with experience and knowledge about grafting techniques and the use of bone substitutes;
- It is recommended to request preoperative tests according to the patient's anamnesis;
- Any product that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item Product Disposal;
- The surgeon must understand and plan all the steps of the procedure and verify the amount he intends to use. Always have the product available in quantities higher than the one you intend to use;
- The surgeon must be aware of the patient's clinical conditions. In the preoperative examination, the site where the Plenum® Bone ortho™ will be placed must be rigorously verified through clinical and imaging examinations. The surgeon must consider the general clinical status of the patient, analyze local or systemic factors that may interfere with the process of bone and soft tissue healing, and the incorporation and resorption of bone substitute;
- Before starting the surgery, the surgeon must submit the patient to a rigorous antisepsis of the external region where the incision will be made to access the region of interest (recipient bed), to prevent the product from coming into contact with contaminant substances;
- The patient should be warned that failure to comply with the instructions of the postoperative period may lead to failure of the surgical procedure, requiring a new surgery to remove the bone graft;
- Patients with clinical conditions with hematological disorders and other uncontrolled systemic alterations have restrictions on the use of this product;
- Forbidden to use any leftovers of Plenum® Bone ortho™. The product that may be left over after the end of a surgery must be disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;
- The wrong choice of the bone graft that will be used, as well as errors in the indication, manipulation, preparation of bone tissue or soft tissues, and installation of the product, can cause damage to the physical structure and contamination of the product and contribute to its failure;
- The use of this bone graft with inappropriate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results in tissue repair or regeneration;
- The main and most probable causes of the lack of incorporation and resorption of the Plenum® Bone ortho™ with the consequent loss of surgery are: infections, localized or systemic diseases, inappropriate use of membranes, meshes or meshes, insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgical site, contamination of the product in the pre and/or intraoperative period or in the postoperative period when there is exposure of the bone substitute in the cavity Oral disease, mechanical instability of the granules, bone necrosis, lack of professional training, technical error in the installation of the product, and patients who do not follow the postoperative instructions. The surgeon should consider all these

aspects in his preoperative planning;

- Plenum® Bone ortho™ is only an auxiliary material for bone replacement. Use alone does not guarantee the success of the treatment. The success of the therapy depends on the intervention technique and the patient's reaction;
- The success of a bone reconstruction/regeneration surgical procedure using Plenum® Bone ortho™ associated with other substitutes, bone grafts or membranes is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/receptor bed incorporation, as these extrinsic factors can be affected by numerous biological factors, clinical and biomechanical. These factors can cause failure and early loss of the product and compromise the surgical procedure. Therefore, faithful obedience to all indications for use, contraindications, precautions, warnings and other observations contained in these instructions for use are essential for the success of the surgical procedure;
- The patient should be informed of all the limitations, risks, complications, warnings, and possible adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques with the use of synthetic bone substitutes;
- The patient should be informed of the risks involved in the application of excessive mechanical load or trauma affecting the operated region, which may imply failure, loss of mechanical stability, exposure, and displacement of the particles of the synthetic bone substitute;
- It is important and fundamental that the surgeon explains all the above topics to his patient, and that he has a complete and clear understanding of all these aspects.

8. WARNINGS

- This product is restricted to medical use. The Plenum® Bone ortho™ should only be used by trained, qualified surgeons with experience and knowledge of the therapeutic techniques of tissue reconstructions and regenerations using regenerative medical devices;
- Plenum® Bone ortho™ is supplied in a sterile condition, sterilized by Gamma Irradiation. Do not use the product if the packaging is violated or expired;
- The product label presents the expiration date, the code and batch number of the product, ANVISA registration number, technical responsible and the company's data, as well as the respective symbology, allowing full traceability of the product;
- SINGLE-USE product, therefore, reprocessing, resterilization and reuse are prohibited and may bring risks of failure and infection to the patient. If the product is exposed to one of these situations, it will lead to contamination and permanent damage to the product and/or cross-infection in the patient;
- The product should be kept in the original packaging until the moment of use, and should only be handled in a sterile surgical drape;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from tissue regeneration surgeries with the use of resorbable membranes.

9. CONTRAINDICATIONS

The surgeon should prudently assess the patient's condition prior to surgery and the use of Plenum® Bone ortho™, and check for any contraindications

listed below. Contraindications may include, but are not limited to:

- Patient who has insufficient bone or soft tissue quantity and quality to provide primary stability of the synthetic bone graft;
- Patient with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the bone graft due to inefficient tissue repair or regeneration;
- Use in places that suffer stress, direct mechanical stress, or the granules are exposed to constant deformation forces and/or tension;
- Patients who do not have sufficient soft tissues for complete coverage of the bone graft;
- Infection near or at the site of surgery;
- Recent history of systemic or localized infection;
- Fever;
- Leukocytosis;
- Patient with allergy or sensitivity to the raw materials of biphasic bioceramics: hydroxyapatite and tricalcium β -phosphate;
- Patient with developmental disorders of the skeletal system;
- Pregnancy and lactation period;
- Calcium metabolism disorder;
- Treatment with steroids and other drugs that interfere with calcium metabolism;
- Use in patients with immunological alterations and systemic diseases that result in poor healing or the absence of integrity of the existing tissue at the surgical site;
- Patient with vascular disorders or coagulation disorders;
- Alcoholic patient or drug user who cannot correctly follow postoperative instructions;
- Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;
- Patient unable to maintain good hygiene at the surgical site;
- Immunosuppressive pathologies;
- Patient submitted to radiotherapy or use of drugs that compromise the vascularization of tissues, bone and soft;

10. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure has risks and possibilities of complications, such as infections, hemorrhages, allergic drug reactions and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Bone ortho™ may include, but are not limited to:



- The surgical procedure for the installation of bone grafts or substitutes can bring complications in the trans and postoperative periods, such as: hemorrhage, pain, edema, hematoma, abscess, fistula, dehiscence, nerve injury or compression (sensory or motor paresthesia, partial or total, transient or permanent), soft tissue injury, local or systemic infection, and other complications typical of a procedure for the placement of bone grafts or substitutes;
- Sensitivity of the tissue at the surgical site;

- Pain or discomfort due to the presence of the bone graft/substitute or due to the surgical procedure;
- Lack of incorporation of the bone graft, or even delay in healing, which can result in product failure;
- Loss of stability, displacement of the granules of the bone graft resulting from non-incorporation and early partial absorption, or even caused by inadequate installation of the product;
- Allergic reaction to the raw material of Plenum® Bone ortho™;
- Intense fibrotic tissue reaction around the surgical area;
- Bone necrosis, insufficient vascularization, and insufficient bone formation, which can cause loss of the grafted area;
- Premature dehiscence of sutures, tension in the coaptation of the surgical edges, or premature exposure of the Plenum® Bone ortho™ and the grafted area can alter the healing processes and result in failure of the surgical procedure;
- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from tissue regeneration surgeries and bone grafting techniques with the use of absorbable synthetic bone grafts/substitutes.

11. SYMBOLS DESCRIPTION

TABLE 2. Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside

	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)
	Consult the instructions for use in electronic format (QR code)
	ANVISA device registration number

NOTE. The graphic symbols related to the labeling of the Plenum® Bone ortho™, referenced above, meet the requirements set forth in the current ISO 15223-1 standards.

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

12. STERILIZATION

- Plenum® Bone ortho™ is supplied in sterile condition, sterilized by Gamma Irradiation;
- Single-use medical product and cannot be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the package and the expiration date. Do not use the product if the packaging is violated or the product has expired;
- Date of manufacture, expiration date, product code, batch number and other information regarding identification and traceability are on the product labeling;

- To maintain sterility, the product should be opened only at the time of use and handled in appropriate and aseptic environments.

13. EXPIRATION DATE

- The shelf life of the product indicated on the label;
- The date of manufacture, expiration date, code, and batch number of the product must be consulted on the product labeling;
- Before use check the expiration date, do not use the product if it has expired.

14. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Bone ortho™ and comply with the requirements of current legislation, the surgeon or someone from his team must keep in the patient's medical record the information regarding the product used. In addition, this information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the surgeon in charge and/or the hospital/health service establishment to ensure the complete traceability of the Plenum® Bone ortho™, by noting in the patient's medical record the code and batch of the product used, or by means of one of the 3 traceability labels that accompany the packaging, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that he has information about the product used in his procedure.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or business model of the product;
- Product code;
- Product batch number;
- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or by the patient himself to the National Health Surveillance Agency - ANVISA and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the product if the packaging is damaged or violated, or with the expiration date expired. If the packaging is damaged, the product should be considered as non-sterile and discarded;
- The product must be stored and transported in a clean environment, protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any

damage or alteration in the packaging or its physical form that makes it impossible to use;

- The storage of the product on shelves or drawers/cabinet should be in such a way that it facilitates the operator's viewing and handling. Do not store the product near lamps, so that the information contained in the package is not affected;
- The product should be opened only at the time of use in a sterile surgical field;
- The handling of the product should only be done by qualified and trained personnel to perform this procedure.

16. DISPOSAL OF DEVICE

- Discard the product if any of the successive packages that ensure sterility are damaged or opened;
- After the use or eventual removal of the product from the patient, the synthetic bone graft (Plenum® Bone ortho™) that is not the object of further analysis or studies must be discarded by health service establishments (hospitals, offices and clinics);
- It is recommended that the Plenum® Bone ortho™ be completely de-characterized, so that it cannot be reused. The decharacterization of the product is the sole responsibility of the hospital institution/health service establishment;
- After decharacterization, the Plenum® Bone ortho™ must be identified as being unfit for use and disposed of in accordance with the standards of waste management of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation;
- Leftover Plenum® Bone ortho™ that is not used during surgery should also be discarded following the guidelines described above.

17. CUSTOMER SERVICES

- If the device presents any adverse event, with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, and batch, and patient data; directly notifying the Plenum Customer Relations Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and properly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Telephone: +55 11 3109- 9000; Email: customerservice@plenum.bio;
- ANVISA also provides a channel available to report incidents, adverse events and technical complaints related to the use of products and services subject to health regulations, which can be found on the website of the National

Health Regulatory Agency - ANVISA at the electronic address www.anvisa.gov.br, using the NOTIVISA link.

18. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health professional and addressed to regulatory authority standards (see the Customer Services section) when it is observed that issues involve a risk to the patient's health.

The warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing, re-sterilization, or reuse of the device by the user;
- Improper use or use other than in a surgical environment;
- Device not used according to instructions for use.

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