PLENUM[®] GUIDE
Português

PLENUM[®] GUIDE
English

PLENUM[®] GUIDE
Español

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Plenum[®] Guide

INSTRUCTIONS FOR USE

TECHNICAL NAME: Regenerative membrane

TRADE NAME: Plenum[®] Guide

MODEL: For a detailed description of the device, see Table 1.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA REGISTRATION N°: 81684340007.

Device code, batch code, and date of manufacture and sterilization: See device labeling.

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	3
3. DEVICE PRESENTATION FORM	3
4. DEVICE RECOMMENDATIONS, PURPOSE OR INTENDED USE	4
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	4
7. PRECAUTIONS, RESTRICTIONS AND WARNINGS	5
8. CONTRAINDICATIONS	6
9. POSSIBLE ADVERSE EFFECTS	7
10. SYMBOLS DESCRIPTION	7
11. STERILIZATION	8
12. EXPIRATION DATE	8
12. IDENTIFICATION AND TRACEABILITY	8
14. STORAGE, TRANSPORT AND HANDLING CONDITIONS	8
15. DIPOSAL OF DEVICE	9
16. CUSTOMER SERVICES	9
17. LIMITED WARRANTY	9

1. DETAILED DEVICE DESCRIPTION

The Plenum® Guide is a white absorbable-synthetic membrane composed of poly(dioxanone) used in the reconstruction and/or guided tissue regeneration of bone and gingival defects in the craniomaxillofacial region, particularly intraosseous and peri-implant soft tissue defects. The microstructure of the Plenum® Guide ensures high swelling efficiency, acts as a mechanical barrier promoting the maintenance of spaces, and offers efficient guided tissue regeneration. The membrane contains sub and micrometric fibers, morphology, and a surface topography similar to the extracellular matrix, the physicochemical and morphological characteristics of which facilitate the diffusion of biological fluids and cell adhesion. In general, resorbable alloplastic materials are the materials of choice for bone repair since they are concomitantly replaced by neoformed tissue. Plenum® Guide can be used as a mechanical barrier to maintain the defect space and as a membrane that increases the volume augmentation/reconstruction of bone and soft tissues. The membrane's physical structure and composition allow for controlled reabsorption within 6 to 12 months.

The Plenum® Guide is a single-use device supplied in sterile conditions, and is available in the models, codes, and dimensions listed below in Table 1.

TABLE 1. Illustrative image, codes, models, and dimensions of Plenum® Guide.

Codes	Models/Dimensions (thickness x width x length)
018-25-B20-80-15	0.25x15x20 mm
018-25-B20-80-25	0.25x25x25 mm
018-25-B20-80-30	0.25x30x40 mm
018-25-B20-81-15	0.5x15x20 mm
018-25-B20-81-25	0.5x25x25 mm
018-25-B20-81-30	0.5x30x40 mm
018-25-B20-82-15	1.0x15x20 mm
018-25-B20-82-25	1.0x25x25 mm
018-25-B20-82-30	1.0x30x40 mm
018-25-B20-83-15	2.0x15x20 mm
018-25-B20-83-25	2.0x25x25 mm
018-25-B20-83-30	2.0x30x40 mm

2. DEVICE COMPOSITION

Plenum® Guide is a 100% synthetic polydioxanone-based membrane.

3. DEVICE PRESENTATION FORM

The Plenum® Guide is not provided with accessories.

The Plenum® Guide, in the models described in the table above, is supplied in a single package in sterile condition and presented in the following manner:

- Primary packaging is a polyethylene terephthalate (PET) blister heat-sealed with surgical-grade paper;
- Secondary packaging is a PET/aluminum/PP envelope containing one (01) unit of the medical device, sterilized using Ethylene Oxide - EtO;
- External packaging is a rigid cardboard box containing the primary and secondary packaging and three (03) identification and traceability labels;

It is the responsibility of the hospital/health care facility to ensure complete traceability of the Plenum® Guide, through the annotation in the patient's chart of the code and batch code of the device used or identification through one of the three (03) traceability labels accompanying the package;

A printed version of the Instructions for Use (including shipping) can be requested at no additional cost at the telephone number +55 11 3109-9000, or through the electronic address (customerservice@plenum.bio), or through the website (www.plenum.bio/ifu);

- Illustrative images of the primary and secondary (Figure 1) and external packaging (Figure 2) of commercialized Plenum® Guide models can be found below.

FIGURE 1: ILLUSTRATIVE IMAGES OF PRIMARY AND SECONDARY PACKAGING.

FIGURE 2: ILLUSTRATIVE IMAGE OF EXTERNAL PACKAGING.

M3 Health Ind. Com. de Prod. Med.
Odont. e Correlatos S.A.

Avenida Ain Ata, 640 - Lote 15 - Quadra B
Jundiaí/SP-Brasil -CEP 13.212-213
PABX +55 11 3109-9000
WWW.plenum.bio
Resp. Tec: Alberto Blay
CRO-SP 50517

4. DEVICE RECOMMENDATIONS, PURPOSE OR INTENDED USE

The Plenum® Guide is a resorbable synthetic membrane composed of polydioxanone, indicated for use in guided tissue regeneration, intraosseous defects resulting from congenital, post-traumatic, post-surgical complications or diseases, such as periodontal and peri-implant lesions that are not intrinsic to bone structure stability. It can also be used associated with bone grafts or synthetic bone substitutes for augmentation/reconstruction of the alveolar ridge used in guided tissue regeneration. This synthetic membrane acts as a barrier that facilitates the maintenance of the bone cavity or the space created for the augmentation of the alveolar ridge. The Plenum® Guide is reabsorbed by the body within 6 to 12 months and replaced by neoformed tissue during the repair or regeneration.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The Plenum® Guide is a synthetic resorbable polydioxanone membrane used in the reconstruction and/or guided regeneration of bone defects and soft tissues in the craniomaxillofacial region. The microstructure membrane

works as a mesh and as a mechanical barrier promoting the maintenance of spaces and an efficient tissue repair/regeneration process. The Plenum® Guide contains sub and micrometric fibers and a morphology and surface topography similar to the extracellular matrix, the physicochemical and morphological characteristics of which facilitate the diffusion of biological fluids and cell adhesion. The Plenum® Guide can be used as a mechanical barrier for defect space maintenance and as a membrane that increases the volume augmentation/reconstruction of bone and soft tissues.


6. INSTRUCTIONS FOR USE

- Carefully read the instructions for use, recommendations, and guidelines before using the Plenum® Guide membrane;
- Verify the device's expiration date before use. Do not use devices if the packaging has been damaged or if the expiration date has expired;
- The Plenum® Guide is supplied ready to use and does not require prior hydration. It is available in different sizes and thicknesses. The selected membrane must be placed directly onto the area that is to be treated. The membrane microstructure offers excellent adaptability near the host bone, and the membrane is hydrated by the patient's blood and biological fluids while also allowing for the use of fixation systems (dental bone tacks and screws). The mechanical barrier generated by the membrane promotes efficient maintenance of the space to be repaired or regenerated;
- The Plenum® Guide is available in different models and sizes, and it is up to the surgeon (physician or dentist) to choose the device model with the most appropriate characteristics and dimensions after the individual clinical analysis of each case;
- The techniques employed for the use and application of Plenum® Guide vary according to the preference of the surgeon (physician or dentist), and they are responsible for choosing the therapeutic approach that is used, the model of the device, the application technique, the use of a fixation system or complementary biomaterials, as well as the criteria for monitoring and evaluating the results of surgery;
- Adequate preparation of the site (host bone), where the Plenum® Guide will be applied, is very important. The membrane must cover the entire area to be treated. The cavity and/or bone or soft tissue and the grafted area must have a homogeneous surface and the final format of the desired recipient bed already prepared, as well as with controlled bleeding, so that the membrane has a perfect adaptation to the host site and achieve the expected results;
- The surgeon shall only use appropriate and sterile surgical instruments during procedures for manipulation and implantation of Plenum® Guide;
- Before the surgical procedure, the surgeon (physician or dentist) shall carry out rigorous pre-operative planning through clinical examinations and imaging (radiographs or computed tomography). They must also consider the need for a detailed history and complementary exams regarding the general health of the patient (complete hemogram, coagulogram, calcium dosage, etc.), since this information directly influences the body's biological response to the installation of the Plenum® Guide and then define a surgical plan in advance, selecting the membrane size and quantity most appropriate to the condition of the patient;
- The surgeon must indicate the surgical technique using their clinic and surgical experience and select the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients. Surgical protocols must be carried out according to references and previous experiences of the surgeon, always considering the most appropriate choice of model and membrane size,

the technique, the installation sequence, and the use of membrane fixation methods (sutures, dental bone tacks, or screws). These decisions must be based on conventional and established therapeutic techniques for guided tissue regeneration and bone grafting;

- Postoperative care and proper hygiene/cleanliness of the operated site are extremely important. The patient must be instructed about existing risks until there is the complete integration of bone and/or soft tissues and membrane reabsorption;
- The operated site's postoperative care and proper hygiene/cleanliness are critical. The patient must be instructed regarding existing risks until complete integration of bone and/or soft tissues and membrane reabsorption;
- Elderly patients, patients with mental health problems, individuals with chemical dependency, or any other difficulty in adhering to postoperative recommendations may represent a higher risk of failure of this device, since these patients may ignore instructions for this procedure;
- The surgeon must adopt the proper clinical/surgical planning described in the specific literature for guided bone/tissue regeneration and bone grafting techniques for immediate and subsequent pre-and postoperative periods;
- The surgeon must guide the patient, who must understand the importance of monitoring during the postoperative period, in which the stability, integrity, and positioning of the membrane, as well as the incorporation and/or reabsorption at the appropriate time, must be verified;
- In order to protect the patient and their surgeon, it is essential that the correct and complete traceability of the Plenum® Guide used to be guaranteed by noting the device code and batch in the patient's records or by attaching it to the identification and traceability labels that come with the device to these records.

Handling and Manipulation of the Plenum® Guide:

ATTENTION: The figures are merely illustrative. They do not represent the actual dimensions and characteristics of the device.

7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

- Surgeons must not use the Plenum® Guide before carefully read the instructions for use and information provided on the packaging. All information on the use of membranes and the warnings, precautions and recommendations included in these instructions for use must be taken into consideration. A failure to observe these points may result in complications during procedures;
- The use of this device is restricted to medical and dental applications. The Plenum® Guide must be exclusively used by surgeons (physicians or dentists) trained, qualified, with experience and knowledge of therapeutic techniques for guided tissue regeneration and use of bone substitutes;
- It is recommended that preoperative examinations be requested according to the knowledge of the patient's condition;
- The Plenum® Guide is supplied in a sterile condition, sterilized using Ethylene Oxide - EtO. Do not use devices if the packaging has been damaged or if the expiration date has expired;
- The Plenum® Guide is a SINGLE-USE device. Reprocessing, reesterilization and reuse of the device are prohibited and may promote risks

of procedure failure and patient infection. Exposure of the device to one of these situations will result in contamination and permanent damage to the device and/or cross-infection in the patient;

- Device labeling include the date of manufacture, device code, batch code, ANVISA (Brazilian National Healthy Regulatory Agency) device registration number, technical product manager and company information, in addition to the respective symbology that allows the device to be correctly traced;
- The device must be kept in the original packaging until it is used and must only be handled in a sterile surgical fields;
- Any device that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item DIPOSAL OF DEVICE;
- The surgeon must understand and plan all steps involved in the surgical procedure and verify if the membrane model they intend to use meets the requirements for preoperative planning. Always ensure that the device is available in quantities that exceed those expected to be used;
- Surgeons must be aware of the patient's clinical conditions. During preoperative examinations, clinical examinations and imaging of the site at which Plenum® Guide will be placed must be rigorously verified. Surgeons must consider the patient's general clinical condition and analyze local or systemic factors that may interfere with the process of bone and soft tissue healing, and membrane integration and resorption;
- Before starting the surgery, the surgeon must submit the patient to a rigorous oral antisepsis process to prevent the device from coming into contact with contaminating substances;
- Patients must be warned that a failure to comply with instructions for the post-operative period may lead to failure of the surgical procedure, requiring additional surgery to adjust or remove the membrane;
- The use of this device by patients with inadequate oral hygiene, clinical conditions involving hematological disorders, and other uncontrolled systemic alterations are restricted;
- The use of any remnant Plenum® Guide materials is prohibited. Any remnant device residuals shall be disposed of after the end of the surgery following health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- Improper selection of the membrane that will be used, as well as mistakes in the indicated use, handling, preparation of bone or soft tissues, and installation of the device, may result in damage to the device's structure, which may contribute to failure;
- The use of inappropriate techniques and procedures with this membrane and inadequate surgical conditions may cause harm to the patient and compromise results for tissue repair and regeneration;
- The primary and most likely causes of the lack of incorporation and reabsorption of the membrane with the consequent compromising of surgery include: infections, localized or systemic diseases, inappropriate use of membranes or meshes, insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgery site, contamination of the device pre- and/or trans-operative or postoperative whenever there is exposure of the membrane in the oral cavity, mechanical instability of granules, bone necrosis, lack of professional training, technical error in installation of the device, and a failure to adhere to postoperative instructions on the part of patients. The surgeon must consider all these aspects during preoperative planning;
- The Plenum® Guide has not been tested on lactating women. Therefore, it is unknown whether the metabolites contained in Plenum® Guide are passed on to breast milk.

- The Plenum® Guide is merely intended for use as an auxiliary material in bone regeneration. Exclusive use of the Plenum® Guide alone does not guarantee the success of the treatment. The success of therapy depends on the intervention technique used and the patient's reaction;

- The success of a surgical procedure of guided tissue regeneration using the Plenum® Guide associated with bone substitutes (alloplastic materials or bone grafts) is related to the correct indication and selection of the device, correct surgical technique, preparation, the stability of the device, the period and quality of the integration of biomaterial(s)/receiving sites, since these extrinsic factors can be affected by numerous biological, chemical and biomechanical factors. These factors may cause failure and early loss of the device and compromise the surgical procedure. Following the recommendations, contraindications, precautions, warnings, and other observations contained in these instructions for use are, therefore, essential for the success of the procedure.

- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques involving the use of resorbable membranes;

- Patients must be informed of the risks involved in procedures involving the use of bone substitutes in which the operation site will be subject to stress and direct mechanical forces, for example, exposure to constant deformation and/or tension;

- Patients must be informed of the risks involved in the application of excessive mechanical load or trauma that affects the operated region, which may imply failure, loss of mechanical stability, exposure and displacement of the membrane;

- It is of the utmost importance that the surgeon explains the entirety of the above topics to their patient, and that patients in turn have a clear and complete understanding of these procedural aspects.

- It is important to inform patients of the need to maintain periodic follow-ups with the surgeon to ensure the integrity and proper positioning of the membrane in relation to the estimated time required for its complete degradation.

8. CONTRAINDICATIONS

Surgeons shall prudently assess the patient's condition before surgery and use of the Plenum® Guide, checking for any of the contraindications reported below. Contraindications may include, but are not limited to:

- Patients with an insufficient quantity and general quality of hard or soft tissue in general to provide stability and primary mechanical attachment to the membrane;
- Patients with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the bone membrane due to inefficient tissue regeneration;
- Use of the membrane in locations subject to stress, direct mechanical forces, or the membrane being exposed to constant deformation and/or tension forces;
- Patients that do not have enough soft tissue to provide complete coverage of membrane;
- Inflammation and/or infection near or at the surgical site;
- Recent history of systemic or localized infection;
- Fever;

- Leukocytosis;
- Patients that present signs of allergy or hypersensitivity to the composition of the raw materials used in manufacturing of Plenum® Guide (see “DEVICE COMPOSITION”);
- Patient with developmental disorders of the skeletal system;
- Pregnancy and lactation period;
- Use in patients with immunological changes and systemic diseases that result in poor healing or in the absence of integrity of existing tissue at the surgical site;
- Calcium metabolism disorders;
- Treatment with steroids and other drugs that have interfered with calcium metabolism;
- Patient with vascular disorders or coagulation disorders;
- Alcoholic patient or drug user who cannot correctly follow postoperative instructions;
- Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;
- Patient unable to maintain proper hygiene at the surgical site;
- Immunosuppressive pathologies;
- Patient undergoing radiotherapy or using medications that compromise the vascularization of hard and soft tissues;
- All contraindications included in the available dental and scientific literature for the use of resorbable synthetic membranes.

9. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved during all surgical procedures. Every surgical procedure presents risks and the possibility of complications, such as infections, hemorrhaging, allergic reactions to drugs and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Guide may include, but are not limited to:

- The surgical procedure involved in installation of resorbable membranes may cause complications during transitory and postoperative periods, such as: hemorrhaging, pain, edema, abscesses, fistula, dehiscence, injury or compression of nerves (partial or total sensory or motor paresthesia that is transient or permanent), soft tissue injury, local or systemic infection, and other typical complications for placement procedures for resorbable membranes;
- Tissue sensitivity at the surgical site;
- Pain or discomfort due to the presence of the membrane or due to the surgical procedure;
- Lack of integration and reabsorption of the membrane, or delays in healing, which may result in device failure;
- Loss of stability, displacement or failure of the membrane from failure in incorporation and early partial reabsorption, as well as inadequate installation of the device;
- Displacement or exposure of the membrane in the oral cavity due to use at sites that may be subject to stress, direct mechanical force, exposure to constant deformation and/or tension, and inadequate fastening of the device;

- Failure of the procedure due to exposure and contamination of the membrane;
- Allergic reaction to the raw materials used in the membrane;
- Intense fibrotic tissue reaction around the surgical site;
- Bone necrosis, insufficient vascularization, reabsorption of the bone substitute, and insufficient bone formation, which may result in loss of the grafted area and membrane;
- Premature dehiscence in sutures, tension in the coaptation of surgical edges, or early exposure of membrane and the grafted area may alter cystic processes and result in failure of the surgical procedure;
- The patient must be informed and understand that complications are possible in all surgical procedures. The surgeon must inform the patient of these complications and the risks associated with guided tissue regeneration surgeries and bone grafting techniques involving the use of resorbable membranes.

10. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 2. Symbols used for labeling to express information related to the identification and use of the device

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Do not re-sterilize
	Do not use if package is damaged
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Temperature limit
	Do not re-use
	Consult electronic instructions for use

UDI	Unique device identifier
MAT	Device manufacturing material
Rx only	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
Proibido reprocessar	Reprocessing prohibited (Brazilian Requirement)

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

11. STERILIZATION

- The Plenum® Guide is supplied in sterile condition, sterilized using Ethylene Oxide - EtO;
- The Plenum® Guide is a single-use medical device and must not be reprocessed, re-sterilized or reused;
- Before use, check the integrity of the packaging and the expiration date. Do not use the device if the package is damaged or if the expiration date has expired;
- Date of manufacture, sterilization expiration date, code, batch code and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments;
- The opening of the packaging for use of the device must only be done by qualified and trained personnel to perform this procedure.

12. EXPIRATION DATE

- Indicated on device labeling;
- Before use, check the expiration date, do not use the device if it has expired.

13. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Guide and comply with health surveillance requirements, the surgeon or other staff must keep the information regarding the membrane used in the patient's medical records. Additionally, this information must also be passed on to the device distributor and patient, in order to complete the traceability cycle. The information necessary for traceability is related to the membrane used, surgical procedure, and patient, as follows:

- Surgeon's name
- Hospital/clinic/office identification
- Manufacturer and distributor identification
- Surgery date

- Device code
- Batch code
- Quantity used
- ANVISA device registration number

The responsible surgeon and their team must make use of the information contained in the three (03) traceability labeling supplied with the membrane, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that they have information about the device used in their procedure.

Labels contain the following information necessary for the traceability of the membrane used.

- Manufacturer's identification
- Device name or commercial model
- Device code
- Batch number
- ANVISA device registration number
- Traceability information is necessary for notification by the health service and/or by the patient to the regulatory authority and to the manufacturer when complaints, technical complaints, and/or adverse events occur for the carrying out of the appropriate investigations.

14. STORAGE, TRANSPORT AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the device if the packaging is damaged or if the expiration date has expired. If the packaging is damaged, the device must be considered as non-sterile and discarded;
- The product must be stored and transported in a clean environment, protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any damage or alteration in the packaging or its physical form that makes it impossible to use;
- It is the sole responsibility of the users of the product, health service establishments (hospitals, clinics and offices) to store, preserve and handle the product in a controlled environment under refrigeration. Temperature limits 2°C to 15 °C (35.6 °F to 59 °F);
- The device must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the device near lamps, so that the information contained in the packaging is not affected;
- The device should only be opened at the time of use in a sterile surgical field;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

15. DIPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After the use of the device or their removal from the patient, the membranes that are not subject to analysis or subsequent studies must be disposed of by the health service establishment (hospitals, offices, and clinics);
- It is recommended that the membrane be completely disassembled to ensure that it cannot be reused. The health service institution/establishment is solely responsible for disassembling the device;
- After disassembly, the Plenum® Guide must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- Remnant residuals of the Plenum® Guide not used during surgery must also be discarded following the guidelines above.

16. CUSTOMER SERVICES

- If the device presents any adverse event with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, batch code, and patient data; directly notifying the Plenum Customer Relations Service;
- If necessary, send the cleaned, decontaminated packaged device, properly identified and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Zip Code 13212-213 - Jundiaí - SP - Brazil. Telephone: +55 11 3109- 9000; Email: customerservice@plenum.bio;
- ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints related to the use of devices and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA.

17. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the Customer Services section above) must be observed in cases involving a risk to the patient's health.

The Warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing, re-sterilization or reuse of the device by the user;
- Improper use or use other than in a surgical setting;
- Device not used according to instructions for use.

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