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PLENUM<sup>®</sup> OSS<sub>hp</sub>  
**English**

PLENUM<sup>®</sup> OSS<sub>hp</sub>  
**Español**

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TECHNICAL NAME: Bone Grafts

TRADE NAME: Plenum<sup>®</sup> Oss<sub>hp</sub>

MODEL: For a detailed description of the device, see Table 1.

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E-mail: [customerservice@plenum.bio](mailto:customerservice@plenum.bio)

Website: [www.plenum.bio](http://www.plenum.bio)

ATTENTION: Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA REGISTRATION N°: 81684340008.

Code, lot, and date of manufacture and sterilization: See device labeling.

Rev. 07

## SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	4
3. DEVICE PRESENTATION FORM	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	5
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	5
6. INSTRUCTIONS FOR USE	5
7. PRECAUTIONS	7
8. WARNINGS	7
9. CONTRAINDICATIONS	8
10. POSSIBLE ADVERSE EFFECTS	8
11. SYMBOLS DESCRIPTION	8
12. STERILIZATION	9
13. EXPIRATION DATE	9
14. IDENTIFICATION AND TRACEABILITY	9
15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	10
16. DISPOSAL OF DEVICE	10
17. CUSTOMER SERVICES	10
18. LIMITED WARRANTY	10




# 1. DETAILED DEVICE DESCRIPTION

The Plenum® Oss<sub>hp</sub> (high porosity, hp) bone graft is a synthetic biphasic bioceramic composed of hydroxyapatite and tricalcium β-phosphate (HA:β-TCP) at a proportion of 70:30%. The macro and microstructure, surface area (approx. 2 m<sup>2</sup>/g), and percentage of porosity (approx. 50%) of Plenum® Oss<sub>hp</sub> provide effective osteoconductivity, high hydrophilic properties, controlled re-absorption of the constituent phases of the bioceramic which are partially reabsorbed by the body and replaced with neoformed tissue during the process of bone tissue repair or regeneration. Plenum® Oss<sub>hp</sub> promotes efficient osteogenic cell growth and bone regeneration at the bone defect site. Similarly, the microporosity and surface roughness of Plenum® Oss<sub>hp</sub> facilitate the diffusion of biological fluids and cell adhesion. Plenum® Oss<sub>hp</sub> is highly hydrophilic and easy to manipulate (using fluids such as saline, platelet-rich plasma or fibrin, or the patient's blood), making it suitable for various applications for bone reconstruction.

In general, resorbable alloplastic materials are the materials of choice for bone repair since these materials are free of possible contaminants or organic molecules from allografts and xenografts, which can induce immune responses. Alloplastic materials, therefore, minimize the risk of infection, inflammation, and other postoperative complications. The bone graft Plenum® Oss<sub>hp</sub> works as a synthetic matrix for the increase in volume/reconstruction of bone tissue and for the maintenance of the defect space and is gradually replaced by the newly formed bone. The controlled re-absorption of the constituent phases of the bioceramic ensures gradual replacement with newly formed bone tissue, providing stability and support to the newly formed bone. While the fastest reabsorption of β-TCP promotes a continuous increase in the porosity of the material, which allows for the integration of the material/tissue due to cell adhesion and proliferation and growth of the blood vessels (subsequent deposition of the newly formed bone), HA provides stability of the newly formed bone volume for a longer period.

Plenum® Oss<sub>hp</sub> is single-use only, supplied in sterile conditions and is available in the device codes, dimensions, and quantities described in Table 1.

**TABLE 1.** Device images, codes, dimensions, and quantities of Plenum® Oss<sub>hp</sub>

Device image	Codes	Models
<b>Granule</b>		<b>(Size / Quantity)</b>
	018-25-B10-25-25	250 - 500 μm / 0,25g
	018-25-B10-25-05	250 - 500 μm / 0.5g
	018-25-B10-25-10	250 - 500 μm / 1.0g
	018-25-B10-25-20	250 - 500 μm / 2.0g
	018-25-B10-50-25	500 - 1000 μm / 0,25 g
	018-25-B10-50-05	500 - 1000 μm / 0.5g
	018-25-B10-50-10	500 - 1000 μm / 1.0g
	018-25-B10-50-20	500 - 1000 μm / 2.0g
	018-25-B10-10-25	1 - 2 mm / 0,25 g
	018-25-B10-10-05	1 - 2 mm / 0.5g
	018-25-B10-10-10	1 - 2 mm / 1.0g
	018-25-B10-10-20	1 - 2 mm / 2.0g
	<b>Block</b>	
	018-25-B10-00-41	BLOCK 7,5X10X03mm
	018-25-B10-00-01	BLOCK 10X10X03mm
	018-25-B10-00-02	BLOCK 10X15X03mm
	018-25-B10-00-42	BLOCK 10X15X05mm
	018-25-B10-00-03	BLOCK 10X10X05mm
	018-25-B10-00-04	BLOCK 20X20X08mm
	018-25-B10-00-05	BLOCK 20X20X10mm
	018-25-B10-00-06	BLOCK 20X30X08mm
018-25-B10-00-07	BLOCK 20X30X10mm	
<b>Wedge</b>		<b>Length (L) x Width (W) x Height (H) - angle</b>
	018-25-B10-00-43	WEDGE 20 x 10 x 6 mm - 9°
	018-25-B10-00-08	WEDGE 20 x 15 x 6 mm - 9°
	018-25-B10-00-09	WEDGE 20 x 15 x 8 mm - 14°
	018-25-B10-00-10	WEDGE 20 x 15 x 10 mm - 19°
	018-25-B10-00-11	WEDGE 20 x 15 x 12 mm - 24°
	018-25-B10-00-12	WEDGE 20 x 15 x 14 mm - 29°
<b>Cylinder</b>		<b>Diameter (Ø) x Height (A)</b>
	018-25-B10-00-44	CYLINDER Ø 08 x 08 mm
	018-25-B10-00-13	CYLINDER Ø 08 x 10 mm
	018-25-B10-00-14	CYLINDER Ø 10 x 10 mm
	018-25-B10-00-15	CYLINDER Ø 12 x 10 mm
	018-25-B10-00-16	CYLINDER Ø 14 x 10 mm
	018-25-B10-00-17	CYLINDER Ø 08 x 20 mm
	018-25-B10-00-18	CYLINDER Ø 10 x 20 mm
	018-25-B10-00-19	CYLINDER Ø 12 x 20 mm
	018-25-B10-00-20	CYLINDER Ø 14 x 20 mm

## 2. DEVICE COMPOSITION

Plenum® Oss<sub>hp</sub> is a resorbable synthetic bone substitute based on biphasic bioceramics composed of hydroxyapatite and tricalcium β-phosphate (HA:β-TCP) at a proportion of 70:30%.

## 3. DEVICE PRESENTATION FORM

The Plenum® Oss<sub>hp</sub> is not provided with accessories.

Plenum® Oss<sub>hp</sub> is supplied packaged in the models described in table 1, in a sterile condition, and presented as follows:

### Plenum® Oss<sub>hp</sub> - Granule:

- Primary packaging is a polyethylene terephthalate (PET) blister thermo-sealed with surgical grade paper;
- Secondary packaging is a case with a transparent crystal polystyrene (GPPS) cover, and white polypropylene (PP) base, containing one (01) blister for the device using one of the presentations described above;
- Tertiary packaging is a PET blister thermo-sealed with surgical grade paper and sterilized using Gamma Irradiation;
- External packaging is a rigid cardboard box containing the primary, secondary, and tertiary packaging, in addition to three (03) identification and traceability labels.

### Plenum® Oss<sub>hp</sub> - Block / Wedge / Cylinder:

- Primary packaging is a case with a transparent crystal polystyrene (GPPS) cover, and white polypropylene (PP) base, containing one (01) blister for the device using one of the presentations described above;
- Secondary packaging is a PET blister thermo-sealed with surgical grade paper and sterilized using Gamma Irradiation;
- External packaging is a rigid cardboard box containing the primary, secondary, and tertiary packaging, in addition to three (03) identification and traceability labels.
- The hospital/health service establishment is responsible for guaranteeing the complete traceability of Plenum® Oss<sub>hp</sub>, the device through means of annotation of the code and batch code of the device used in the patient's medical records, or identification through one of the three (03) traceability labels that accompany the packaging;
- A printed version of the Instructions for Use (including shipping) can be requested at no additional cost at the telephone number +55 11 3109-9000 or through the electronic address (customerservice@plenum.bio) or the website (www.plenum.bio);
- See below for illustrative images of the primary, secondary , tertiary and external packaging (Figures 1, 2 and 3) of the Plenum® Oss<sub>hp</sub> models, in the forms in which they will be marketed.


Figure 1: Illustrative images of Primary, Secondary, and Tertiary Packaging - Plenum® Oss<sub>hp</sub> - Granule


Figure 2: Illustrative Images of External Packaging - Plenum® Oss<sub>hp</sub> - Granule


Figure 3: Illustrative images of the primary, secondary and outer packaging of Plenum® Oss<sub>hp</sub> - Block / Wedge / Cylinder

## 4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Plenum® Oss<sub>hp</sub> is a synthetic bone graft for filling and/or reconstruct bone cavities or defects in the skeletal system that are not intrinsic to the stability of the bone structure. It is indicated for guided tissue regeneration, augmentation/reconstruction and filling of bone defects such as atrophic alveolar ridge, long bones, spine, empty spaces (gaps), bone cavities created surgically or through traumatic injuries, bone deformities, and for facial bone augmentation and/or remodeling needs.

Plenum® Oss<sub>hp</sub> can also be used for alveolar ridge augmentation/reconstruction associated with resorbable or non-resorbable membranes or meshes for guided tissue regeneration. This synthetic bone substitute works like a bone graft that maintains the space of the bone defect, repairing the defect or even augmenting the alveolar ridge. Plenum® Oss<sub>hp</sub> is partially resorbed by the body and replaced by newly formed tissue during the process of bone tissue repair or regeneration of bone tissue.

## 5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The macro and microstructure, surface area, porosity percentage, and physicochemical characteristics of Plenum® Oss<sub>hp</sub> promote efficient osteoconductive properties, high hydrophilicity, and adaptability near the host bone, allowing the diffusion of biological fluids and cell adhesion to the product.

As a resorbable synthetic bone graft, Plenum® Oss<sub>hp</sub> is designed to maintain the space of the bone defect while promoting the growth of new bone tissue. The constituent bioceramic phases are partially reabsorbed by the body and replaced by neofomed tissue during bone tissue repair or regeneration.

Overall, Plenum® Oss<sub>hp</sub> offers a safe and effective solution for the filling and/or reconstructing of bone defects, including bone defect space maintenance, defect repair, alveolar ridge augmentation or maxillary sinus elevation.

## 6. INSTRUCTIONS FOR USE

- Proper preparation of the site (host bone) at which the Plenum® Oss<sub>hp</sub> will be applied is of the utmost importance. The device must fill the entire area to be treated, and any bleeding must be controlled so that the bone substitute adapts perfectly to the host bone and provides the expected results;
- The surgeon shall use only appropriate and sterile surgical instruments during procedures for manipulation and implantation of Plenum® Oss<sub>hp</sub>;
- Before the surgical procedure, the surgeon (physician or dentist) shall carry out rigorous pre-operative planning through clinical examinations and imaging (radiographs or computed tomography). They shall also consider the need for a detailed history and complementary exams regarding the general health of the patient (complete hemogram, coagulogram, calcium dosage, etc.) since this information directly influences the body's biological response to the use of the bone substitute, Plenum® Oss<sub>hp</sub>, and then define a surgical plan in advance, selecting the granules, blocks, wedges, cylinders sizes and quantity most appropriate to the condition of the patient;
- The surgical technique to be used must be indicated by the surgeon using their clinic and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique should not be applied to all patients. Surgical protocols must be carried out according to references and previous experiences of the surgeon, always considering the most appropriate choice of granule size, blocks, wedges, cylinders models and quantity of the device, the technique, the grafting sequence, and the use of membranes or meshes whenever necessary. These decisions must be based on conventional and established therapeutic techniques for guided tissue regeneration and bone grafting;
- Postoperative care and proper hygiene/cleanliness of the operated site are extremely important. The patient must be instructed regarding existing risks until there is complete osseointegration of the implant;
- The patient must be properly oriented about care during the postoperative phase. The patient's ability and willingness to follow the instructions of the surgeon is one of the most important aspects of a surgical procedure;
- Elderly patients, patients with mental health problems, individuals with chemical dependency, or any other difficulty in adhering to postoperative recommendations may represent a higher risk of failure of this device since these patients may ignore instructions for this procedure;
- The surgeon must adopt the proper clinical/surgical planning described in the specific literature for guided bone regeneration and bone grafting techniques for immediate and subsequent pre- and postoperative periods;
- The surgeon must clearly guide the patient, who must understand the importance of monitoring during the postoperative period, in which the stability and integrity of the device must be verified, as well as the incorporation and/or reabsorption at the appropriate period;
- To protect the patient and their surgeon, it is important that the correct and complete traceability of the Plenum® Oss<sub>hp</sub> used to be guaranteed by noting the device code and batch code in the patient's records or by attaching to the identification and traceability labels that come with the device to these records.
- Carefully read the instructions for use, recommendations, and guidelines before using the Plenum® Oss<sub>hp</sub> synthetic bone substitute;
- Verify the device's expiration date before use. Do not use devices if the packaging has been damaged or if the expiration date has expired;
- Plenum® Oss<sub>hp</sub> is supplied ready for use. It is available in granule format, in different quantities and in three particle size ranges (250 - 500 µm, 500 - 1000 µm and 1 -2 mm). It is also available in the form of blocks, wedges and cylinders of different sizes. Plenum® Oss<sub>hp</sub> does not require prior hydration, but hydration in saline solution or the patient's own blood is recommended to improve the efficiency of the product's handling.
- The surgeon (physician or dentist) is responsible for choosing the device model with the most appropriate characteristics and dimensions after an individual clinical analysis of each case;
- The techniques employed for the use and application of Plenum® Oss<sub>hp</sub> depend on the preference of the surgeon (physician or dentist), and they are responsible for choosing the therapeutic approach that is used, model and dimensions of the device, application technique, as well as the criteria for monitoring and evaluating the results of surgery;

## Handling and Manipulation of Plenum® Oss<sub>hp</sub> - Granule

## Handling and Manipulation of Plenum® Oss<sub>hp</sub> - Block / Wedge / Cylinder

1. Remove the blister containing the product from the outer carton.

2. Open the blister on a sterile surgical drape and place the primary packaging (case) plus product on the operating table.

3. Open the sterile glove case by turning the lid anticlockwise

4. Once the case is open, remove the block using tweezers or your thumb and forefinger.

**ATTENTION:** The figures are merely illustrative. They do not represent the actual dimensions and characteristics of the device.

## Customization of the Plenum® Oss<sub>hp</sub> Block / Wedge / Cylinder and adaptation to the receiving bed:

The Plenum® Oss<sub>hp</sub> Blocks / Wedges / Cylinders are medical devices that can be adapted and adjusted, but the mechanical conformation limits must be respected according to the table below.

Minimum thickness for screw fixing	3 mm	
Minimum graft edge thickness for screw fixation	3 mm	
Maximum number of fixing screws per commercial model	<p>Minimum: 1 screw (for fixings in the center of the bone graft);</p> <p>Maximum: 2 screws (for fixation at the ends of the bone graft); Note: The fixation screws must be positioned in such a way that they respect the minimum distance of 3 mm between them and between the ends of the graft (edge of the graft).</p>	

Customize the Plenum® Oss<sub>hp</sub> Blocks / Wedges / Cylinders according to the anatomy of the recipient bed with the aid of sterile straight-piece surgical/prosthetic drills and discs (examples: 702 drills, Maxicut line models and diamond discs);

Finish properly, leaving no sharp edges or corners to avoid perforating the flap;

Adapt the block to the recipient bed without leaving any steps, leveling all sides of the graft to the bone defect.

Maximum recommended torque limit for screw fixation: 10 N.cm.

5. Fitting Plenum® Oss<sub>hp</sub> over the surgical bed and stabilizing with a screw.

6. Final appearance of the Plenum® Oss<sub>hp</sub> adapted and fixed to the surgical bed.

**ATTENTION:** The figures are merely illustrative. They do not represent the actual dimensions and characteristics of the device.

## 7. PRECAUTIONS

- Surgeons must not use the Plenum® Oss<sub>hp</sub> until they have carefully read the instructions for use and the information provided on the packaging. All information on the use of bone grafts and the warnings, precautions, and recommendations included in these instructions for use must be taken into consideration. A failure to observe these points may result in complications during procedures;
- The Plenum® Oss<sub>hp</sub> must be exclusively handled by properly trained surgeons (physicians or dentists) qualified for this procedure, with experience and knowledge of the techniques used in guided tissue regeneration and use of bone substitutes;
- It is recommended that preoperative examinations be requested according to the knowledge of the patient's condition;
- Device labels include the date of manufacture, device code, batch code, ANVISA (Brazilian National Healthy Regulatory Agency) device registration number, technical device manager, and company information, in addition to the respective symbology that allows the device to be correctly traced;
- The device shall be kept in the original packaging until it is used and shall only be handled in a sterile surgical field;
- Any device that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item Disposal of Device;
- The surgeon should understand and plan all steps involved in the procedure and verify the number of devices they intend to use. Always ensure that the device is available in quantities that exceed those expected to be used;
- Surgeons must be aware of the patient's clinical conditions. During preoperative examinations, clinical examinations and imaging of the site at which Plenum® Oss<sub>hp</sub> will be placed must be rigorously verified. Surgeons must consider the patient's general clinical condition and analyze local or systemic factors that may interfere with the process of healing hard and soft tissue, and bone substitute incorporation and reabsorption;
- Before starting the surgery, the surgeon must submit the patient to a rigorous oral antisepsis process to prevent the device from coming into contact with contaminating substances;
- Patients must be warned that a failure to comply with instructions for the post-operative period may lead to failure of the surgical procedure, requiring additional surgery to adjust or remove substitute bone;
- The use of this device by patients with inadequate oral hygiene, clinical conditions involving hematological disorders, and other uncontrolled systemic alterations are restricted.
- The use of any remnant Plenum® Oss<sub>hp</sub> materials is prohibited. Any remnant device residuals shall be disposed of after the end of the surgery, following health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- Improper selection of the bone substitute that will be used, as well as mistakes in the indicated use, handling, preparation of bone or soft tissues, and installation of the device, may result in damage to the device's structure, which may contribute to failure;
- The use of inappropriate techniques and procedures with this bone substitute and inadequate surgical conditions may cause injuries to the patients and compromise results for tissue repair and regeneration;
- The primary and most likely causes of the lack of incorporation and reabsorption of Plenum® Oss<sub>hp</sub> with the consequent compromising of

surgery include infections, localized or systemic diseases, inappropriate use of membranes or meshes insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgery site, contamination of the device pre- and/or trans-operative or postoperative whenever there is an exposure of the bone substitute in the oral cavity, mechanical instability of granules, blocks, wedges, cylinders, bone necrosis, lack of professional training, technical error in installation of the device, and a failure to adhere to postoperative instructions on the part of patients. The surgeon must consider all these aspects during preoperative planning;

- The Plenum® Oss<sub>hp</sub> has not been tested on lactating women. Therefore, it is unknown whether the metabolites contained in Plenum® Oss<sub>hp</sub> are passed on to breast milk.
- The Plenum® Oss<sub>hp</sub> is merely intended for use as an auxiliary material in bone replacement. Exclusive use of the Plenum® Oss<sub>hp</sub> alone does not guarantee the success of the treatment. The success of therapy depends on the intervention technique used and the patient's reaction;
- The success of a guided tissue regeneration or bone reconstruction surgical procedure using Plenum® Oss<sub>hp</sub> associated with other substitutes, bone grafts, or membranes is related to the correct indication and selection of the device, use of the correct surgical technique, preparation, the stability of the device, the period and quality of the incorporation of biomaterial(s)/ host since these extrinsic factors can be affected by numerous biological, chemical and biomechanical factors. These factors may cause failure and early loss of the device and compromise the surgical procedure. Following the recommendations, contraindications, precautions, warnings, and other observations contained in these instructions for use are, therefore, essential for the success of the procedure.
- Patients must be informed of all limitations, risks, complications, warnings, and potential adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques involving the use of synthetic bone substitutes;
- Patients must be informed of the risks involved in procedures using bone substitutes in which the operation site will be subject to stress and direct mechanical forces, for example, exposure to constant deformation and/or tension;
- Patients must be informed of the risks involved in the application of excessive mechanical load or trauma that affects the operated region, which may imply failure, loss of mechanical stability, exposure, and displacement of the synthetic bone substitute's parts;
- It is of the utmost importance that the surgeon explains the entirety of the above topics to their patient and that patients, in turn, have a clear and complete understanding of these procedural aspects.

## 8. WARNINGS

- The Plenum® Oss<sub>hp</sub> is supplied in sterile condition, sterilized using Gamma Irradiation. Do not use if the package is opened or damaged or if expiration date has been exceeded.
- This is a SINGLE USE product, so reprocessing, resterilization and reuse are prohibited and may pose risks of failure and infection to the patient. If the product is exposed to any of these situations, it will result in contamination and permanent damage to the product and/or cross-infection of the patient;
- Do not overfill defects.
- Do not leave defect open.
- Do not compromise blood supply to the defect area.
- The device should be secured to prevent motion and migration, use in areas where the graft can be adequately contained.

## 9. CONTRAINDICATIONS

Surgeons shall prudently assess the patient's condition before the surgery and use of the Plenum® Oss<sub>hp</sub>, checking for any contraindications reported below. Contraindications may include, but are not limited to:

- Patients with an insufficient quantity and general quality of bone or soft tissue to provide primary stability of the synthetic bone substitute;
  - Patients with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the bone substitute due to inefficient tissue repair or regeneration;
  - Use of bone substitute in locations subject to stress, direct mechanical forces, or granules, blocks, wedges and cylinders exposed to constant deformation and/or tension forces;
  - Patients that do not have enough soft tissue to provide complete coverage of the bone substitute;
  - Inflammation and/or infection near or at the surgical site;
  - Recent history of systemic or localized infection;
  - Fever;
  - Leukocytosis;
  - Patients that present signs of allergy or hypersensitivity to the composition of the raw materials used in the manufacturing of Plenum® Oss<sub>hp</sub> (see "DEVICE COMPOSITION");
  - Patients with developmental disorders of the skeletal system;
  - Pregnancy and lactation period;
  - Use in patients with immunological changes and systemic diseases that result in poor healing or the absence of integrity of existing tissue at the surgical site;
  - Patients with vascular disorders or coagulation disorders;
  - Calcium metabolism disorders;
  - Treatment with steroids and other drugs that have interfered with calcium metabolism;
  - Alcoholic patient or drug user who cannot correctly follow postoperative instructions;
  - Patients with mental illness or unable to follow the recommendations and care of the postoperative phase;
  - Patients unable to maintain proper hygiene at the surgical site;
- Immunosuppressive pathologies;
- Patients undergoing radiotherapy or using medications that compromise the vascularization of bone and soft tissues;
  - All contraindications included available dental and scientific literature for resorbable bone substitutes.

## 10. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved during all surgical procedures. Every surgical procedure presents risks and the possibility of complications, such as infections, hemorrhaging, allergic reactions to drugs, and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum®

Oss<sub>hp</sub> may include, but are not limited to:


- The surgical procedure involved in the installation of bone grafts or substitutes may cause complications during transitory and postoperative periods, such as hemorrhaging, pain, edema, abscesses, fistula, dehiscence, injury or compression of nerves (partial or total sensory or motor paresthesia that is transient or permanent), soft tissue injury, local or systemic infection, and other typical complications for placement procedures for a bone graft or substitutes;
- Tissue sensitivity at the surgical site;
- Pain or discomfort due to the presence of the bone substitute or depending on the surgical procedure;
- Lack of incorporation of the bone substitute or even delayed healing, which may result in device failure;
- Loss of stability, displacement of bone substitute granules, blocks, wedges and cylinders resulting from failure in the incorporation and early partial reabsorption, as well as inadequate installation of the device;
- Displacement or exposure of the bone graft granules, blocks, wedges and cylinders in the oral cavity due to use at sites that may be subject to stress, direct mechanical force, exposure to constant deformation and/or tension, and inadequate fixation of the device;
- Allergic reaction to the raw materials used in components;
- Intense fibrotic tissue reaction around the surgical site;
- Bone necrosis, insufficient vascularization, reabsorption of the bone graft, and inadequate bone formation, which may result in loss of the grafted area;
- Premature dehiscence in sutures, the tension in the coaptation of surgical edges, or early exposure of Plenum® Oss<sub>hp</sub> and the grafted area may alter cystic processes and fail the surgical procedure;
- The patient must be informed and understand that there is a possibility of complications in all surgical. The surgeon must notify the patient of these complications and the risks associated with guided tissue regeneration surgeries and bone grafting techniques involving resorbable synthetic bone grafts.

## 11. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. **Please check the physical labeling for the symbols applied to the device.** The symbols express information related to the identification and use of the device.

**TABLE 2.** Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using ethylene oxide

	Sterilized using Irradiation
	Sterilized using vaporized hydrogen peroxide
	Non-sterile
	Do not re-sterilize
	Do not use if package is damaged
	Single sterile barrier system
	Double sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Keep away from sunlight
	Keep dry
	Lower limit of temperature
	Upper limit of temperature
	Temperature limit
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist.  (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)
	Consult the instructions for use in electronic format (QR code)
	ANVISA device registration number

**NOTE:** The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

## 12. STERILIZATION

The Plenum® Oss<sub>hp</sub> is provided in sterile conditions and sterilized using Gamma Irradiation;

The Plenum® Oss<sub>hp</sub> is single-use only and must not be reprocessed, resterilized, or reused;

- Before use, check the integrity of the packaging and the expiration date. Do not use the device if the package is damaged or if the expiration date has expired;
- Date of manufacture, sterilization expiration date, code, batch number, and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in appropriate and aseptic environments;
- The packaging opening for the device's use must only be done by qualified and trained personnel to perform this procedure.

## 13. EXPIRATION DATE

- Indicated on device labeling.
- Before use, check the expiration date, do not use the device if it has expired.

## 14. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Oss<sub>hp</sub> and comply with health surveillance requirements, the surgeon or other staff must keep the information regarding the device used in the patient's medical records. Additionally, this information must also be passed on to the device distributor and patient to complete the traceability cycle. The required information for traceability is related to the device used, surgical procedure, and patient, as follows:

- Surgeon's name;
- hospital/clinic/office identification;
- manufacturer and distributor identification;
- Date of surgery;
- Device code used;
- Batch code;
- Quantity used;
- ANVISA device registration number.

The responsible surgeon and their team must use the information in the three (03) traceability labels supplied with the Plenum® Oss<sub>hp</sub>, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that they have information about the device used in their procedure.

Labels contain the following information necessary for the traceability of the device used.

- Manufacturer identification;
- Device name or commercial model;
- Device code;
- Batch code;
- ANVISA device registration number.

Traceability information is necessary for notification by the health service and/or by the patient to the regulatory authority, and to the manufacturer, when complaints, technical complaints and/or adverse events occur, for the carrying out of the appropriate investigations.

## 15. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the device if the packaging is damaged or if the expiration date has expired. If the packaging is damaged, the device must be considered non-sterile and discarded;
- The device must be stored and transported in a clean environment (upper limit temperature of 60°C / 140°F), away from sources of heat or humidity, sunlight and the action of the weather, to prevent any damage or change in the packaging or its physical form that makes it impossible to use;
- The device must be stored on shelving or drawers/cabinets to allow it to be easily viewed and handled by the operator. Do not store the device near lamps so that the information contained in the packaging is not affected;
- The device shall only be opened at the time of use in a sterile surgical field;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

## 16. DISPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After using the device or removing it from the patient, the synthetic bone graft (Plenum® Oss<sub>hp</sub>), that are not subject to analysis or subsequent studies, must be disposed of by the health service establishment (hospitals, offices, and clinics);
- It is recommended that the Plenum® Oss<sub>hp</sub> be completely disassembled to ensure that it cannot be reused. The health service institution/establishment is solely responsible for disassembling the device;
- After disassembly, the Plenum® Oss<sub>hp</sub> must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards (hospitals, offices, and clinics) and applicable environmental legislation;
- Remnant residuals of the Plenum® Oss<sub>hp</sub> not used during surgery must also be discarded following the guidelines above.

## 17. CUSTOMER SERVICES

- If the device presents any adverse event, with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, and batch, and patient data; directly notifying the Plenum Customer Relations Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and properly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Telephone: +55 11 3109- 9000; Email: customerservice@plenum.bio;
- ANVISA also provides a channel available to report incidents, adverse events and technical complaints related to the use of products and services subject to health regulations, which can be found on the website of the National Health Regulatory Agency - ANVISA at the electronic address [www.anvisa.gov.br](http://www.anvisa.gov.br), using the NOTIVISA link.
- FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

## 18. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health professional and addressed to regulatory authority standards (see the Customer Services section) when it is observed that issues involve a risk to the patient's health.

The warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing, re-sterilization, or reuse of the device by the user;
- Improper use or use other than in a surgical environment;
- Device not used according to instructions for use.

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